



Constant temperature incubator shaker OD Monitor

| | |
|---|-----|
| ■ Model Selection Guide..... | 010 |
| ■ Constant temperature incubator shaker with NewMax drive | |
| DWMMax VBR-104..... | 012 |
| DWMMaxDwMax MBR-034P/034/032P/032/104P..... | 013 |
| ■ Compact size constant temperature incubator shaker for Well plate and Microtube | |
| Maximizer MBR-022UP..... | 015 |
| Bioshaker MBR-024..... | 016 |
| ■ Medium size constant temperature incubator shaker for Well plate | |
| Bioshaker MBR-420FL..... | 017 |
| ■ Small size constant temperature incubator shaker | |
| Bioshaker BR-21/22/23 Series..... | 018 |
| ■ Medium size constant temperature incubator shaker | |
| Bioshaker BR-53FP..... | 020 |
| Bioshaker BR-40/41/42/43Series..... | 022 |
| ■ Large size constant temperature incubator shaker | |
| Bioshaker BR-300LF/BR-300LF LED irradiation..... | 024 |
| Bioshaker Double platform BR-300LF..... | 026 |
| Bioshaker BR-180LF/BR-180LF-70RT..... | 028 |
| Bioshaker GBR-200/300..... | 030 |
| Bioshaker BR-3000LF/Double platformBR-3000LF..... | 032 |

| | |
|--|-----|
| ■ Non-contact turbidimeter for Shaking culture | |
| OD-Monitor A&S/B&L/C&T..... | 034 |
| OD-Monitor Turbidostat System A/B..... | 035 |
| ■ Digital colorimeter | |
| Easy operation OD Monitor miniphoto 518R..... | 036 |
| ■ LED irradiation unit | |
| LED irradiation unit LC-LED450W..... | 037 |

Selection guide

Constant temperature incubator shaker

CO₂ incubator shaker

Shaker

Mixer Rotator Stirrer

Bead beater homogenizer Ultrasonic homogenizer

Aluminum block bath Minisize bath

Water bath Shaking water bath Immersion cooler

Hybridization oven Constant temperature chamber

Centrifugal concentrator Cold trap

Freeze dryer

Electrophoresis and blotting apparatus

Constant temperature water circulating system [Chiller]

Appendix

NewMax drive

DWMMax MBRSeries



DWMMax VBR-104



Maximizer MBR-022UP



Bioshaker MBR-024



Maximizer MBR-420FL



for Deepwell plate

Small size

Low temp.
FP TypeMiddle temp.
UM TypeHigh temp.
FH Type

Medium size

Low temp.
LF/FL/FP TypeMiddle temp.
FM TypeHigh temp.
FH Type

Large size

Bioshaker BR-300LF Series



Bioshaker BR-180LF Series



Bioshaker GBR Series



Bioshaker BR-3000LF Series



Constant temperature incubator shaker series

Widely used for Cultivation of microbes such as Yeast, E.Coli. etc. that require temp. control.

With NewMax drive

| Page | Model | Door type | Cooling method | Temperature range |
|-------|-------------|----------------|-----------------|--------------------------------------|
| P.013 | ① M-BR-034P | • lift-up door | Peltier element | • 7°C below Room temperature to 60°C |
| | ① M-BR-034 | | N/A | • 7°C above Room temperature to 60°C |
| | ② M-BR-032P | | Peltier element | • 15°C to 60°C |
| | ② M-BR-032 | | N/A | • 7°C above Room temperature to 60°C |
| P.012 | ① M-BR-104P | | Peltier element | • 7°C below Room temperature to 60°C |
| | ③ V-BR-104 | | N/A | • 8°C above Room temperature to 60°C |

For Well plate and Microtube (Compact size & Medium size)

| Page | Model | Door type | Cooling method | Temperature range |
|-------|--------------|---------------------|-----------------|--------------------------------------|
| P.015 | ④ M-BR-022UP | • Lift-up door | Peltier element | • 15°C to 60°C |
| P.016 | ⑤ M-BR-024 | | N/A | • 5°C above Room temperature to 65°C |
| P.017 | ⑥ M-BR-420FL | • Single swing door | ✓ | • 4°C to 70°C |

Bioshaker (Small size)

| Page | Model | Door type | Cooling method | Temperature range |
|-------|--------------|-------------------------------------|-----------------|---------------------------------------|
| P.018 | ⑦ BR-21FP-MR | • Single swing door (Leftward open) | Peltier element | • 15°C to 60°C |
| | ⑦ BR-22FP-MR | | | |
| | ⑦ BR-23FP-MR | | | |
| P.018 | ③ BR-21UM-MR | • lift-up door | N/A | • 5°C above Room temperature to 70°C |
| | ③ BR-22UM-MR | | | |
| | ③ BR-23UM-MR | | | |
| P.018 | ⑦ BR-21FH-MR | • Single swing door (Leftward open) | | • 5°C above Room temperature to 100°C |
| | ⑦ BR-22FH-MR | | | |
| | ⑦ BR-23FH-MR | | | |

Bioshaker (Medium size)

| Page | Model | Door type | Cooling method | Temperature range |
|-------|--------------|--|-----------------|---------------------------------------|
| P.022 | ⑨ BR-40LF | • Swing lift-up door | ✓ | • 4°C to 70°C |
| | ⑥ BR-41FL-MR | | | |
| | ⑥ BR-42FL-MR | | | |
| | ⑥ BR-43FL-MR | | | |
| P.020 | ⑥ BR-53FP | • Double hinged door (Rightward/Leftward open) | Peltier element | • 15°C to 55°C (*) |
| P.022 | ⑥ BR-41FM-MR | | N/A | • 5°C above Room temperature to 70°C |
| | ⑥ BR-42FM-MR | | | |
| | ⑥ BR-43FM-MR | | | |
| P.022 | ⑥ BR-43FH-MR | | | • 5°C above Room temperature to 100°C |

Bioshaker (Large size & Double platform & LED irradiation)

| Page | Model | Door type | Cooling method | Temperature range |
|-------|---------------------------|--|----------------|-------------------|
| P.024 | ⑩BR-300LF | ●Single swing door (Leftward open) | ✓ | ●4°C to 70°C |
| P.026 | ⑩LED irradiationBR-300LF | | | |
| | ⑩Double platformBR-300LF | | | |
| P.028 | ⑪BR-180LF | ●Clamshell split door (Divided into Upper and Lower, Upward opening and Front opening) | | ●4°C to 80°C |
| | ⑪BR-180LF-70RT | | | |
| P.030 | ⑫G-BR-200 | ●Swing lift-up door | | ●4°C to 80°C |
| | ⑫G-BR-300 | | | |
| P.032 | ⑬BR-3000LF | ●Double swing door | ●4°C to 70°C | |
| | ⑬Double platformBR-3000LF | | | |

LED irradiation unit

LED irradiation unit that is optimum used with Middle sized Bioshaker in combination.

-->P.037




Compact size

Small size

Medium size

Large size








150cm



| | Shaking motion | Speed range | Features | Applications | Page |
|---|-------------------------------|----------------|--|---|-------|
| ● | •Horizontal eccentric shaking | •200~1600r/min | •Well plate×4 •Microtube rack×4 | •Cold-Shock expression for E.coli | P.013 |
| ● | | | •Well plate×2 •Microtube rack×2 | •Shaking Culture with Deep well plate or Microtube. | |
| ● | •Orbital | •50~500r/min | •15mL disposable tube ×16 •50ml disposable tube ×8 | •Cold-Shock expression for E.coli | P.012 |
| ● | | | •100m Erlenmeyer flask×6 (Each requires optional accessories.) | •Shaking Culture with Deep well plate or Microtube. | |
| ● | | | | •Enzyme digestion, Biodegradation test. | |
| ● | | | | •Cultivation of Microbes such as E.coli. | |

| | Shaking motion | Speed range | Features | Applications | Page |
|---|-------------------------------|----------------|---------------------------------|---|-------|
| ● | •Horizontal eccentric shaking | •300~2500r/min | •Well plate×2 •Microtube rack×2 | •Shaking culture with Deep-well plate and Microtube rack. | P.015 |
| ● | | •100~1200r/min | •Well plate×4 •Microtube rack×2 | •ELISA and Phage display method. | P.016 |
| ● | | •200~1500r/min | •Well plate ×12 •Microplate ×24 | •Possible for shaking culture with 96-deep well plate. | P.017 |
| | | | | •ELISA and Phage display method. | |

| | Shaking motion | Speed range | Shaking width | Features | Applications | Page |
|---|-----------------------|--------------|---------------|---|--|-------|
| ● | •Reciprocal | •20~300r/min | •25mm | •500m Erlenmeyer flask×4 •2 Erlenmeyer flask×1 | •Cultivation of yeast and insect cell. | P.018 |
| ● | •Orbital | | | | | |
| ● | •Reciprocal / Orbital | | | | | |
| ● | •Reciprocal | | | | •Cultivation of Microbes such as E.coli. | P.018 |
| ● | •Orbital | | | | | |
| ● | •Reciprocal / Orbital | | | | | |
| ● | •Reciprocal | •20~300r/min | •25mm | •500m Erlenmeyer flask×8 •5 Erlenmeyer flask×1 | •Cultivation of Thermophile and Hybridization. | P.018 |
| ● | •Orbital | | | | | |
| ● | •Reciprocal / Orbital | | | | | |

| | Shaking motion | Speed range | Shaking width | Features | Applications | Page |
|---|-----------------------|--------------|---------------|---|--|-------|
|  | ●Reciprocal / Orbital | ●20~200r/min | ●25mm | ●500m Erlenmeyer flask×6 ●2 Erlenmeyer flask×2 | ●Cultivation of Microbes such as E.coli. ●Cultivation of yeast and insect cell. | P.022 |
|  | ●Reciprocal | ●20~300r/min | | ●500m Erlenmeyer flask×8 ●5 Erlenmeyer flask×1 | | |
|  | ●Orbital | | | | ●Cultivation of Microbes such as E.coli. | P.022 |
|  | ●Reciprocal / Orbital | | | | ●Cultivation of Thermophile and Hybridization. | P.022 |
|  | ●Reciprocal | | | | | |
|  | ●Orbital | | | | | |
|  | ●Reciprocal / Orbital | | | | | |

(*) The value at 25°C of the rear surface intake temp.

| | Shaking motion | Speed range | Shaking width | Features | Applications | Page |
|---|-----------------------|------------------|---------------|--|--|-------|
| ● | •Reciprocal / Orbital | •25~250r/min | •10~50mm | •500m Erlenmeyer flask×15 •5 Erlenmeyer flask×3 | •Cultivation of microbes such as E.Coli. etc. | P.024 |
| ● | | •25~160r/min | | •500m Erlenmeyer flask×30 •1 Erlenmeyer flask×18 | •AGP test | P.026 |
| ● | | •25~250r/min (*) | | •500m Erlenmeyer flask×24 •5 Erlenmeyer flask×4 | •Cultivation of microbes such as E.Coli. etc. | P.028 |
| ● | •Reciprocal | •25~250r/min | •70mm | | •Cultivation of actinomycete. | |
| ● | •Reciprocal / Orbital | •25~200r/min | •25/50mm | •For as a laboratory bench. Stackable up to Two stage. | •Cultivation of E. coli, yeast and insect cells. | P.030 |
| ● | | •25~300r/min | | •500m Erlenmeyer flask×24 •5 Erlenmeyer flask×4 | •Cultivation of Thermophile and Hybridization. | |
| ● | | •25~250r/min | •10~50mm | •500m Erlenmeyer flask×32 •5 Erlenmeyer flask×4 | | P.032 |
| ● | | •25~160r/min | | •500m Erlenmeyer flask×64 •1 Erlenmeyer flask×36 | •Cultivation of E. coli, yeast and insect cells. | |

(*) When the width is 25mm, only orbital can be 25 to 400r/min.

OD-Monitor Series

Non-contact turbidimeter that realizes to measure the turbidity of sample at OD600 while shaking in noncontact.
Available 3 models for Large and Small Erlenmeyer flask and Test tube.

-->P.034



● Protuberances not included in Dimensions. ● Vessels of photo not included.

We contribute to the development of research and industry.
[2019-2020 General Catalog] **TATEC**

Selection guide

Constant temperature incubator shaker

CO₂ incubator shaker

Shaker

Mixer Rotator Stirrer

Bead beater homogenizer Ultrasonic homogenizer

Aluminum block bath Minisize bath

Water bath Shaking water bath Immersion cooler

Hybridization oven Constant temperature chamber

Centrifugal concentrator Cold trap

Freeze dryer

Electrophoresis and Blotting apparatus

Constant temperature circulating system [Chiller]

Appendix

DwMax VBR-104

**Our proprietary technology "NewMax drive" enables excellent aeration effectiveness Max. 500r/min.
Optimum for bacteriophage culture and small-scale cultivation.**




Our proprietary technology "NewMax drive" enables a great stirring effectiveness to preventing the wobble.

The shaker generates synergistically the wobble by its vibration and the wobbling of platform by inertia. The wobbling sometimes spoils the effectiveness of shaking and stirring. NewMax drive realizes the shaker less wobbled. For example, we obtained the effectiveness dramatically the variation between wells in the cultivation using deep wells improved (see the page 014).

| Model | VBR-104 |
|-----------------------------|---|
| Temperature range (*1) | 8°C above RT (*2) to 60°C |
| Temp. control accuracy (*3) | ±0.5°C to ±1.0°C |
| Shaking motion | Orbital shaking |
| Speed range | 50 to 500r/min |
| Platform size | 240 x 220mm (Max. inside height: 168mm) |
| Ambient temp. range | 5°C to 30°C |
| Heating method | Heater 80W |
| Safety devices/functions | Fuse, Self-diagnosis function for temp. control error, Braking when hood open, Motor overload/Over-voltage protection |
| Dimensions (W x D x H) | 375 x 445 x 260mm |
| Maximum load | 1kg (Platform included) |
| Weight | Approx. 20kg |
| Power supply | AC100 to 240V/2A |
| Standard accessories | 1 x Stick sheet |

(*1)It may not reach the maximum temperature by the surroundings. (*2)RT (Room temperature) has been set 25°C as the standard. (*3)The value at the ambient temperature 25°C.

Dimensions



Features

- Excellent aeration effectiveness.
- Suitable for small Erlenmeyer flask, Centrifuge tube and Microtube.
- Compact bench-top unit is easy to handle.

Applications

- Cultivation of Bacteriophage
- Cultivation of Microbes such as E.coli
- Solubility test for Refractory


Easy to detach the platform. Plenty variations.

Comes with Stick sheet. Spring net shaking platform, Universal shaking platform and others are available as an option.

Spring net shaking platform
MR-2420 (Option)



Optional parts : Platforms



| Product Name / Model | Capacity of Vessels |
|--|--------------------------|
| Shaking platform for 1.5mL Microtube V-BR-1140 | 1.5/2mL Microtube x40 |
| Horizontal shaking platform for 15mL centrifuge tube V-BR-1816 | 15mL disposable tube x16 |
| Diagonal shaking platform table for 15mL centrifuge tube V-BR-1816S | 15mL disposable tube x16 |
| Horizontal shaking platform for 50mL centrifuge tube V-BR-3508 | 50mL disposable tube x8 |
| Diagonal shaking platform table for 50mL centrifuge tube V-BR-3508S | 50mL disposable tube x8 |
| Shaking platform for 50mL Erlenmeyer flask V-BR-0050 | 50mL Erlenmeyer flaskx6 |
| Shaking platform for 100mL Erlenmeyer flask V-BR-0100 | 100mL Erlenmeyer flaskx6 |
| Spring net shaking platform MR-2420 | 100mL Erlenmeyer flaskx6 |

DwMax MBR-034P/034/032P/032/104P

Powerful, Uniform and Silent! Our proprietary technology "NewMax drive" enables Excellent aeration effectiveness of 96-well plate, microtube and disposable centrifuge tube.

Experimental data --> P.014



MBR-034P
Deep well plate



MBR-032P
Microtube rack (option)



MBR-104P
Shaking platform (option)

Features


- Optimum for using 96-well plate and stirring for microtube [Model 03]
- Optimum for stirring of Microtube and/or Centrifuge tube [Model 10]
- Temperature control below R.T. by peltier cooling [P Type]

Our proprietary technology "NewMax drive™" enables a great stirring effectiveness to preventing the wobble.

The shaker generates synergistically the wobble by its vibration and the wobbling of platform by inertia. The wobble sometimes weak the effects for shaking and stirring. NewMax drive realizes Shaker less wobbled. For example, the effect of dramatically the variation between wells in the cultivation using well plate rack improved has been obtained (see the next page).

Applications

- Shaking Culture with Deep well plate or Microtube
- Cold-Shock expression for E.coli [Model-032], ELISA etc.
- Enzyme digestion, Biodegradation test, Metabolomic analysis etc.



| Capacity | 4 x Well plate (*1) rack or 4 x Microtube rack | | 2 x Well plate or 2 x Microtube rack | | 1.5 or 2 mL microtube rack and 15 or 50mL disposable centrifuge tube by changing the platform. |
|----------------------------|---|---------------------------|--------------------------------------|---------------------------|--|
| Model | MBR-034P | MBR-034 | MBR-032P | MBR-032 | MBR-104P |
| Temperature range | 7°C below RT (*2) to 60°C | 7°C above RT (*2) to 60°C | 15°C to 60°C | 7°C above RT (*2) to 60°C | 7°C below RT (*2) to 60°C |
| Temp. control accuracy | ±0.5°C to ±1.0°C (*3) | | | | |
| Shaking motion/Speed range | Horizontal eccentric shaking, 200 to 1600r/min | | | | Orbital, 50 to 500r/min |
| Platform size | - | | | | 240 x 220mm (Max. inside height : 100mm) |
| Ambient temp. range | 5°C to 30°C | | | | |
| Heating method | Peltier element | Heater 70W | Peltier element | Heater 70W | Peltier element |
| Cooling method | Peltier element | - | Peltier element | - | Peltier element |
| Other functions | Digital display of temperature and speed, 1 x Memory function of temperature and speed. | | | | |
| Safety devices/functions | Fuse, Self-diagnosis function for temp. control error, Braking when hood open, Motor overload/Over-voltage protection | | | | |
| Dimensions (W x D x H) | 405 x 445 x 230mm | 405 x 402 x 230mm | 365 x 400 x 230mm | 365 x 347 x 230mm | 405 x 445 x 230mm |
| Weight | Approx. 18kg | Approx. 16kg | Approx. 15kg | Approx. 13kg | Approx. 18kg |
| Power supply | AC100 to 240V/1.5A | | | | |
| Standard accessories | - | | | | |

(*1) Equipped with the platform for Well plate. Since designed for 96-well plate, the agitating effectiveness will be different by using different well sizes such as 24-well and 384-well plate.

(*2) RT (Room temperature) has been set 25°C as the standard. (*3) The value at the ambient temperature 25°C.

Optional parts

| Product Name / Model | Capacity of Vessels |
|-----------------------------|-----------------------|
| Microtube rack EM-0524 | 0.5mL Microtube x24 |
| Microtube rack EM-1524 | 1.5/2mL Microtube x24 |
| Microtube rack EM-1515 (*1) | 5mL tube x12 |

| Product Name / Model | Capacity of Vessels |
|--|--------------------------|
| Shaking platform for 1.5mL Microtube V-BR-1140 | 1.5/2mL Microtube x40 |
| Horizontal shaking platform for 15mL centrifuge tube V-BR-1816 | 15mL disposable tube x16 |
| Horizontal shaking platform for 50mL centrifuge tube V-BR-3508 | 50mL disposable tube x8 |

(*1) The suitable 5mL Microtubes are WATSON's (Japan) and Argos's (USA) not Eppendorf's (Germany). The Microtube racks are for 03 models (034P/034/032P/032) and the shaking platform is for 10 model (104P).

DWMax experimental data

History for shaking culture with Deep well plate

Deep well plate (hereinafter, DWP) is a container originally used to store samples; However, when the method of High-throughput screening (HTS) was developed around 2000, it was started being used for not only Storage but also other applications such as Cultivation and Reaction. At that time, we had a customer's inquiry "Do you have any ideas how to perform a shaking culture with DWP effectively?" Then, we managed to developed "Bioshaker® MBR-024 for DWP" that was the predecessor of current Bioshaker®.

USER'S VOICE

Really amazed an effective shaking culture of heavy yeast in such a small well was realized!




The verification for mounting position and the variation among the wells and the effectiveness of shaking culture of yeast with 96-DWP.

Validation for NewMax drive. To reduce the variations due to the position of vessels and wells.


We verified the variation in turbidity (OD₆₀₀) for the yeast in four wells out of 96-deep well plate (squared well). Since yeast is a larger cell than E. coli, it is inclined to sink to the bottom of well in particular shaking culture with DWP if the stirring effectiveness is poor, which causes the variation in turbidity. At the same time, we also compared the popular shaking culture method using a Erlenmeyer flask to 96-deep well plate in turbidity (OD₆₀₀). The variation in among the position of deep wells are indicated with the bar graph. The deviation of each well plate in the same DWP is indicated with error bars in black (Fig.1). As a result, each one was no big difference in deviation. We have concluded that the stirring effectiveness at any positions in DWP made uniformly. Compared it with the Erlenmeyer flask in turbidity, finally the value doubled (OD₆₀₀ = above 8) of the upper value in the Erlenmeyer flask was obtained (Figure.2).

[Measurement position of OD₆₀₀]




The six wells per group was respectively allocated at the four corner and center of DWP. Measured the five deep wells allocated among from respectively five groups at the four corner and center of DWP. *The deep well circled in red were measured for the variation in Front left; Also same with Left back, Right front and Right back.

1) Variations were measured 9 hours since the shaking culture started.




Comparison between DWP and Erlenmeyer flasks in OD 600.



Comparison with the conventional model: Shaking culture of E.coli with 96-well plate

We made shaking culture of E. coli using M-BR-034P and the conventional model with 96-well plate to compare. 4 hours after the started we measured the turbidity of theirs; As a result, we discovered that OD₆₀₀ of M-BR-034P reached "OD₆₀₀=1.0" in half the time early that of former model even the shaking speed slower.

Turbidity 4 hours after the start of shaking culture



[Conditions] *All models in common

- 96-deep well plate: Square hole (conical bottom), Capacity approx. 2mL per well
- Sealing Gas Permeable Adhesive Seals
- Cultivation volume: 1 mL (approx. 50% of the well volume is optimal for aeration). In the case 1 mL, the solution does not reach the sealing up to 1500r/min while that reaches it at 1600r/min.

Maximizer MBR-022UP

Optimum for such as 96-well plate, Microtube rack and Vial. Possible for Temp control below R.T. For Multi-specimen processing such as Plasmid purified that requires a high-speed micro vibrations.


High-end model (improved uniformity and stability) --> P.013 Without temp. control "Micro mixer E-36" --> P.061

Features


- Optimum for stirring of 96-well plate and/or Microtube
- Enables Shaking culture with 96-deep well plate
- Peltier cooling enables temp. control from 15°C

Applications

- Shaking culture with Deep-well plate and Microtube rack
- ELISA and Phage display method
- Enzymatic digestion, Biodegradability tests, Metabolome analysis etc.



Optional parts : Microtube racks/Vial rack



| Product Name / Model | Microtube rack | | Vial rack |
|----------------------|----------------|----------------|-----------|
| | EM-0524 | EM-1524 | |
| Capacity of Vessels | 0.5mLx24 | 1.5/2.0mLx24 | 4mLx35 |
| Dimensions | 126x86x28.5Hmm | 126x86x40.5Hmm | — |

The suitable Microtubes are Eppendorf's.



96-well plate



USER'S VOICE


This device made me start Cultivation with Microtube rack and Deep well plate.

| Model | MBR-022UP |
|-----------------------------|---|
| Temperature range | 15°C to 60°C |
| Temp. control accuracy (*1) | ±0.3°C to ±0.8°C |
| Shaking motion | Horizontal eccentric shaking |
| Speed range | 300 to 2500r/min |
| Capacity | 2 x Well plate, 2 x Microtube rack (24 pcs) or 1 x Up to 4mL Vial rack (35 pcs) |
| Ambient temp. range | 5°C to 35°C |
| Heating/cooling method | Peltier element |
| Display (Temp./Speed) | Digital display |
| Program functions | 1 program, 9 segments (Temperature and Time) |
| Safety devices/functions | Fuse, Self-diagnosis function for temp. control error, Braking when hood open, Motor overload protection, Overload/Overvoltage protection |
| Dimensions (W x D x H) | 270 x 377 x 245mm |
| Weight | Appox. 12kg |
| Power supply | AC100 to 240V/1.5A |
| Standard accessories | 1x Plate adapter, 1 x Water-reception tray |

(*1)The value at the ambient temperature 25°C as the standard.

•Since designed for 96-well plate, the agitating effectiveness will be different by using different well sizes such as 24-well and 384-well plate.

Dimensions



Bioshaker MBR-024

Optimum for such as 96-well plate and Microtube. The capacity doubled as MBR-022UP. For Multi-specimen processing such as Plasmid purified.

High-end model (improved uniformity and stability) --> P.013 Invitro shaker Mix-EVR --> P.056



96-deep well plate

| Model | MBR-024 |
|-----------------------------|--|
| Temperature range (*1) | 5°C above RT (*2) to 65°C |
| Temp. control accuracy (*2) | ±0.3°C to ±1.0°C |
| Shaking motion | Horizontal eccentric shaking |
| Speed range | 100 to 1200r/min |
| Capacity | 4 x Well plate or 2 x Microtube rack (36 pcs) |
| Ambient temp. range | 15°C to 35°C |
| Heating method | Heater 80W |
| Display (Temp./Speed) | Digital display (Temp./Speed switchable) |
| Timer | 60min for operation |
| Safety devices/functions | Alarm of abnormal sample temp. (Set as desired), Braking when hood open, Motor overload protection, Self-diagnosis function for sensor error |
| Dimensions (W x D x H) | 300 x 377 x 318mm |
| Weight | Apporx. 12kg |
| Power supply | AC100V/2A |
| Standard accessories | 1 x 2 types adapter for plate |

(*1)It may not reach the maximum temperature by the surroundings.

(*2)The value at the ambient temperature 25°C.

•Since designed for 96-well plate, the agitating effectiveness will be different by using different well sizes such as 24-well and 384-well plate.

Features

- Optimum for stirring of 96-well plate or Microtube rack.
- Possible for shaking culture with 96-deep well plate.
- Microtube racks available as an option.

Applications

- Shaking culture with Well plate or Microtube rack
- ELISA and Phage display method
- Enzymatic digestion, Biodegradability tests etc.




96-deep well plate



1.5mm Microtube rack (option)

Optional parts : Microtube racks




EM-1036

| Model | EM-0836 | EM-1036 |
|---------------------|--------------|----------|
| Capacity of Vessels | 0.5mLx36 | 1.5mLx36 |
| Dimensions | 262x82x32Hmm | |

•Two tube racks must to hold at a time on the shaking rack.

The shape of Deep well plate (compatible with MBR series)

Common shape of Deep well plate is shown below left. Another shape Deep well plate with rims is shown below right and it requires adapter (standard or option depending models) for use.




Common shape DWP



DWP with rims

Dimensions



Bioshaker MBR-420FL

Large capacity model is optimum for 96-well plate.

Possible for Temp control from 4°C.

For Multi-specimen processing such as Plasmid purified.

Features

- Possible for shaking culture with 96-deep well plate
- Inverting function of shaking (Up to 1000 r/min)
- Possible for Temp. control from 4°C

Applications

- Shaking culture with Deep well plates
- ELISA and Phage display method
- Enzymatic digestion, Biodegradability tests etc.



Shaking platform for Deep well plates (Option)

Optional parts: Shaking platform

Shaking platform for Deep well plates →

6 plates can be placed on the platform. 2 sets of platform can be attached to the unit.



Shaking platform for Microplates →

6 plates can be placed on the platform. 4 sets of platform can be attached to the unit.



Example of use

| Model | MBR-420FL |
|-------------------------------|---|
| Temperature range (*1) | 4°C to 70°C |
| Temp. control accuracy | ±0.3°C to ±1.0°C |
| Shaking motion | Horizontal eccentric shaking (Inverting function) |
| Speed range | 200 to 1500r/min (Up to 1000 r/min at inverting shaking) |
| Capacity (*2) | 12 x Well plate or 24 x Microplate (*3) |
| Ambient temp. range | 5°C to 35°C |
| Shaking inverting function | The time and speed for inverting can be set up. |
| Heating/Cooling method | Heater 800W/Compressor 140W |
| Display (Temp./Speed) | Digital display |
| Program functions | 4 programs, 9 segments: Temp., Shaking speed and time/per segment (Setting range: 00h00min to 99h59min), Wait function (*4) |
| Safety devices/functions | Leakage/Overcurrent breaker, Fuse, Alarm of abnormal sample temp. (Set as desired), Detecting the door opening and closing, Self-diagnosis function for sensor error, Motor current limiter, Compressor overheating cut-off circuit |
| Dimensions (W x D x H)/Weight | 600 x 732 x 643mm, Approx. 110kg (Platform not included) |
| Power supply | AC100V/12A |

(*1)The temp. of defrost function that just starting is not included.


(*2)Shaking platform (Option) is required to fix the plates.

(*3)Varies depending the number of shaking platform (option).

(*4)Shaking will start when reaching the preset temp.

| Product Name / Model | Capacity of Vessels / Remarks | Mounting qty on MBR-420FL | Capacity of vessel per platform | |
|---|--|--------------------------------|---------------------------------|----|
| Shaking platform for Deep well plates DWP-2412 | Standard deep well plate | 2 | 6 | 12 |
| Adaptor ADP-2412 | Required for mounting Deep Well plate with rims to DWP-2412, 2 sets. | 2 sets of adaptor per platform | — | — |
| Shaking platform for Microplates MTP-2412 | Standard microplate (Do not use lid, Sealing films recommended) | 4 | 6 | 24 |
| Stick sheet shaking platform SR-2412 | Flat bottoms such as Petri dish can be held simply by placing on the platform. Max shakin speed 1000r/min. | 4 | — | — |

Dimensions



Bioshaker BR-21/22/23 Series

*Widely used for Cultivation of microbes etc.
that require temp. control. The compact sized can be
put under and on the Lab bench.*

Dimensions --> P.038 OD-Monitor Turbidimeter while shaking --> P.034
Shaker NR-1/2/3 --> P.062



[FP type]



BR-23FP-MR
The compact sized can be put
under and on the Lab bench.



BR-23UM-MR
Lift-up door allows good visibility
inside the chamber.

USER'S VOICE

FR type is CFC-free that relieves
us from the gas leakage and the
procedures for disposal.



Features

- "9 models" can be selected by temp. range, shaking motion and door (right or left-opening)
- Stackable up to Two levels. The platform can be changed.
- Capacity : Erlenmeyer flask 500mL x 4 pcs, 2L x 1 pc

Applications


- Cultivation of Microbes such as E.Coli.
- Cultivation of Yeast and Insect cell (FR)
- Cultivation of Thermophile and Hybridization (FH).

Compact design and Stackable

The compact sized can be placed under the Lab bench. Stackable up to Two levels with "Stackable base (Option)." Excellent performance for microbial cultures thanks to the shaking speed Max. 300r/min even 200r/min at the upper of stacked units.

Sliding retractable control panel

Sliding retractable control panel allows space saving and excellent in operability. Setting for Temperature, Shaking speed and Operation program with simply.



Drawer type control panel




| Temperature range | FP type: 15°C to 60°C (*1) | | | UM type: 5°C above RT to 70°C | | | FH type: 5°C above RT to 100°C | | |
|-----------------------------|--|--------------|------------------------|-------------------------------|----------------|------------------------|-----------------------------------|--------------|------------------------|
| Model | BR-21FP-MR | BR-22FP-MR | BR-23FP-MR | BR-21UM-MR | BR-22UM-MR | BR-23UM-MR | BR-21FH-MR | BR-22FH-MR | BR-23FH-MR |
| Shaking motion | Reciprocal | Orbital | Reciprocal/ Orbital | Reciprocal | Orbital | Reciprocal/ Orbital | Reciprocal | Orbital | Reciprocal/ Orbital |
| Temp. control accuracy (*1) | ±0.2°C | | | ±0.2°C to ±0.4°C | | | ±0.1°C to ±0.5°C | | |
| Speed range/Shaking width | 20 to 300r/min (200r/min at the upper of stacked units), 25mm | | | | | | | | |
| Platform size/Capacity | 290 x 250mm (Max. inside height: FP/FH = 270mm, UM = 315mm), Erlenmeyer flask 500mL x 4, 2L x 1 | | | | | | | | |
| Door type | Single swing door (Leftward open) | | | Swing lift-up door | | | Single swing door (Leftward open) | | |
| Ambient temp. range | 5°C to 35°C | | | | | | | | |
| Maximum load | Approx. 3.5kg (Platform and Vessel holder included) | | | | | | | | |
| Heateing method | Peltier element (4 x 70W) | | | Heater (350W) | | | | | |
| Cooling method | | | | - | | | | | |
| Display (Temp./Speed) | Digital display | | | | | | | | |
| Program functions | 1 program, 9 segments: Temp., Shaking speed and time/per segment (Setting range: 00h00min to 99h59min), Temp. gradient function (*2) | | | | | | | | |
| Safety devices/functions | Leakage/Overcurrent breaker, Fuse, Alarm of abnormal sample temp. (Set as desired), Detecting the door opening and closing, Self-diagnosis function for sensor error, Overcurrent protection of Peltier (FR) , Motor current limiter | | | | | | | | |
| Dimensions (W x D x H) | 437 x 585 x 545mm | | | 422 x 550 x 540mm | | | 437 x 550 x 545mm | | |
| Weight | Approx. 48kg | Approx. 49kg | Approx. 49.5kg | Approx. 41.5kg | Approx. 42.5kg | Approx. 43kg | Approx. 42kg | Approx. 43kg | Approx. 43.5kg |
| Power supply (*3) | AC100V/5A | | | | | | | | |
| Standard accessories | 1 x Spring net shaking platform (MR-2925), 1 x Drain nozzle, 1 x Drain hose | | | | | | | | |

(*1)The value at the ambient temperature 25°C. (*2)The transfer time of temperature can be set up. (*3)Voltage variation rate is ±5% for AC100V.

Optional accessories and Capacity/Qty of Vessels

Spring net shaking platform and Capacity/Qty of Vessels


| Model | MR-2925 | | | |
|---|--|---------|----------------------------|---|
|  | Vessels can be easily mounted on regardless its shapes. It can be used by removing some springs and changing the height of the spring net as necessary. Universal shaking platform (w/o Clamp) is option recommended if vessel should be tightly held. | | | |
| Spring net Pitch | 20mm | | | |
| Spring net stage | 1 stage (Height is changeable.) | | | |
| External dimension | 290x250x108Hmm | | | |
| Weight | Apporx.1.1kg | | | |
| Capacity of Vessel | φ16 mm test tubes (*1) | 50 (*2) | Erlenmeyer flask 300mL | 5 |
| | Disposable Centrifuge tube for 50 mL | 20 (*2) | Erlenmeyer flask 500mL | 4 |
| | Erlenmeyer flask 50mL | 20 | Erlenmeyer flask 1L | 2 |
| | Erlenmeyer flask 100mL | 10 | Erlenmeyer flask 2L | 1 |
| | Erlenmeyer flask 200mL | 9 | Sakaguchi flask 500mL (*3) | 2 |
| | Erlenmeyer flask 250mL | 5 | | |

(*1) φ16 mm or less should be with the stand to be mounted on in Spring net.

(*2) The qty when Angle rack tilted 45°.

(*3) Be careful when using Shake flask as it is fragile.

Universal shaking platform and Capacity/Qty of Flasks



| Model | MT-2925 | | | |
|---|--|----|------------------------|---|
|  | It is used together Clamps, Angle racks and "OD-Monitor" series (Clamps of Erlenmeyer flask and Shake flask are available). It can hold flasks tightly than Spring shaking platform. | | | |
| External dimension | 290x250x23Hmm | | | |
| Weight | Apporx.0.9kg | | | |
| Capacity of Flask | Erlenmeyer flask 50mL | 20 | Erlenmeyer flask 500mL | 4 |
| | Erlenmeyer flask 100mL | 10 | Erlenmeyer flask 1L | 2 |
| | Erlenmeyer flask 200mL | 9 | Erlenmeyer flask 2L | 1 |
| | Erlenmeyer flask 250mL | 5 | Sakaguchi flask 500mL | 3 |
| | Erlenmeyer flask 300mL | 5 | | |
| | | | | |

Angle rack for Test tubes and Disposable Centrifuge tubes (Mounted on MT-2925 for use)


| Model | AT-1250 | AT-1620 | AT-1650 | AT-3518 (*1) |
|-------------------------|----------------------|------------------------|------------------------|--|
| Vessels / Capacity | φ12mm test tubes x50 | φ16.5mm test tubes x20 | φ16.5mm test tubes x50 | Disposable Centrifuge tube for 50 mL x18 |
| Mounting qty on MT-2925 | 1 (*2) | 2 (*2) | 1 (*2) | 1 |
| External dimension (mm) | 250x100x95H | 285x70x80H | 286x110x100H | 285x110x102H |
| Weight | Apporx.0.8kg | Apporx.0.6kg | Apporx.0.9kg | Apporx.0.8kg |

(*1) Be sure to be used below 50°C. (*2) The qty when Angle rack with 160mm test tube tilted 45°.

Other optional accessories

| Product Name / Model | Stackable base STB-4842 | Inner rack RSB-3430 |
|----------------------|--|---|
| Shape |  |  |
| Description | The fixing tool that prevents slippage when the stacked in use (Required when stacked) | It can be installed to the upper of chamber. The height of space that the sample placed differs depending on the mounting position. |

Stick sheet shaking platform

| Model | SR-2925 |
|--|--|
|  | Flat bottoms such as Erlenmeyer flask and Petri dish can be held simply by placing on the platform. The decreasing an adhesive force by dirt can be recovered by washing with water and drying it. Even Erlenmeyer flasks can be used up to 200r/min by using Optional Clamp belt. |
| Max. shaking speed | 100r/min (200r/min when using Clamp belt) |
| External dimension / Weight | 290x250x30Hmm / Apporx.1.1kg |
| Standard accessories | Stick sheetx1 |

Stick sheet, Clamp belt



Stick sheet (ST-2925) above also available as sold separately.

Clamp belt (CB-0205/1000) is for securely holding Erlenmeyer flask on Stick sheet shaking platform.

| Model | ST-2925 | Model | CB-0205 | CB-1000 |
|--------------------|-----------|------------------------------|---------------|--------------|
| External dimension | 290x250mm | Capacity of Erlenmeyer flask | 200 to 500 mL | 1 to 3 L (*) |

(*) Up to 2L of Erlenmeyer flask can be used with the BR-21/22/23 series.


Clamps (Mounted on MT-2925 and AF-2070 for use.)

| Shape | Vessels / Capacity | Model |
|--|--------------------|---------------|
|  Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet. | 50mL | CF-0050 |
| | 100mL | CF-0100 |
| | 200mL | CF-0200 |
| | 250mL | CF-0250 |
| | 300mL | CF-0300 |
| | 500mL | CF-0500 |
| | 1L | CF-1000 |
| | 2L | CF-2000 |
|  MT-2925 with Multiple clamps mounted. | Sakaguchi flask | 500mL SF-0500 |

Angle rack for Erlenmeyer flasks (Mounted on MT-2925 *Clamps required for use)

| Model | AF-2070 |
|-----------------------------|---|
| Capacity | 50mL:4, 100mL:3, 200 to 300mL:2, 500mL:1 |
| Mounting qty on MT-2925 | 2 sets (The qty when Angle rack tilted 30°) |
| External dimension / Weight | 250x110x50Hmm, Apporx.0.8kg |

Tube holders (Mounted on MT-2925 for use)

| | |
|---|--|
|  M-0150 M-0500 | Tube holder M-0015 3 sets x 8 racks for 1.5/2.0 mL Microtube and Micro Vial |
| Tube holder M-0150 Disposable Centrifuge tube for 15 mL 3 sets x 6 racks (*) | Tube holder M-0500 Disposable Centrifuge tube for 50 mL 3 sets x 4 racks (*) |

(*) Up to 2 sets can be mounted on MT-2925 in BR-21/22/23 series.

The mounting direction of Angle rack. AT-3518

The tilt direction of Angle rack should be mounted orthogonal to the shaking direction when reciprocating due to the inside dimension of chamber.



Bioshaker BR-53FP

Peltier cooling with CFC-free. Temperature range ensures 15°C to 55°C. Cut power consumption by 80% with the conventional one.



USER'S VOICE

CFC-free relieves us from the gas leakage and the procedures for disposal. We should apply for "Subsidy scheme for energy efficiency".



Features

- Peltier cooling with CFC-free. Battery friendly.
- Capacity : Erlenmeyer flask 500mL x 8 pcs, 5L x 1 pc
- Stackable to two levels. Wide variation platforms available as an option.

Applications

- Cultivation of microbes such as Yeast, Insect cell, E.Coli. etc.
- Cultivation of Protein expression vector for E.coli.
- Food testing system and inspection for QC.

Peltier cooling with CFC-free.

Cut power consumption by 80% compared with the conventional one (100r/min at 37°C). Energy saving with 7A. CFC-free relieves you from the gas leakage and the procedures for disposal.

| Model | Size | Power supply | Power consumption |
|----------------|---------------|--------------|----------------------------------|
| BR-23FP | Small | 5A | 80Wh (at +37°C, 100r/min) |
| BR-53FP | Medium | 7A | 90Wh (at +37°C, 100r/min) |
| BR-43FL | Medium | 12A | 510Wh (at +37°C, 100r/min) |

Open rightward and leftward "Double hinged door"

Open rightward and leftward "Double hinged door" without caring about a place of installation.

Detachable shaking platform. Wide variation available.

Comes with Spring net shaking platform. Others such as Universal shaking platforms available as an option.


Stackable to two levels

Different temperature can be set respectively upper unit and lower unit for instance, which contribute a space saving.

The front panel in black has improved the high-visibility on the display.

The temp. range 5°C to 55°C works well for general cultivation. The program function enables to correspond to the various protocols such as the temp. range can be changed to another one after shaking. Temp. control is stable thanks to Peltier cooling.

Enables the culture at normal temp.



The left figure shows a test result for cultivation of E. coli measured by Non-invasive turbidity meter. Cell growth reached OD value 2.5 (Upper limit) in 55 hours.

| Model | BR-53FP |
|-------------------------------|---|
| Temperature range (*1) | 15°C to 55°C |
| Shaking motion | Switchable Reciprocal/Orbital shaking |
| Temp. control accuracy | ±0.3°C |
| Speed range/Shaking width | 20 to 300r/min (200r/min at the upper of stacked units) / 25mm |
| Capacity | Universal shaking platform MT-4430: Erlenmeyer flask 500mLx8 pcs Universal shaking platform MT-4030: Erlenmeyer flask 500mLx6 or 5L x 1 |
| Door type | Double hinged door (Rightward/Leftward open) |
| Maximum load | Approx. 7 kg (Platform included) |
| Ambient temp. range (*2) | 5°C to 30°C |
| Heating method | Peltier element |
| Program functions | 1 program, 9 segments: Temp., Shaking speed and time/per segment (Setting range: 00h00min to 99h59min), Temp. gradient function (*3) |
| Safety devices/functions | Leakage/Overcurrent breaker, Fuse, Alarm of abnormal sample temp. (Set as desired), Detecting the door opening and closing, Self-diagnosis function for sensor error, Motor current limiter |
| Dimensions (W x D x H)/Weight | 600 x 686 x 643mm, Approx. 86 kg (Platform included) |
| Power supply | AC100V to 240V/7A |
| Standard accessories | 1 x Power cable, 1 x Spring net shaking platform (MR-4030) |

(*1) The value at 25°C of the rear surface intake temp. (*2) Dew condensation (*3)The transfer time of temperature can be set up.

External output compatible with Pharmaceutical validation (Customization).

Customized external output for Pharmaceutical validation available.

- Temperature output •Speed output
- Collectively signal output when safety device activated.

Optional parts: Inner rack

RSB-5032

It can be installed to the upper of chamber. The Maximum load 1kg.




Optional parts: Stackable base

STB-6070


The fixing tool for Two levels.

Shaking platform and Capacity/Qty of Flasks

| Model | MT-4430 | |
|--|--|----|
|  | <p>The universal platform for 250mL to 1L. The capacity (Num. of flasks) is larger than that of shaking platform MT-4030. It can hold flasks tightly than shaking platform MR-4030</p> | |
| External dimension / Weight | 440x300x23Hmn, Apporx.1.4kg | |
| Capacity of Flask | Erlenmeyer flask 250mL | 12 |
| | Erlenmeyer flask 300mL | 12 |
| | Erlenmeyer flask 500mL | 8 |
| | Erlenmeyer flask 1L | 6 |

Connection holes included as standard equipment.

Equipped with Connection holes for installing such as the sensor of recorder inside the chamber. It is used with Junction pipe (Option).



Optional parts: Junction pipe

JCP-3657



The pipe (w/o Lid) for connection to the outside. Be sure be used below 50°C. (*)

(*)Electrical equipment should not be placed inside the chamber.



Optional parts: Clamps

Mounted on Platforms MT-4030/4430 and Angle racks AF-2070/2570 for use.

| Vessels / Capacity | | Model |
|--|---|----------------|
|  Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet. | 50mL | CF-0050 |
| | 100mL | CF-0100 |
| | 200mL | CF-0200 |
| | 250mL | CF-0250 |
| | 300mL | CF-0300 |
| | 500mL | CF-0500 |
| | 1L | CF-1000 |
| | 2L | CF-2000 |
| | 3L | CF-3000 |
| | 5L | CF-5000 |
|  MT-4030 with Multiple clamps mounted. | Up to 2L sets of Erlenmeyer flask can be used in BR-40LF. | |
| | Sakaguchi flask | |
| | Medium bottle | |
| | 500mL | SF-0500 |
| | 1L | SF-1000 |
| | 2L | SF-2000 |
| | 250mL | MF-0250 |

Other optional accessories are compatible with BR-4X series. -->P.023

Example of customization for BP-53FR Customized model "Equipped with HEPA (High efficiency particulate air) filter for CPC (Cell Processing Center)".



Stackable up to Two levels

Supply record of EU-GMP (European Union Good Manufacturing Practice) Grade B cleanrooms.

HEPA filter removes the particle in exhaust wind. Supplied to EU-GMP (European Union Good Manufacturing Practice) Grade B clean rooms that requires Strict cleanliness.

| Model | BR-53FP-CM |
|-------------------------------|--|
| Temp. control accuracy (*1) | 18°C to 50°C |
| Shaking motion | Switchable Reciprocal/Orbital shaking |
| Temp. control accuracy | ±0.5°C |
| Speed range/Shaking width | 20 to 300r/min, 25mm |
| HEPA Filter | Accumulation rate: Above 99.99% for above 0.3µm particle |
| Door type | Double hinged door (Rightward/Leftward open) |
| Maximum load | Approx. 7 kg (Platform included) |
| Ambient temp. range (*2) | 18°C to 25°C |
| Dimensions (W x D x H)/Weight | 625 x 890 x 650mm, Approx. 108 kg |
| Power supply | AC100V or 220 to 240V/7A |

(*1) The value at 25°C of the rear surface intake temp. (*2) Dew condensation.

Bioshaker BR-40/41/42/43 Series

Widely used for Cultivation of microbes such as Yeast, E.Coli. etc. that require temp. control. Sufficient capacity, Space saving and High versatility.

Dimensions --> P.038 OD-Monitor Turbidimeter while shaking --> P.034 LED irradiation unit --> P.037 Shaker NR-10/20/30 --> P.064

Lift-up door offers allows good visibility inside chamber (BR-40LF)




Open rightward and leftward "Double hinged door" (BR-41/42/43 series)

Features

- "8 models" can be selected by temp. range and shaking motion.
- Capacity : Erlenmeyer flask 500mL x 6 pcs, 5L x 1 pc
- Stackable up to Two levels. Platform can be changed.

Applications

- Cultivation of microbes such as E. Coli. etc.
- Cultivation of Yeast and Insect cell (LF/FL)
- Cultivation of Thermophile and Hybridization (FH).



Open rightward and leftward "Double hinged door" and "Stackable"

Open rightward and leftward "Double hinged door". Stackable up to Two levels with "Stackable base (Option)." Excellent performance for Cultivation of microbe thanks to the shaking speed Max. 300r/min even 200r/min at the upper of stacked units.

LED irradiation unit is available as an option.

LED (Light emitting diode) irradiation unit available as an option. Good effect on photosynthesis. * Recommended to use with Small and Medium sized unit (See ? page for details)



| Temperature range | 4°C to 70°C | | | | 5°C above RT to 70°C | | | 5°C above RT to 100°C |
|-----------------------------|--|---|-------------------|------------------------|----------------------|------------|------------------------|-----------------------|
| Model | BR-40LF | BR-41FL-MR | BR-42FL-MR | BR-43FL-MR | BR-41FM-MR | BR-42FM-MR | BR-43FM-MR | BR-43FH-MR |
| Shaking motion | Reciprocal/ Orbital | Reciprocal | Orbital | Reciprocal/ Orbital | Reciprocal | Orbital | Reciprocal/ Orbital | Reciprocal/Orbital |
| Temp. control accuracy (*1) | ±0.3 to ±1.0°C | | | | | | | |
| Speed range | 20 to 200r/min | 20 to 300r/min (200r/min at the upper of stacked units) | | | | | | |
| Shaking width | 10 to 40mm | 25mm | | | | | | |
| Capacity | Universal shaking platform MT-4430: Erlenmeyer flask 500mLx8, Universal shaking platform MT-4030: Erlenmeyer flask 5Lx1 *BR-40LF: 2L x 2 (*2) | | | | | | | |
| Door type | Swing lift-up door | Double hinged door (Rightward/Leftward open) | | | | | | |
| Maximum load | approx. 7 kg (Platform and Vessel holder included) | | | | | | | |
| Ambient temp. range | 5°C to 30°C | 5°C to 35°C | | | | | | |
| Heating method | Heater 500W | Heater 800W | | | | | | |
| Cooling method | Compressor 75W | Compressor 140W | | | - | | | |
| Timer/Program functions | BR-40LF (Timer): Temp. transition time and ON-OFF changeover (Non-functional by only temp. transition time *e.g.: Shaking at 37°C → Shaking stop at 4°C.) BR-41/42/43 (Program): 1 program, 9 segments: Temp., Shaking speed and time/per segment (Setting range: 00h00min to 99h59min), Temp. gradient function (*3) | | | | | | | |
| Safety devices/functions | Leakage/Overcurrent breaker, Fuse, High temp, Alarm of abnormal sample temp. (Set as desired), Detecting the door opening and closing, Self-diagnosis function for sensor error, Motor current limiter, Compressor overheating cut-off circuit (BR-40LF and FL) | | | | | | | |
| Dimensions (W x D x H) | 586 x 630 x 679mm | | 600 x 686 x 643mm | | | | | |
| Weight | Approx. 70 kg | Approx.95 kg | | | Approx.77 kg | | Approx.80 kg | |
| Power supply | AC100V/9A | | AC100V/12A | | | AC100V/9A | | |
| Standard accessories | 1 x Spring net shaking platform (MR-4030), 1 x Drain-reception plastic bottle with hose (BR-40LF) | | | | | | | |

(*1)The temp. of defrost function that just starting is not included.


(*2)Maximum inside height are 331mm (40LF) and 350mm (41/42/43).

(*3)The transfer time of temperature can be set up.

•Temp. range and Temp. stability may not be meet the specs when frost by continuous operation. Required to defrost it regularly.

Optional accessories and Capacity/Qty of Vessels

Spring net shaking platform and Capacity/Qty of Vessels


| Model | MR-4030 | | | |
|---|---|---------|--------------------------|---|
|  | Vessels can be easily mounted on regardless its shapes. It can be used by removing some springs and changing the height of the spring net as necessary. Universal shaking platform recommended for 5 L Erlenmeyer flask or 1 to 2L Shake flask. | | | |
| Spring net Pitch | 20mm | | | |
| Spring net stage | 2 stage (Height is changeable.) | | | |
| External dimension / Weight | 400x300x108Hmm, Apporx.3.3kg | | | |
| Capacity of Vessel | φ16 mm test tubes (*1) | 96 (*2) | Erlenmeyer flask 300mL | 9 |
| | Disposable Centrifuge tube for 50 mL | 48 (*2) | Erlenmeyer flask 500mL | 6 |
| | Erlenmeyer flask 50mL | 24 | Erlenmeyer flask 1L | 3 |
| | Erlenmeyer flask 100mL | 20 | Erlenmeyer flask 2L | 2 |
| | Erlenmeyer flask 200mL | 12 | Erlenmeyer flask 3L (*3) | 1 |
| | Erlenmeyer flask 250mL | 9 | Sakaguchi flask 500mL | 6 |

(*1) φ16 mm or less should be with the stand to be mounted on in Spring net.

(*2) The qty when Angle rack tilted 45°.

(*3) It cannot be used with BR-40LF.

Shaking platforms and Capacity of Flasks

| Model | MT-4030/MT-4430 | | | |
|--|---|--------|-------------------------------|-------|
|  | MT-4030 is used together Clamps, Angle racks and "OD-Monitor" series. It can hold flasks tightly than Spring shaking platform. In MT-4430 the mounting position for Flasks are fixed and the capacity of (Num. of flasks) is larger than that of MT-4030. | | | |
| External dimension / Weight | 400/440x300x23Hmm, Apporx.1.4kg | | | |
| Capacity of Flask | Erlenmeyer flask 50mL | 35 / - | Erlenmeyer flask 2L | 2 / - |
| | Erlenmeyer flask 100mL | 18 / - | Erlenmeyer flask 3L (*1) | 2 / - |
| | Erlenmeyer flask 200mL | 12 / - | Erlenmeyer flask 5L (*1) (*2) | 1 / - |
| | Erlenmeyer flask 250mL | 9/12 | Sakaguchi flask 500mL | 6 / - |
| | Erlenmeyer flask 300mL | 9/12 | Sakaguchi flask 1L | 4 / - |
| | Erlenmeyer flask 500mL | 6 / 8 | Sakaguchi flask 2L | 2 / - |
| | Erlenmeyer flask 1L | 4 / 6 | | |

(*1) It cannot be used with BR-40LF. (*2) The height is up to 350mm height.

Angle rack for Test tubess (Mounted on MT-4030 for use)

| Model | AT-1250 | AT-1620 | AT-1650 | AT-1820 |
|-------------------------|--------------|--------------|--------------|--------------|
| Capacity | φ12mm x50 | φ16.5mm x20 | φ16.5mm x50 | φ18mm x20 |
| Mounting qty on MT-4030 | 2 (*1) | 3 (*1) | 2 (*1) | 3 (*1) |
| External dimension (mm) | 250x100x95H | 285x70x80H | 286x110x100H | 320x70x80H |
| Weight | Apporx.0.8kg | Apporx.0.6kg | Apporx.0.9kg | Apporx.0.7kg |

(*1) The qty when Angle rack with 160mm test tube tilted 45°.

Other optional accessories

Inner rack RSB-5032

It can be installed to the upper of chamber. The Maximum load 1kg.

| Stackable base | Junction pipe |
|---|---|
| The fixing tool for Two stages. For BR-41/42/43/53: STB-6070 For BR-40LF: Stackable baseA | The pipe (w/o Lid) for connection to the outside. (JCP-3657) Be sure to be used below 50°C. (*1) |


(*1) Electrical equipment should not be placed inside the chamber.



Angle rack for Disposable Centrifuge tubes

| Model | AT-3518 |
|-------------------------|--------------------------------|
| Capacity | 50mLx18 |
| Mounting qty on MT-4030 | 2 (When Angle rack tilted 45°) |
| External dimension | 285x110x102Hmm |
| Weight | Apporx.0.8kg |

Stick sheet shaking platform

| Model | SR-4030 |
|--|--|
|  | Flat bottoms such as Erlenmeyer flask and Petri dish can be held simply by placing on the platform. The decreasing an adhesive force by dirt can be recovered by washing with water and drying it. Even Erlenmeyer flasks can be used up to 200r/min by using Optional Clamp belt. |
| Max. shaking speed | 100r/min (200r/min when using Clamp belt) |
| External dimension / Weight | 400x300x30Hmm, Apporx.1.7kg |
| Standard accessories | Stick sheetx1 |

Stick sheet, Clamp belt


Stick sheet (ST-4030) above also available as sold separately. Clamp belt (CB-0205/1000) is for securely holding Erlenmeyer flask on Stick sheet shaking platform.

| Model | ST-4030 | Model | CB-0205 | CB-1000 |
|--------------------|-----------|------------------------------|---------------|---------------|
| External dimension | 400x300mm | Capacity of Erlenmeyer flask | 200 to 500 mL | 1 to 3 L (*1) |


(*1) Up to 2L of Erlenmeyer flask can be used with the BR-40LF.

Clamps (Mounted on MT-4030 and AF-2070/2570 for use)

| Vessels / Capacity | | Model |
|---|-------|---------|
| Up to 2L sets of Erlenmeyer flask can be used in BR-40LF. | 50mL | CF-0050 |
| | 100mL | CF-0100 |
| | 200mL | CF-0200 |
| | 250mL | CF-0250 |
| | 300mL | CF-0300 |
| | 500mL | CF-0500 |
| | 1L | CF-1000 |
| | 2L | CF-2000 |
| | 3L | CF-3000 |
| | 5L | CF-5000 |
| Sakaguchi flask | 500mL | SF-0500 |
| | 1L | SF-1000 |
| | 2L | SF-2000 |



Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet.



MT-4030 with Multiple clamps mounted.



MT-4030 with Multiple clamps mounted.

Angle rack for Erlenmeyer flasks (Mounted on MT-4030 *Clamps required)

| Model | AF-2070 | AF-2570 |
|-------------------------|--|--|
| Capacity | 50mL:4, 100mL:3, 200 to 300mL:2, 500mL:1 | 50mL:5, 100mL:4, 200mL:3, 250 to 500mL:2 |
| Mounting qty on MT-4030 | 2 (*1) | 2 (*1) |
| External dimension | 250x110x50Hmm | 300x110x50Hmm |
| Weight | Apporx.0.8kg | Apporx.0.9kg |

(*1) The qty when Angle rack tilted 30°. *Clamps are available as an option.

Tube holders (Mounted on MT-4030 for use)

| Tube holder M-0015 | | Tube holder M-0500 | |
|--|--|--|--|
|  | | 3 sets x 8 racks for 1.5/2.0 mL Microtube and Micro Vial | |
| Tube holder M-0150 | | Tube holder M-0500 | |
|  | | Disposable Centrifuge tube for 50 mL 3 sets x 4 racks | |



Disposable Centrifuge tube for 15 mL 3 sets x 6 racks

Disposable Centrifuge tube for 50 mL 3 sets x 4 racks

Bioshaker BR-300LF/BR-300LF LED irradiation

Widely used for Cultivation of microbes such as Yeast, E.Coli. etc. that require temp. control. The power saving LED irradiation model and a model without compressor available.

Dimensions --> P.038 OD-Monitor Turbidimeter while shaking --> P.034 Double Shaker NR-150 --> P.066





BR-300LF

BR-300LF LED irradiation

USER'S VOICE

Newly adopted LED irradiation unit for Power saving.




Providing a line up of incorporated LED irradiation unit.

Large Bioshakers are long seller by continued improvements such as Shaking platform drawer lever. Equipped with Switchable shaking motion and adjustment of shaking width in stepless between 10mm and 50mm. LED irradiation unit has been adopted to save power more than conventional one. It has been proved comparable to Fluorescent lamp in the cultivation of algae.



Power code
(Lower left back)



7x Straight tube LED

Equipped with the remote terminal for the temp. control.

Connecting with Program unit PU-5 (see the right page) enables the remote control setting and the program operation.



| Model | BR-300LF | BR-300LF-LED |
|-----------------------------------|---|--|
| Temperature range | 4°C to 70°C | |
| Temp. control accuracy (*1) | ±0.3°C to ±1.0°C | |
| Shaking motion/Speed range/Width | Switchable Reciprocal/Orbital shaking, 25 to 250r/min, 10 to 50mm (Stepless) | |
| Platform size/Capacity | 600 x 400mm (Max. inside height: 575mm), Erlenmeyer flask 500mLx15, 5Lx3 | |
| Light source/Irradiation capacity | - | 7 x Straight tube LED, 8000 lux (*2), Stepless dimming |
| Irradiation time setting | - | 24-hour timer setting method (Switchable ON/OFF per 15min) |
| Door type | Single swing door (Leftward open) | |
| Maximum load | Approx. 10kg (Platform not included) | |
| Ambient temp. range (*1) | 5°C to 35°C | |
| Display (Temp./Speed) | Digital display (Switchable Set temp./Chamber temp.) | |
| Heating/Cooling method | Heater 700W/Compressor 125W | |
| Other functions | Shaking duration integrator (0.1h to 999.9h, Automatic resetting function *Required Program unit PU-5 or 6/Option to use), Chamber lamp | |
| Safety devices/functions | Leakage/Overcurrent breaker, Compressor overheating cut-off circuit, Abnormal sample temp. protection, High/Low temp., Display of deviation from set value (upper / lower limit), Self-diagnosis function for temperature sensor error, | |
| Dimensions (W x D x H)/Weight | 1110 x 682 x 1120mm, Approx. 210kg | 1110 x 682 x 1218mm, Approx. 230kg |
| Power supply | AC220 to 240V/11A | AC220 to 240V/12A |
| Standard accessories | 1 x Universal shaking platform MT-6040E (*3), 4 x Vibration absorbing pad, 4 x Spacer, 1 x Drain-reception plastic bottle (With hose), 2 x Safety pin for the mounting stay (BR-300LF-LED), 1 x Light shielding plate (BR-300LF-LED) | |

(*1)The temp. of defrost function that just starting is not included.


(*2)The value was measured at the position of 230mm from the center of platform. The value is 5500 lux at the position of 30mm from that.

(*3)Clamps not included.


•Temp. range and Temp. stability may not meet the specs when frost by continuous operation. Required to defrost it regularly.

Optional accessories and Capacity/Qty of Vessels

Shaking platform and Capacity/Qty of Flasks

| Model | | MT-6040E | | |
|---|------------------------|--|-------------------------------------|----|
|  | | It is used together Clamps, Angle racks and "OD-Monitor" series (Clamps of Erlenmeyer flask and Sakaguchi flask available). It can hold flasks tightly than Spring shaking platform. | | |
| External dimension / Weight | | 600x400x42Hmm, Apporx.4kg | | |
| Capacity of Flask | Erlenmeyer flask 50mL | 70 | Erlenmeyer flask 2L | 6 |
| | Erlenmeyer flask 100mL | 35 | Erlenmeyer flask 3L | 4 |
| | Erlenmeyer flask 200mL | 28 | Erlenmeyer flask 5L/Disposable 3&5L | 3 |
| | Erlenmeyer flask 250mL | 21 | Sakaguchi flask 500mL | 15 |
| | Erlenmeyer flask 300mL | 21 | Sakaguchi flask 1L | 9 |
| | Erlenmeyer flask 500mL | 15 | Sakaguchi flask 2L | 6 |
| | Erlenmeyer flask 1L | 9 | Sakaguchi flask 3L | 4 |

Spring net shaking platform and Capacity/Qty of Vessels (Mounted on MT-6040E for use)



| Model | MR-4030L | | | |
|--|---|--|--|--|
|  | | The usable Shaking platform that Vessels can be easily mounted on regardless its shapes. Fixed with the screws to Platform MT-6040E for using. Up to Two sets can be mounted. (See Left figure.) | | |
| Spring net Pitch | 20mm | | | |
| Spring net stage | 2 stage (Height is changeable.) | | | |
| External dimension / Weight | 400x300x108Hmm, Apporx.2kg | | | |
| Capacity of Vessels (≈) | φ16 mm test tubes (*1) | 96pcs (When Two sets mounted: 2 x 96 = 192pcs). | | |
| | Disposable Centrifuge tube for 50 mL (*1) | 48pcs (When Two sets mounted: 2 x 48 = 96pcs). | | |
| | Erlenmeyer flask | 50mL | 24pcs (When Two sets mounted: 2 x 24 = 48pcs). | |
| | | 100mL | 20pcs (When Two sets mounted: 2 x 20 = 40pcs). | |
| | | 200mL | 12pcs (When Two sets mounted: 2 x 12 = 24pcs). | |
| | | 250mL | 9pcs (When Two sets mounted: 2 x 9 = 18pcs). | |
| | | 300mL | 9pcs (When Two sets mounted: 2 x 9 = 18pcs). | |
| | | 500mL | 6pcs (When Two sets mounted: 2 x 6 = 12pcs). | |
| | | 1L | 3pcs (When Two sets mounted: 2 x 3 = 6pcs). | |
| | | 2L | 2pcs (When Two sets mounted: 2 x 2 = 4pcs). | |
| | | 3L | 1pcs (When Two sets mounted: 2 x 1 = 2pcs). | |
| | Sakaguchi flask 500mL | 6pcs (When Two sets mounted: 2 x 6 = 12pcs). | | |

(*1)The qty when Angle rack tilted 45°. φ16mm or less test tube should fixed with the stand per net.



(*2)The Capacity/Qty of vessels are less than those of Universal platform.

Program units

Connecting to the remote temperature setting terminal of BR-300LF series enables remote setting and program setting of temperature (The modification of BR is required for the shaking control).

| | | | |
|---|------------|---|--|
|  | Model | PU-5 | |
| | Functions | Remote program control of Temperature. | |
|  | Main Spec. | Storage programs: 2 program (8 segments per program) | |
| | Model | PU-6 | |
| | Functions | Remote program control of Temperature and Shaking ON/OFF. | |
| | Main Spec. | Storage programs: 2 program (8 segments per program) Time signal (ON/OFF): Each 2 programs | |

Clamps (Mounted on MT-6040E and AF-2070/2570 for use)

| | Vessels / Capacity | Model |
|---|--------------------|------------|
| | | |
|  Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet. | 50mL | CF-0050 |
| | 100mL | CF-0100 |
| | 200mL | CF-0200 |
| | 250mL | CF-0250 |
| | 300mL | CF-0300 |
| | 500mL | CF-0500 |
| | 1L | CF-1000 |
| | 2L | CF-2000 |
| | 3L | CF-3000 |
| | 3&5L (*) | CF-3000DSP |
|  Clamps for Sakaguchi flask and Sakaguchi flask | 500mL | SF-0500 |
| | 1L | SF-1000 |
| | 2L | SF-2000 |
| | 3L | SF-3000 |

(*)For Disposable Erlenmeyer flasks (Corning 3L and Thomson 5L).

Angle rack for Test tubess (Mounted on MT-6040E for use)

| Model | AT-1250 | AT-1620 | AT-1650 | AT-1820 |
|--------------------------|--------------|--------------|--------------|--------------|
| Capacity | φ12mm x50 | φ16.5mm x20 | φ16.5mm x50 | φ18mm x20 |
| Mounting qty on MT-6040E | 4 (*) | 6 (*) | 4 (*) | 5 (*) |
| External dimension (mm) | 250x100x95H | 285x70x80H | 286x110x100H | 320x70x80H |
| Weight | Apporx.0.8kg | Apporx.0.6kg | Apporx.0.9kg | Apporx.0.7kg |

(*)The qty when Angle rack with 160mm test tube tilted 45°.

Angle rack for Disposable Centrifuge tubess (Mounted on MT-6040E for use)

| Model | AT-3518 |
|-----------------------------|--------------------------------|
| Capacity | 50mLx18 |
| Mounting qty on MT-6040E | 4 (When Angle rack tilted 45°) |
| External dimension / Weight | 285x110x102Hmm, Apporx.0.8kg |
| Remarks | Be sure be used below 50°C. |

Angle rack for Erlenmeyer flasks (Mounted on MT-6040E *Clamps required)

| Model | AF-2070 | AF-2570 |
|--------------------------|--|--|
| Capacity | 50mL:4, 100mL:3, 200 to 300mL:2, 500mL:1 | 50mL:5, 100mL:4, 200mL:3, 250 to 500mL:2 |
| Mounting qty on MT-6040E | 4 (*) | 4 (*) |
| External dimension | 250x110x50Hmm | 300x110x50Hmm |
| Weight | Apporx.0.8kg | Apporx.0.9kg |

(*)The qty when Angle rack tilted 30°. *Clamps are available as an option.

Bioshaker Double platform BR-300LF

Enables to use widely for Cultivation of microbes such as Yeast, E.Coli. etc. that require temp. control. Equipped with double platform, the capacity is doubled as BR-300LF (Single platform).

Dimensions --> P.038 OD-Monitor Turbidimeter while shaking --> P.034 Double Shaker NR-150 --> P.066



Double platform BR-300LF



Double platform realizes the capacity twice as BR-300LF (Single platform).

Large Bioshakars are long seller by the repeated improvements such as Shaking platform drawer lever. The capacity (Up to 1 liter Erlenmeyer flask) doubled as BR-300LF (Single platform). Note: The upper platform should be removed when above 2L Erlenmeyer flask and 1L Shake flask are placed. (Ask us for details).



Shaking platform drawer lever

Features


- Wide temperature range for 4°C to 70°C
- Capacity : Erlenmeyer flask 500mL x 30 pcs, 1Lx18 pcs
- Shaking width adjustable in stepless between 10mm and 50mm.

Applications

- Cultivation of microbes such as E.Coli. etc.
- Cultivation of Yeast, Psychrophile and Insect cell
- The other samples that require shaking and temp. control.

Equipped with the remote terminal for the temp. control.

Conectting with Program unit PU-5 (see the right page) enables the remote control setting and the program operation.



| Model | Double platform BR-300LF |
|-------------------------------|---|
| Temperature range | 4°C to 70°C |
| Temp. control accuracy (*1) | ±0.3°C to ±1.0°C |
| Shaking motion | Switchable Reciprocal/Orbital shaking |
| Speed range | 25 to 160r/min |
| Shaking width | 10 to 50mm (stepless) |
| Platform size/Capacity | 600 x 400mm (Capable/Max. vessel size: Up to 1L Erlenmeyer flask and 500mL Shake flask), Erlenmeyer flask 500mLx30, 1Lx18 |
| Door type | Single swing door (Leftward open) |
| Maximum load | Approx. 10kg/per platform (Platform not included) |
| Ambient temp. range (*1) | 5°C to 35°C |
| Display (Temp./Speed) | Digital display (Switchable Set temp./Chamber temp.) |
| Heating/Cooling method | Heater 700W/Compressor 125W |
| Other functions | Shaking duration integrator (0.1h to 999.9h, Automatic resetting function *Required Program unit PU-5 or 6/Option to use), Chamber lamp |
| Safety devices/functions | Leakage/Overcurrent breaker, Compressor overheating cut-off circuit, Alarm of abnormal sample temp. (Set as desired), High/Low temp., Display of deviation from set value (upper / lower limit), Self-diagnosis function for temperature sensor error |
| Dimensions (W x D x H)/Weight | 1100 x 682 x 1120mm, Approx. 220kg |
| Power supply | AC220 to 240V/11A |
| Standard accessories | 2 x Universal shaking platform MT-6040 (*2), 4 x Vibration absorbing pad, 4 x Spacer, 1 x Drain-reception plastic bottle (With hose) |

(*1)The temp. of defrost function just run is not included.

(*2)Clamps not included.

•Temp. range and Temp. stability may not meet the specs when frost by continuous operation. Required to defrost it regularly.

Bioshaker BR-180LF/BR-180LF-70RT

Split doors enable to take large vessels in and out easily. Can be modified to Forced ventilation system of carbon dioxide. Shaking width 70mm has an effective on the cultivation of Actinomycete (BR-180LF-70RT).

Dimensions --> P.038 OD-Monitor Turbidimeter while shaking --> P.034 Forced ventilation system of carbon dioxide --> P.046



BR-180LF



BR-180LF-70RT

Features

- Max. speed 400r/min.
- Split doors enable easy to take large vessels.
- Capacity : Erlenmeyer flask 500mL x 24 pcs, 5L x 4 pcs



Easy to take large vessels in and out.

"Clamshell split doors" enable to take large vessels in and out easily. Equipped casters enable easy to move (BR-180LF). The Max. speed 400r/min. 70RT (orbital shaking/shaking width 70mm) can be shaken up to 250r/min that is optimum for especially Cultivation of microbe with Erlenmeyer flask.

Applications

- Cultivation of E.Coli and yeast
- Cultivation of Actinomycete (70RT)

Sliding retractable control panel

Setting for Temperature, Shaking speed and Operation program with simply.



Drawer type control panel

| Model | BR-180LF | BR-180LF-70RT |
|-------------------------------|--|--|
| Shaking motion | Switchable Reciprocal/Orbital shaking | Orbital shaking |
| Temperature range | 4°C to 70°C | 4°C to 70°C |
| Temp. control accuracy (*1) | ±0.3°C to ±1.0°C | |
| Speed range (*2) | 25 to 250 r/min, 25 to 400r/min at 25mm/orbital | 25 to 250r/min |
| Shaking width | Switchable 25/50mm (Default: 25mm) | 70mm |
| Platform size/Capacity | 700 x 500mm (Max. inside height: 420mm), Erlenmeyer flask 500mLx24, 5Lx4 | 680 x 480mm (Max. inside height: 420mm), Erlenmeyer flask 500mLx24, 5Lx4 |
| Door type | Clamshell split door (Divided into Upper and Lower, Upward opening and Front opening) | |
| Maximum load | approx. 10kg (Platform not included) | |
| Ambient temp. range (*1) | 5°C to 30°C | |
| Heating/Cooling method | Heater 700W/Compressor 125W | |
| Program functions | 1 program, 9 segments (4 combinations can be saved): Temp., Shaking speed and Time/per segment (Setting range: 00h00min to 99h59min), Temp. gradient function (*3) | |
| Safety devices/functions | Leakage/Overcurrent breaker, Compressor overheating cut-off circuit, Alarm of abnormal sample temp. (Set as desired), Overload/Overvoltage protection, Fuse, Current limiter, Detecting the door opening and closing, Self-diagnosis function for temperature sensor error, Nonvolatile memory error, Breaking device when motor runaway | |
| Dimensions (W x D x H)/Weight | 1110 x 716 x 990mm, Approx. 180kg | 1110 x 716 x 930mm, Approx. 190kg |
| Power supply | AC220 to 240V/12A | AC220 to 240V/12A |
| Standard accessories | 1 x Universal shaking platform MT-7050 (LF)/MT-6848 (LF-70RT) (*4), 1 x Drip tray (LF)/Drain-reception plastic bottle with hose (LF-70RT), 4 x Vibration absorbing pad, 4 x Spacer (LF-70RT) | |

(*1)The temp. of defrost function just run is not included.



(*2)There is the limitation depending the volume of vessels and the shaking diameter. See the figure "Note" on the right page.

(*3)The transfer time of temperature can be set up. (*4)Clamps not included.


•Temp. range and Temp. stability may not be meet the specs when frost by continuous operation. Required to defrost it regularly.

Optional accessories and Capacity/Qty of Vessels

Shaking platform and Capacity/Qty of Flasks

| Model | | MT-7050 (Refer MT-6848 below for 70RT.) | | |
|---|------------------------|--|---------------------------------------|-----|
|  | | It is used together Clamps, Angle racks and "OD-Monitor" series (Clamps of Erlenmeyer flask and Sakaguchi flask available). It can hold flasks tightly than Spring shaking platform. | | |
| External dimension / Weight | | 700x500x22Hmm, Apporx.6kg | | |
| Capacity of Flask | Erlenmeyer flask 50mL | 96 | Erlenmeyer flask 2L | 8 |
| | Erlenmeyer flask 100mL | 48 | Erlenmeyer flask 3L | 6 |
| | Erlenmeyer flask 200mL | 35 | Erlenmeyer flask 5L / Disposable 3&5L | 4/5 |
| | Erlenmeyer flask 250mL | 28 | Sakaguchi flask 500mL | 18 |
| | Erlenmeyer flask 300mL | 28 | Sakaguchi flask 1L | 12 |
| | Erlenmeyer flask 500mL | 24 | Sakaguchi flask 2L | 6 |
| | Erlenmeyer flask 1L | 13 | Sakaguchi flask 3L | 5 |
| Model | | MT-6848 (for 70RT, Pre-assembled) | | |
|  | | Designed for BR-180LF-70RT. Fixed with Optional clamps. Not equipped with Drawer mechanism due to designed for High load. | | |
| External dimension / Weight | | 680x480x22Hmm, Apporx.5.2kg | | |
| Capacity | | same as MT-7050 (see above) | | |

Spring net shaking platform and Capacity/Qty of Vessels (Mounted on MT-7050 for use)

| Model | MR-4030L (Not used with 70RT) | | |
|---|---|---|--|
|  | The usable Shaking platform that Vessels can be easily mounted on regardless its shapes. Fixed with the screws to Platform MT-7050 for using. Up to Two sets can be mounted. (See Left figure.) | | |
| Spring net Pitch | 20mm | | |
| Spring net stage | 2 stage (Height is changeable.) | | |
| External dimension / Weight | 400x300x108Hmm, Apporx.2kg | | |
| Capacity of Vessels (°v) | φ16 mm test tubes (*1) | 96pcs (When Two sets mounted: 2 x 96 = 192pcs). | |
| | Disposable Centrifuge tube for 50 mL (*1) | 48pcs (When Two sets mounted: 2 x 48 = 96pcs). | |
| | 50mL | 24pcs (When Two sets mounted: 2 x 24 = 48pcs). | |
| | 100mL | 20pcs (When Two sets mounted: 2 x 20 = 40pcs). | |
| | 200mL | 12pcs (When Two sets mounted: 2 x 12 = 24pcs). | |
| | 250mL | 9pcs (When Two sets mounted: 2 x 9 = 18pcs). | |
| | 300mL | 9pcs (When Two sets mounted: 2 x 9 = 18pcs). | |
| | 500mL | 6pcs (When Two sets mounted: 2 x 6 = 12pcs). | |
| | 1L | 3pcs (When Two sets mounted: 2 x 3 = 6pcs). | |
| | 2L | 2pcs (When Two sets mounted: 2 x 2 = 4pcs). | |
| | 3L | 1pcs (When Two sets mounted: 2 x 1 = 2pcs). | |
| | Sakaguchi flask 500mL | 6pcs (When Two sets mounted: 2 x 6 = 12pcs). | |

(*1)The qty when Angle rack tilted 45°. φ16mm or less test tube should fixed with the stand per net.
(*2)The Capacity/Qty of vessels are less than those of Universal platform.

Machine height of BR-180LF and BR-180LF-70RT





BR-180LF-70RT equipped with Vibration absorbing pad but not casters and adjusters unlike BR-180RT. Due to it the height is slightly lower than BR-180LF.

[Note] limitation of Shaking speed by Shaking width and Vessel sizes. (70RT)

| Shaking width / Shaking motion | 70mm, Orbital |
|--------------------------------|---------------|
| ~500mL | 250r/min |
| 1L | 250r/min |
| 2L | 200r/min |
| 3L | 160r/min |
| 5L | 150r/min |

Clamps (Mounted on MT-7050/6848 and AF-2070/2570 for use.)

| | Vessels / Capacity | | Model |
|--|--------------------|-----------------|------------|
| | Erlenmeyer flask | Sakaguchi flask | |
|  | 50mL | | CF-0050 |
| | 100mL | | CF-0100 |
| | 200mL | | CF-0200 |
| | 250mL | | CF-0250 |
| | 300mL | | CF-0300 |
| | 500mL | | CF-0500 |
| | 1L | | CF-1000 |
| | 2L | | CF-2000 |
| | 3L | | CF-3000 |
| | 5L | | CF-5000 |
|  | 3&5L (*) | | CF-3000DSP |
| | 500mL | | SF-0500 |
| | 1L | | SF-1000 |
| | 2L | | SF-2000 |
| | 3L | | SF-3000 |

Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet.



MT-7050 with multiple clamps mounted.

(*)For Disposable Erlenmeyer flasks (Coming 3L and Thomson 5L).

Angle rack for Test tubes (Mounted on MT-7050 for use)

| Model | AT-1250 | AT-1620 | AT-1650 | AT-1820 |
|-------------------------|--------------|--------------|--------------|--------------|
| Capacity | φ12mm x50 | φ16.5mm x20 | φ16.5mm x50 | φ18mm x20 |
| Mounting qty on MT-7050 | 4 (*) | 6 (*) | 4 (*) | 6 (*) |
| External dimension (mm) | 250x100x95H | 285x70x80H | 286x110x100H | 320x70x80H |
| Weight | Apporx.0.8kg | Apporx.0.6kg | Apporx.0.9kg | Apporx.0.7kg |

(*)The qty when Angle rack with 160mm test tube tilted 45°.
•It cannot be used with BR-180LF-70RT.

Angle rack for Disposable Centrifuge tubes (Mounted on MT-7050 for use)

| Model | AT-3518 |
|-----------------------------|--------------------------------|
| Capacity | 50mLx18 |
| Mounting qty on MT-7050 | 4 (When Angle rack tilted 45°) |
| External dimension / Weight | 285x110x102Hmm, Apporx.0.8kg |
| Remarks | Be sure to be used below 50°C. |

•It cannot be used with BR-180LF-70RT.

Angle rack for Erlenmeyer flask (Mounted on MT-7050 / Need Clamps)

| Model | AF-2070 | AF-2570 |
|-------------------------|--|--|
| Capacity | 50mL:4, 100mL:3, 200 to 300mL:2, 500mL:1 | 50mL:5, 100mL:4, 200mL:3, 250 to 500mL:2 |
| Mounting qty on MT-7050 | 4 (*) | 4 (*) |
| External dimension | 250x110x50Hmm | 300x110x50Hmm |
| Weight | Apporx.0.8kg | Apporx.0.9kg |

(*)The qty when Angle rack tilted 30°.
•Clamps are available as an option.
•It cannot be used with BR-180LF-70RT.

[Note] limitation of Shaking speed by Shaking width and Vessel sizes. (BR-180LF)

| Shaking width | 25mm | | | 50mm | | | |
|----------------|-----------------|---------|-----|------------|-----|----------|-----|
| Shaking motion | Recip- rocal | Orbital | | Reciprocal | | Orbital | |
| r/min | 250 | 250 | 400 | 200 | 250 | 200 | 250 |
| ~500mL | ○ | ○ | ○ | ○ | ○ | ○ | ○ |
| 1L | ○ | ○ | x | ○ | ○ | ○ | ○ |
| 2L | ○ | ○ | x | ○ | x | ○ | ○ |
| 3L | ○ | ○ | x | ○ | x | ○ | x |
| 5L | ○ | ○ | x | ○ | x | 170r/min | x |


Bioshaker GBR-200/GBR-300

Enables to use widely for Cultivation of microbes such as Yeast, E.Coli. etc. that require temp. control. High-end model that features for stackable, wide range of temp. range and speed. Cut power consumption by 75% with the conventional one!

Dimensions --> P.038 OD-Monitor Turbidimeter while shaking --> P.034 LED Light irradiation unit --> P.037



Stackable up to Two levels



Swing lift-up door and Stackable up to Two levels.

Swing lift-up door and Stackable up to Two levels contribute space-saving while keeping large capacity. It can be shaken up to 300r/min (200r/min at the upper of stacked units) that is optimum for shaking culture.



Features

- Max. speed 200r/min (GBR-200) and 300r/min (GBR-300)
- Capacity : Erlenmeyer flask 500mL x 24 pcs, 5L x 4 pcs
- For as a laboratory bench. Stackable up to Two stage.

Applications

- Cultivation of microbes such as E.Coli. etc.
- Cultivation of Yeast, Psychrophile and Insect cell
- Cultivation of Thermophile and Hybridization (FH).

The top board can be practically used for Lab bench.

Adopted the same material of Lab bench top board. The top board height is almost the same that of Lab bench.



Contributes save energy.

When set above 10°C, the compressor not run, while the heater runs for temp. control that is energy saving. It cuts power consumption by 75% with the conventional one!

| Model | GBR-200 | GBR-300 |
|-------------------------------|---|---|
| Shaking motion | Switchable Reciprocal/Orbital shaking | |
| Temperature range | 4°C to 80°C | |
| Temp. control accuracy (*1) | ±0.3°C to ±1.0°C | |
| Speed range (*2) | 25 to 200r/min | 25 to 300r/min (200r/min at the upper of stacked units) |
| Shaking width | Switchable 25/50mm (Default: 25mm) | |
| Platform size/Capacity | 700 x 500mm (Max. inside height: 385mm), Erlenmeyer flask 500mLx24, 5Lx4 | |
| Door type | Swing lift-up door | |
| Maximum load | approx. 10kg (Platform not included) | |
| Ambient temp. range (*1) | 5°C to 30°C | |
| Heating/Cooling method | Heater 700W/Compressor 125W | |
| Other functions | 1 pattern, 9 segments (4 combinations can be saved): Temp., Shaking speed and Time/per segment (Setting range: 00h00min to 99h59min), Wait function (*3) | |
| Safety devices/functions | Leakage/Overcurrent breaker, Compressor overheating cut-off circuit, Alarm of abnormal sample temp. (Set as desired), Overload/Overvoltage protection, Fuse, Current limiter, Detecting the door opening and closing, Self-diagnosis function for temperature sensor error, Nonvolatile memory error, Breaking device when motor runaway, Self-check/lamp of filter | |
| Dimensions (W x D x H)/Weight | 1200 x 800 x 810mm (Including Top board: 1200 x 750 x 25mm), Approx.235kg | |
| Power supply | AC100V/12A | |
| Standard accessories | 1 x Universal shaking platform MT-7050 (*4), 2 x 2 types Spacer, 2 x Fixing bracket for stackable (*4) | |

(*1)The temp. of defrost function just run is not included.

(*2)There is the limitation depending the volume of vessels and the shaking width. See the figure "Note" on the right page.


(*3)Shaking will start when reaching the preset temp. (*4)It comes with by the request at order.

(*4)Clamps not included.


•Temp. range and Temp. stability may not be meet the specs when frost by continuous operation. Required to defrost it regularly.

Optional accessories and Capacity/Qty of Vessels

Shaking platform and Capacity/Qty of Flasks


| Model | | MT-7050 | | |
|---|------------------------|--|---------------------------------------|-------|
|  | | It is used together Clamps, Angle racks and "OD-Monitor" series (Clamps of Erlenmeyer flask and Sakaguchi flask are available). It can hold flasks tightly than Spring shaking platform. | | |
| External dimension / Weight | | 700x500x22Hmm, Apporx.6kg | | |
| Capacity of Flask | Erlenmeyer flask 50mL | 96 | Erlenmeyer flask 2L | 8 |
| | Erlenmeyer flask 100mL | 48 | Erlenmeyer flask 3L | 6 |
| | Erlenmeyer flask 200mL | 35 | Erlenmeyer flask 5L / Disposable 3&5L | 4 / 5 |
| | Erlenmeyer flask 250mL | 28 | Sakaguchi flask 500mL | 18 |
| | Erlenmeyer flask 300mL | 28 | Sakaguchi flask 1L | 12 |
| | Erlenmeyer flask 500mL | 24 | Sakaguchi flask 2L | 6 |
| | Erlenmeyer flask 1L | 13 | Sakaguchi flask 3L | 5 |

Spring net shaking platform and Capacity/Qty of Vessels (Mounted on MT-7050 for use)

| Model | MR-4030L | | |
|---|---|--|---|
|  | <p>The usable Shaking platform that Vessels can be easily mounted on regardless its shapes. Fixed with the screws to Platform MT-7050 for using. <u>Up to Two sets can be mounted.</u> (See Left figure.)</p> | | |
| Spring net Pitch | 20mm | | |
| Spring net stage | 2 stage (Height is changeable.) | | |
| External dimensions / Weight | 400x300x108Hmm, Apporx.2kg | | |
| Capacity of Vessels (*) | Erlenmeyer flask | φ16 mm test tubes (*) | 96pcs (When Two sets mounted: $2 \times 96 = 192$ pcs). |
| | | Disposable Centrifuge tube for 50 mL (*) | 48pcs (When Two sets mounted: $2 \times 48 = 96$ pcs). |
| | | 50mL | 24pcs (When Two sets mounted: $2 \times 24 = 48$ pcs). |
| | | 100mL | 20pcs (When Two sets mounted: $2 \times 20 = 40$ pcs). |
| | | 200mL | 12pcs (When Two sets mounted: $2 \times 12 = 24$ pcs). |
| | | 250mL | 9pcs (When Two sets mounted: $2 \times 9 = 18$ pcs). |
| | | 300mL | 9pcs (When Two sets mounted: $2 \times 9 = 18$ pcs). |
| | | 500mL | 6pcs (When Two sets mounted: $2 \times 6 = 12$ pcs). |
| | | 1L | 3pcs (When Two sets mounted: $2 \times 3 = 6$ pcs). |
| | | 2L | 2pcs (When Two sets mounted: $2 \times 2 = 4$ pcs). |
| | | 3L | 1pcs (When Two sets mounted: $2 \times 1 = 2$ pcs). |
| | Sakaguchi flask 500mL | 6pcs (When Two sets mounted: $2 \times 6 = 12$ pcs). | |

(*)1) The qty when Angle rack tilted 45°. φ16mm or less test tube should fixed with the stand per net.
 (*)2) The Capacity/Qty of vessels are less than those of Universal platform.


Junction pipe




| Model | JCP-3687 | | | |
|-------------|---|--|--|--|
| Description | The pipe (w/o Lid) for connection to the outside. It is for installing such as the sensor of recorder inside the chamber. Electrical equipment should not be placed inside the chamber. | | | |
| Remarks | Be sure be used below 50°C. | | | |

Clamps (Mounted on MT-7050 and AF-2070/2570 for use.)

| Vessels / Capacity | | Model |
|--------------------|------------|---------|
| Erlenmeyer flask | 50mL | CF-0050 |
| | 100mL | CF-0100 |
| | 200mL | CF-0200 |
| | 250mL | CF-0250 |
| | 300mL | CF-0300 |
| | 500mL | CF-0500 |
| | 1L | CF-1000 |
| | 2L | CF-2000 |
| | 3L | CF-3000 |
| | 5L | CF-5000 |
| 3&5L (*) | CF-3000DSP | |
| Sakaguchi flask | 500mL | SF-0500 |
| | 1L | SF-1000 |
| | 2L | SF-2000 |
| | 3L | SF-3000 |



Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet.



MT-7050 with multiple clamps mounted.

(*)For Disposable Erlenmeyer flasks (Corning 3L and Thomson 5L).

Angle rack for Test tubes (Mounted on MT-7050 for use)

| Model | AT-1250 | AT-1620 | AT-1650 | AT-1820 |
|-------------------------|--------------|--------------|--------------|--------------|
| Capacity | φ12mm x50 | φ16.5mm x20 | φ16.5mm x50 | φ18mm x20 |
| Mounting qty on MT-7050 | 4 (*) | 6 (*) | 4 (*) | 6 (*) |
| External dimension (mm) | 250x100x95H | 285x70x80H | 286x110x100H | 320x70x80H |
| Weight | Apporx.0.8kg | Apporx.0.6kg | Apporx.0.9kg | Apporx.0.7kg |

(*)The qty when Angle rack with 160mm test tube tilted 45°.

Angle rack for Disposable Centrifuge tubes (Mounted on MT-7050 for use)

| Model | AT-3518 | | | |
|-----------------------------|--------------------------------|--|--|--|
| Capacity | 50mLx18 | | | |
| Num. of mounting to MT-7050 | 4 (When Angle rack tilted 45°) | | | |
| External dimension / Weight | 285x110x102Hmm, Apporx.0.8kg | | | |
| Remarks | Be sure be used below 50°C. | | | |

Angle rack for Erlenmeyer flasks (Mounted on MT-7050 for use)

| Model | AF-2070 | AF-2570 |
|-------------------------|--|--|
| Capacity | 50mL:4, 100mL:3, 200 to 300mL:2, 500mL:1 | 50mL:5, 100mL:4, 200mL:3, 250 to 500mL:2 |
| Mounting qty on MT-7050 | 4 (*) | 4 (*) |
| External dimension | 250x110x50Hmm | 300x110x50Hmm |
| Weight | Apporx.0.8kg | Apporx.0.9kg |

(*)1) The qty when Angle rack tilted 30°. Clamps are available as an option.

[Note] limitation of Shaking speed by Shaking width and Vessel sizes.

| Shaking width | 25mm | | | | 50mm | | | |
|---------------|------------|-----|---------|-----|------------|-----|----------|-----|
| | Reciprocal | | Orbital | | Reciprocal | | Orbital | |
| r/min | 200 | 300 | 200 | 300 | 200 | 300 | 200 | 300 |
| ~500mL | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ |
| 1L | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ |
| 2L | ○ | ○ | ○ | ○ | ○ | x | ○ | ○ |
| 3L | ○ | ○ | ○ | ○ | ○ | x | ○ | x |
| 5L | ○ | ○ | ○ | ○ | 170r/min | x | 170r/min | x |

Bioshaker BR-3000LF Single platform /BR-3000LF Double platform

Enables to use widely for Cultivation of microbes such as Yeast, E.Coli. etc. that require temp. control. Double platform realizes Double capacity.

Dimensions --> P.038 OD-Monitor Turbidimeter while shaking --> P.034 LED irradiation unit --> P.037



BR-3000LF Double platform

Features

- Single platform and Double platform
- Capacity : Erlenmeyer flask 500mL x 32 pcs (Single), 64 pcs (Double)
- Shaking width adjustable in stepless between 10mm and 50mm.

Applications

- Cultivation of microbes such as E.Coli. etc.
- Cultivation of Yeast, Psychrophile and Insect cell
- The other samples that require shaking and temp. control.

The magnificent max. capacity model with the double platform


This is the magnificent max. capacity model with the double platform among TAITEC incubator shakers by the repeated improvements such as Shaking platform drawer lever. Up to 1 L flask can be placed on Double platform that vessel qty doubled as Single platform. Note: The upper platform should be removed when above 2L Erlenmeyer flask and 1L Shake flask are placed. (Ask us for details). Since Switchable shaking motion and Shaking width can be adjusted in stepless the optimum efficiency of cultivation can be obtained depending a vessel.



Shaking platform drawer lever

Equipped with the remote terminal for the temp. control.

Connecting with Program unit PU-5 (see the right page) enables the remote control setting and the program operation.



[Lower left back]


| Model | BR-3000 Single platform | BR-3000 Double platform |
|---------------------------------|---|---|
| Temperature range | 4°C to 70°C | |
| Temp. control accuracy (*1) | ±0.3°C to ±1.0°C | |
| Shaking motion | Switchable Reciprocal/Orbital shaking | |
| Speed range | 25 to 250r/min | 25 to 160r/min |
| Shaking width | 10 to 50mm (stepless) | |
| Platform size/Capacity | 800 x 600mm (Max. inside height: 605mm), Erlenmeyer flask 500mLx32, 5Lx6 | 800 x 600mm (Limitation flask volume per platform: 1L Erlenmeyer flask or 500mL Sakaguchi flask) Erlenmeyer flask 500mLx64, 1Lx36 |
| Door | Double swing door | |
| Maximum load | Approx. 17kg (Platform not included) | Approx. 17kg/per platform (Platform not included) |
| Ambient temp. range (*1) | 5°C to 35°C | |
| Display (Temp./Speed) | Digital display (Switchable Set temp./Chamber temp.) | |
| Heating/Cooling method | Heater 800W/Compressor 125W | |
| Other functions | Shaking duration integrator (0.1h to 999.9h, Automatic resetting function *Required Program unit PU-5 or 6/Option to use), Chamber lamp | |
| Safety devices/functions | Leakage/Overcurrent breaker, Compressor overheating cut-off circuit, Alarm of abnormal sample temp. (Set as desired), Display of deviation from set value (upper / lower limit), Self-diagnosis function for sensor error | |
| External Dimensions (W x D x H) | 1325 x 862 x 1234mm | |
| Weight | Approx. 310kg | Approx. 330kg |
| Power supply | AC220 to 240V / 12A | |
| Standard accessories | 1 x Universal shaking platform MT-8060A (Single type *Double type comes with 2 pcs of platform) (*2), 4 x Vibration absorbing pad, 4 x Spacer, 1 x Drain-reception plastic bottle (With hose) | |

(*1)The temp. of defrost function just run is not included. (*2)Clamps not included.

•Temp. range and Temp. stability may not be meet the specs when frost by continuous operation. Required to defrost it regularly.


Optional accessories and Capacity/Qty of Vessels

Shaking platform and Capacity/Qty of Flasks

| Model | | MT-8060A | | |
|---|------------------------|--|---|-----|
|  | | It is used together Optional Clamps, Angle racks and "OD-Monitor" series. 2 sets used with Double shaking platform type (included in BR-3000LF). | | |
| External dimension / Weight | | 800x600x45Hmm, Apporx.9kg | | |
| Capacity of Flask | Erlenmeyer flask 50mL | 98 | Erlenmeyer flask 2L | 12 |
| | Erlenmeyer flask 100mL | 72 | Erlenmeyer flask 3L | 9 |
| | Erlenmeyer flask 200mL | 63 | Erlenmeyer flask 5L/ Disposable3&5L | 6/7 |
| | Erlenmeyer flask 250mL | 46 | Sakaguchi flask 500mL | 27 |
| | Erlenmeyer flask 300mL | 46 | Sakaguchi flask 1L | 18 |
| | Erlenmeyer flask 500mL | 32 | Sakaguchi flask 2L | 12 |
| | Erlenmeyer flask 1L | 18 | Sakaguchi flask 3L | 9 |
| Capacity/Qty in Double shaking platform type. | | | | |
| Capacity of Flasks (Two stages in total) | Erlenmeyer flask 50mL | 196 | Erlenmeyer flask 2L (*) | 12 |
| | Erlenmeyer flask 100mL | 144 | Erlenmeyer flask 3L (*) | 9 |
| | Erlenmeyer flask 200mL | 126 | Erlenmeyer flask 5L/ Disposable3&5L (*) | 6/7 |
| | Erlenmeyer flask 250mL | 92 | Sakaguchi flask 500mL | 54 |
| | Erlenmeyer flask 300mL | 92 | Sakaguchi flask 1L (*) | 18 |
| | Erlenmeyer flask 500mL | 64 | Sakaguchi flask 2L (*) | 12 |
| | Erlenmeyer flask 1L | 36 | Sakaguchi flask 3L (*) | 9 |

(*) It can be used on the lower platform when that of upper removed.

Spring net shaking platform and Capacity/Qty of Vessels (Mounted on MT-8060A for use)

| Model | MR-4030L | | |
|---|---|--|--|
|  | The usable Shaking platform that Vessels can be easily mounted on regardless its shapes. Fixed with the screws to Platform MT-8060A for using. Up to Four sets can be mounted. (See Left figure.) 8 sets in total can be mounted on Double platforms. | | |
| Spring net Pitch | 20mm | | |
| Spring net stage | 2 stage (Height is changeable.) | | |
| External dimension / Weight | 400x300x108Hmm, Apporx.2kg | | |
| Capacity of Vessels (*) | φ16 mm test tubes (*) | 96pcs (When Four sets mounted: 4 x 96 = 384pcs) x 2 stages | |
| | Disposable Centrifuge tube for 50 mL (*) | 48pcs (When Four sets mounted: 4 x 48 = 192pcs) x 2 stage | |
| | 50mL | 24pcs (When Four sets mounted: 4 x 24 = 96pcs) x 2 stage | |
| | 100mL | 20pcs (When Four sets mounted: 4 x 20 = 80pcs) x 2 stage | |
| | 200mL | 12pcs (When Four sets mounted: 4 x 12 = 48pcs) x 2 stage | |
| | 250mL | 9pcs (When Four sets mounted: 4 x 9 = 36pcs) x 2 stage | |
| | 300mL | 9pcs (When Four sets mounted: 4 x 9 = 36pcs) x 2 stage | |
| | 500mL | 6pcs (When Four sets mounted: 4 x 6 = 24pcs) x 2 stage | |
| | 1L | 3pcs (When Four sets mounted: 4 x 3 = 12pcs) x 2 stage | |
| | 2L (*) | 2pcs (When Four sets mounted: 4 x 2 = 8pcs) | |
| | 3L (*) | 1pcs (When Four sets mounted: 4 x 1 = 4pcs) | |
| | Sakaguchi flask 500mL | 6pcs (When Four sets mounted: 4 x 6 = 24pcs) x 2 stage | |

(*)1) The qty when Angle rack tilted 45° φ16mm or less test tube should fixed with the stand per net.

(*)2) The Capacity/Qty of vessels are less than those of Universal platform.

(*)3) It can be used on the lower platform when that of upper removed.



LED irradiation unit (Mounted to MR-4030 for use)

| Model | LC-LED450W | |
|----------------------------|------------|--|
| LED Color | White | |
| Peak wavelength | 450nm | |
| Dimensions of Light source | 340x240mm | |

• See details on 47 page.

● Protuberances not included in Dimensions. ● Vessels of photo not included.

Clamps (Mounted on MT-8060A and AF-2070/2570 for use)

| | Vessels / Capacity | Model |
|--|--------------------|------------|
|  Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet. | 50mL | CF-0050 |
| | 100mL | CF-0100 |
| | 200mL | CF-0200 |
| | 250mL | CF-0250 |
| | 300mL | CF-0300 |
| | 500mL | CF-0500 |
| | 1L | CF-1000 |
| | 2L (*) | CF-2000 |
| | 3L (*) | CF-3000 |
| | 5L (*) | CF-5000 |
|  MT-8060A with multiple clamps mounted. | 3&5L (*) | CF-3000DSP |
| | 500mL | SF-0500 |
| | 1L (*) | SF-1000 |
| | 2L (*) | SF-2000 |
| | 3L (*) | SF-3000 |

(*)1) It cannot be used in the upper platform and can be used on the lower platform when that of upper removed.

(*)2) For Disposable Erlenmeyer flasks (Corning 3L and Thomson 5L).

Angle rack for Test tubes (Mounted on MT-8060A for use)

| Model | AT-1250 | AT-1620 | AT-1650 | AT-1820 |
|--------------------------|--------------|--------------|--------------|--------------|
| Capacity | φ12mm x50 | φ16.5mm x20 | φ16.5mm x50 | φ18mm x20 |
| Mounting qty on MT-8060A | 10 (*) | 16 (*) | 8 (*) | 6 (*) |
| External dimension (mm) | 250x100x95H | 285x70x80H | 286x110x100H | 320x70x80H |
| Weight | Apporx.0.8kg | Apporx.0.6kg | Apporx.0.9kg | Apporx.0.7kg |

(*) The qty when Angle rack with 160mm test tube tilted 45°.

Angle rack for Disposable Centrifuge tubes (Mounted on MT-8060A for use)

| Model | AT-3518 |
|------------------------------|--------------------------------|
| Capacity | 50mLx18 |
| Num. of mounting to MT-8060A | 8 (When Angle rack tilted 45°) |
| External dimension / Weight | 285x110x102Hmm, Apporx.0.8kg |
| Remarks | Be sure be used below 50°C. |

Angle rack for Erlenmeyer flask (Mounted on MT-8060A for use)

| Model | AF-2070 | AF-2570 |
|--------------------------|--|--|
| Capacity | 50mL:4, 100mL:3, 200 to 300mL:2, 500mL:1 | 50mL:5, 100mL:4, 200mL:3, 250 to 500mL:2 |
| Mounting qty on MT-8060A | 10 (*) | 5 (*) |
| External dimension | 250x110x50Hmm | 300x110x50Hmm |
| Weight | Apporx.0.8kg | Apporx.0.9kg |

(*)1) The qty when Angle rack tilted 30°. • Clamps are available as an option.

Program units

Connecting to the remote temperature setting terminal of BR-300LF series enables remote setting and program setting of temperature (The modification of BR is required for the shaking control).

| Model | PU-5 | PU-6 |
|-------------|--|--|
| Description | Remote program control of Temperature. | Remote program control of Temperature and Shaking ON/OFF. |
| Main Spec. | Storage programs: 2 program (8 segments per program) | Storage programs: 2 program (8 segments per program) Time signal (ON/OFF): Each 2 programs |

OD-Monitor OD-Monitor A&S/B&L/C&T

The industry's first Non-contact turbidimeter realizes to measure the turbidity of sample at OD₆₀₀ while shaking. Adopted among the hundreds of R&D institutions. Operated with Bioshaker (Incubator shaker) in combination.

"Bioshaker BR series" --> P.018-032 Shaker "NR series" --> P.062-066



Features

- The turbidity can be measured at OD₆₀₀ noncontact while shaking.
- Three models for Large / Small Erlenmeyer flask and Test tube.
- Operated with Bioshaker (Incubator shaker) in combination

Applications

- Labor-saving acquisition of growth curve.
- Screening for Microbes such as beneficial bacteria that are sampled difficulty.
- Automation of IPTG induction [B&L with options].

Released from the cumbersome turbidity measurement.

Enables you to perform the shaking culture and the measurement of turbidity efficiently at the same time and released from the cumbersome turbidity measurement. Not has a harmful effect on the growth and metabolism of fungus thanks to no-stop stop shaking.

Operated with Bioshaker in combination.

Fix ODSensor with the clamp on the platform inside Bioshaker BR-series (w/o BR-180LF-70RT) and put the vessel on it. It can be also used in Shaker NR series.

Wide variety of customization.

We can customize it widely and flexibly per your request such as for Strictly anaerobic bacteria, Thermophile, Ultra-high concentration, Water bath and so on. Feel free to ask us.

| Model | OD-Monitor A&S | OD-Monitor B&L | OD-Monitor C&T |
|-----------------------------------|--|---|---|
| Measurement method/range | Light transmission, 0.00 to 2.55 (OD600) | Light transmission/Reflected light, 0.00 to 9.99 (OD600) | Light transmission, 0.00 to 2.55 (OD600) |
| Vessel and Fluid volume | Vessel: Glass Erlenmeyer flask (*2) Size: 100/200/300/500mL Fluid volume: 20/40/60/100mL | Vessel: Glass Erlenmeyer flask (*2) Size: 1/2/3L Fluid volume: 350/700/1000mL | Vessel: Glass Test tube Size: φ16.5/18.0mL Fluid volume: 5mL or 10mL |
| Shaking conditions | Orbital shaking 80 to 400r/min (Shaking width 25mm) 80 to 200r/min (Shaking width 30/40mm) 80 to 250r/min (Shaking width 50mm) | Orbital shaking 100 to 250r/min (Shaking width 25/50mm) | Reciprocal shaking 100 to 250r/min (Shaking width 25/50mm) Angle of Test tube 90/45/30° |
| Measurement intervals | 1/5/10/30/60/120/240min | | |
| Measurement channel | 1 to 8 (*1) | 1 to 4 (*1) | 1 to 8 (*1) |
| Ambient temp. range | 15 to 50°C | | |
| Other functions | Alarm for specific value at OD, Transformation to Spectrophotometer in measurement value, USB flash drive (CSV file), Operation for pump unit (B&L) when above specified OD value. | | |
| Standard accessories for ODBox | 1 x Flat cable (2m) * For A&S, 1 x Power code | | |
| Standard accessories for ODSensor | 1xVibration absorbing pad of Flask (large and small), 1 x Flat cable (0.5m) | 1 x Flask stand (1L), 1 x Flat cable (2m) | 1 x Magnet adapter (2 types), 4 x Flat cable (2m) |

(*1)A&S/B&L equipped with 1 channel per 1 set of ODSensor. Some number of ODSensor should be required depending the number of samples. Although C&T equipped with 8 channels per 1 set of ODSensor it can be only connected with 1 set of ODBox. ODSensor with maximum number of channel cannot be put in Bioshaker depending its size. Ask us for details.
(*2)Flask clamp not come with. It is as an option.


The capacity (number of installation and measurement) of ODSensor in Bioshaker BR-series.

A&S and B&L: In BR-series the number of measurement of the unit equal to that of installation of the unit.

C&T: The number of measurement to the number of installation is by 8 times.

| BR series | ODSensor-S | ODSensor-L | ODSensor-T (at 45°) |
|----------------|-------------------|------------------|---------------------|
| BR-21/22/23 | 3 sets | 1set (Up to 2L) | 2 sets = 2 x C&T |
| BR-40/41/42/43 | 5 sets | 2 sets (*1) | 3 sets = 3 x C&T |
| BR-300LF | 8 sets | 4 sets | 5 sets = 5 x C&T |
| BR-180LF/GBR | 16 sets = 2 x A&S | 5 sets = 2 x B&L | 12 sets = 12 x C&T |
| BR-3000LF | 20 sets = 3 x A&S | 8 sets = 2 x B&L | 14 sets = 14 x C&T |

(*1)Up to 2L of Erlenmeyer flask can be used with the BR-40LF.



Combination of OD-Monitor C&T and BR-23FP

Combination of OD-Monitor A&S and BR-43FL

OD-Monitor Turbidostat System A / System B

*Enables Continuous culture with OD-Monitor.
Two types for large and small capacity available.*

"Bioshaker BR series" --> P.018-032 Shaker "NR series" --> P.062-066

Features

- OD-Monitor for Continuous culture system.
- Realized Small scale continuous culture compared to Mini jar fermentor.
- Used with Bioshaker BR-series in combination.

Applications

- Seed fungus not harmed during the experiment.
- Continuous culture (Producing substance)
- Study for Mutation by Continuous culture.

Continuous culture with Erlenmeyer flask.


The fungus concentration is kept constant at arbitrary OD value while measuring the turbidity at OD600 noncontact by adding medium and discharging liquid.

Operated with Bioshaker in combination.

Fix ODSensor with the clamp on the platform inside Bioshaker BR-series and put the vessel on it. It can be also used in TAITEC Shaker NR series in combination.

System configuration

The component units/parts are highlighted in the green squared. The photo of composition is for "1 Channel" of System B and also System A.



OD-Monitor Turbidostat system A



USER'S VOICE

Turbidostat system enables continuous culture at low cost maintenance.

Realizes Small scale continuous culture at low cost.

Small scale continuous culture (Min 20mL volume) that even Mini jar fermentor cannot. Drastically save the cost of culture medium and effort of washing operation.

| Model | OD-Monitor Turbidostat system A | OD-Monitor Turbidostat system B |
|--------------------------|---|---|
| Measurement method | Light transmission | Light transmission/ Reflected light |
| Measurement range | 0.00 to 2.55 (OD600) | 0.00 to 9.99 (OD600) |
| Vessel and Fluid volume | Glass Erlenmeyer flask Size: 100/200/300/500mL Fluid volume: 20/40/60/100mL | Glass Erlenmeyer flask Size: 1/2/3L Fluid volume: 350/700/1000mL |
| Shaking conditions | Orbital shaking 80 to 400r/min (Shaking width 25mm) 80 to 200r/min (Shaking width 30/40mm) 80 to 250r/min (Shaking width 50mm) | Orbital shaking 100 to 250r/min (Shaking width 25/50mm) |
| Measurement channel | 1 to 4 | |
| Ambient temp. range | 15 to 50°C | |
| Composition at 1 Channel | 1x Control unit, 1 x Measurement unit, 2 x Pump unit, 1 x Turbidostat cap (with Pre-assembled Fluid feeding tube and Needle) | |

•Flask clamps required. •This unit should be fixed to the universal platform of Bioshaker with Flask clamps (Option). See the below the below options for details •Ask us when Fraction Collector required.

Optional accesories for OD-Monitor and OD-Monitor Turbidostat System

| Product Name / Model | Remarks |
|---|---|
| Universal shaking platform MT-2925 | Required for using BR-21/22/23 series with OD-Monitor series and Turbidostat system |
| Universal shaking platform MT-4030 | Required for using BR-40/41/42/43 series with OD-Monitor series and Turbidostat system |
| Auto pump unit APU-01 | For B&L. Starts to operate when a preset OD value is reached in ODBox-B. It can be used for automatically addition of IPTG. |
| Fiat cable OD-CAV0200 | Mainly for A&S. 2m in length, 1 set required when connecting 3 to 4 sets of ODSensor-S, 2 sets required when connecting 5 to 6 sets of that, 3 sets required when connecting 7 to 8 sets of that. |
| Turbidostat Cap | Same as that of Turbidostat system. |
| Needle for OD-CC cyctem | The needle used when inserting the liquid feeding tube of the automatic pump unit into the reservoir. |

•Optional Clamps are required to fix Erlenmeyer flasks on ODSensor (except for C&T). [For 50mL Erlenmeyer flask : CF-0050] [100mL : CF-0100] [200mL : CF-0200] [300mL : CF-0300] [500mL : CF-0500] [1L : CF-1000] [2L : CF-200] [3L : CF-3000]

Easy operation OD Monitor miniphoto 518R

Enables to measure the absorbance (turbidity) by just inserting Test tube and L-shaped test tube used cultivation and reaction. Wavelength can be changed with replacement of filter.

"Bioshaker BR series" --> P.018-032 Shaker "NR series" --> P.062-66 Seesaw shaking for L shaped test tube --> P.120/124/134




Insert of Test tube

| Model | miniphoto 518R |
|------------------------|--|
| Measurement range | Absorbance: 0.000 to 2.000A, Transmittance: 1-100%T |
| Measurement wavelength | 660nm (*1) |
| Vessels | φ16.5 to 18mm Glass test tube or L-shaped test tube (*2) |
| Light source/Detector | High-Brightness LED, Photodiode |
| Display | Optical Density/Transmittance simultaneously, HOLD function. |
| Interface | Printer (Option) |
| Dimensions (W x D x H) | 187 x 106 x 95mm |
| Weight | Apporx. 880g |
| Power supply | AC100 to 240V |
| Standard accessories | 1 x Glass filter (660nm), 1 x Light shielding cap (*3), 1 x AC adapter |

(*1)Wavelength can be changed with replacement of the filter. See the below option for details. (*2)Test tube without scratches should be used. Ask us when a different diameter of Test tube required. (*3)It cannot be used with L shaped test tube. Aluminum foil or something can be substituted if light shielding required.

Option



Interference
Filter



Printer

| Product Name | Remarks |
|----------------------------|--|
| interference Filter | Band-pass filter approx. 10mm FWHM, 4 types R420/R530/R620/R660, Dedicated for miniphoto 518R. |
| L shaped test tube | φ18x120x70Hmm, 10sets |
| Printer | Printer for miniphoto 518R, included with Connecting cable and Roll paper |

Features

- For φ16.5 to 18mm Test tube and L-shaped test tube. Liquid volume 4mL
- Comes with 660nm Glass filter (Changes wavelength by replacement)
- Displays the value just inserting test tube. The display can be pinned with HOLD key.

Applications

- Measurement turbidity of Microbe culture fluid
- Absorbance measurement of Colorimetric samples
- Applications requiring measuring Absorbance/ Transmittance of sample



Measurement by inserting Test tube that has been used cultivation and reaction.



Minimum volume of the measurement is 4mL .


USER'S VOICE

Lightweight compact size and Easy Operation come in handy to me.



Light shielding cap (Standard ACC) and Printer (Option) in combination.

Dimensions



[Top view]

[Right view]

LED irradiation unit LC-LED450W LC-LED450W

Optimum used with Middle and Large sized Bioshaker in combination. The white LED contains wavelength that has an effect of photosynthesis.

Middle sized "Bioshaker BR-40/41/42/43" series --> P.022 Large sized "Bioshaker GBR-200/300" --> P.030

Features


- LED contains wavelength that has an effect of photosynthesis.
- Optimum used with Middle sized Bioshaker mainly.
- The irradiation intensity can be adjusted per 1%.

Applications

- Cultivation of Cyanobacteria and Algae
- Culture of Plant callus
- Cultivation of Plants in the laboratory

Optimum for a substitute for white fluorescent light as an upgrade.

It contains 450 nm that is the main absorption wavelength of chlorophyll as the primary peak and also contains the required intensity around 660 to 740 nm. Optimum for a substitute for white fluorescent light as an upgrade.




USER'S VOICE

It has a peak sufficiently overlapping with those of absorption of Chlorophyll.



Combination of
LC-LED450W+MR-4030

| Product Name / Model | Remarks |
|--|---|
| Spring net shaking platform MR-4030 | Remove the bottom plate and fix the LC-LED plate with screws. See p23 for vessels capacity/qty. |
| Transparent Stick sheet for LC-LED | Transparent stick sheet for protecting LC-LED plate and holding flat bottom vessels. It can be used up to 100r/min. |




LC-LED450W

Optimum for using with Middle sized Bioshaker BR and Middle sized Shaker combination.


Optimum for using with Middle sized Bioshaker BR-40LF/41FL/42FL/43FL/53FP and Middle sized Shaker NR-10/20/30/80 combination. The irradiation from the bottom can be done effectively. * Ask us for how to use it with Large sized Bioshaker in combination.



Installation of Large sized
Bioshaker



Algae cultivated with
LC-LED450



Basic system of CO₂-BR-43FL for cultivation of
Algae combined with CO₂-BR.


| Model | LC-LED450W |
|--------------------------------|---|
| LED Color | White |
| Peak wavelength | 450nm |
| Half width | ±15mm |
| Photon flux density (*) | approx. 215μmol |
| Irradiated surface average | ±15% |
| Dimensions of LED light source | 340 x 240mm |
| Dimming system | PWM control (can be adjusted per 1%) |
| Use conditions | Temperature range: 4 to 50°C, Shaking speed: Up to 200r/min |
| Weight/Power supply | Apporx. 2.5kg, AC100V/2A |
| Standard accessories | 1 x Controller |

(*)The value may be varied depending a light quantum meter.

Dimensions for Bioshaker BR-series.


BR-21/22/23 series (1) --> P.018

The drawing shows the units stacked.
Upper : BR-21UM/22UM/23UM
Lower : BR-21FP/22FP/23FP




BR-21/22/23 series (2) --> P.018

The drawing shows the units stacked.
Upper : BR-21FH/22FH/23FH
Lower : BR-21FP/22FP/23FP




BR-40LF --> P.022

The drawing shows the units stacked.
Upper : BR-40LF
Lower : BR-41/42/43series




BR-41/42/43 series --> P.022


The drawing shows the units stacked.




BR-300LF --> P.024 Double platform BR-300LF --> P.026




BR-300LF-LED --> P.024




BR-180LF --> P.028



BR-180LF-70RT --> P.028



BR-3000LF --> P.032 Double platform BR-3000LF --> P.032



G-BR-200/300 --> P.030

The drawing shows the units stacked.

