



# Water bath Shaking Water bath Immersion cooler

■ Model selection guide .....	110
■ Correspondence table of Unit Water bath	
"Thermminder" series from First-generation models to	
Latest models.....	112
■ Unit Water bath	
Thermominder SDminiN.....	114
Thermominder SDN/EXN-B.....	115
Thermominder SM-05N/SJ-07N/SX-10N/SH-10N/SP-12N.....	116
■ Bench-top Shaking Water bath	
Personal-11•SDN/EXN/SM set.....	120
■ Bench-top Shaking Water bath(Low Temp. / High Temp.)	
Personal Lt-10F•SX set.....	122
Personal H-10•SH set.....	122
■ Shaking Water bath (Large size)	
Water bath Shaker MM-10.....	124
Cool bath Shaker ML-10F.....	124

■ Combination type Shaking Water bath	
Plus Shaker EP-1.....	126
■ Program Setting Device	
Program Unit PU-5/PU-6.....	127
■ Immersion cooler	
Cool Pipe 80LF/150LF/400R.....	128
Cool Pipe 250DF/800R.....	129
■ Immersion Cooler (Related products)	
Temperature control Unit TU-100N/TU-200N.....	130

## Selection guide

Constant temperature incubator shaker

CO<sub>2</sub> incubator shaker

Shaker

Mixer Rotator Stirrer

Bead beater homogenizer Ultrasonic homogenizer

Aluminum block bath Minimize bath

Water bath Shaking water bath Immersion cooler

Hybridization oven Constant temperature chamber

Centrifugal concentrator Cold trap

Freeze dryer

Electrophoresis and blotting apparatus

Constant temperature water circulating system [Chiller]

Appendix

## Unit Water bath

## Thermominder SDminiN



## Thermominder SDN-B EXN-B



## Thermominder SM-05N/SJ-07N SX-10N/SH-10N SP-12N



## Personal-11-SDN set EXN set SM set



•Optional accessories  
Monode kit



## Personal Lt10F-SX set Personal H-10-SH set



## Shaking Water bath (large size)

## Water bath Shaker MM-10 Cool bath Shaker ML-10F



## Combination type Shaking Water bath

## Plus Shaker EP-1



## Immersion cooler

## Cool pipe series



## Unit Water bath

## Combination and Temperature range

For example, up to 70°C can actually be used in this combination. Up to approx 95°C can be raised when the water bath changed to stainless steel insulated water bath. Moreover, nearly 100°C can be raised when it can be changed to the model SX-10N (Thermominder set up to 105°C).

Heat-resistant Plastic Hood  
Heat-resistant temp. up to +100°C

Thermominder Models that can be used up to 100 °C

Tap water Below 70 °C  
Below 70°C (Heat medium for high temp. recommended for temp. above 70°C)

Plastic Water bath

Heat-resistant temp. up to 70°C

Page	Model	Features	Applications
P.114	SDminiN	•Low price	<ul style="list-style-type: none"> <li>•Various incubations such as Enzyme reaction</li> <li>•Temporarily Incubation of Culture medium</li> <li>•Inactivation of Serum</li> </ul>
P.115	SDN-B	•Low water level	
P.115	EXN-B	•Low water level, High Temp.accuracy	
P.116	SM-05N	•Economy	
P.116	SJ-07N	•Standard	
P.116	SX-10N	•High Specifications	
P.116	SH-10N	•High Temperature	
P.116	SP-12N	•External Circulation	•Circulation to Capillary of Evaporator

## Bench-top Shaking Water bath that can be used for many application

Page	Model	Features	Applications
P.120	Personal-11-SDN set	•Standard	<ul style="list-style-type: none"> <li>•Small-scale Culture of Microbe such as E. coli</li> <li>•Various incubations such as Enzyme reaction</li> <li>•Hybridization</li> </ul>
P.120	Personal-11-EXN set	•High Temp. accuracy	
P.120	Personal-11-SM set	•Power saving, Low price	
P.121	MD-1218	•Monode kit for Personal-11	•Shake culture using L-shaped Test tube.
P.122	Personal Lt-10F-SX set	•Low-temperature type	<ul style="list-style-type: none"> <li>•Cultivation of microbe such as Psychrophilic bacteria</li> <li>•Various incubations such as Enzyme reaction</li> </ul>
P.122	Personal H-10-SH set	•High-temperature type	<ul style="list-style-type: none"> <li>•Cultivation of microbe such as Thermophilic bacteria</li> <li>•Various incubations such as Enzyme reaction</li> </ul>

## "The long seller" Large Integrated type shaking constant temp. bath

Page	Model	Features	Applications
P.124	MM-10	•Standard	<ul style="list-style-type: none"> <li>•Cultivation of Microbe such as E. coli</li> <li>•Various incubations such as Enzyme reaction</li> </ul>
P.124	ML-10F	•Low-temperature type	<ul style="list-style-type: none"> <li>•Cultivation of Microbe such as E. coli</li> <li>•Various incubations such as Enzyme reaction</li> <li>•Ames test [ML-10F with PU-6, Some Modification required]</li> </ul>

## Add Shake to Unit Water bath at a low price

Page	Model	Features	Applications
P.126	EP-1	•Shaking unit	<ul style="list-style-type: none"> <li>•Cultivation of Microbe such as E. coli</li> <li>•Various incubations such as Enzyme reaction</li> <li>•Hybridization</li> </ul>

## Simple cooling singly. Low temp. constant temp. bath by combination.

Page	Model	Features	Applications
P.128	80LF	<ul style="list-style-type: none"> <li>•Just throwing the tip of cooling pipe in the object for use.</li> <li>•Cooling pipe made of stainless steel and has movable flexibility.</li> <li>•Only Cooling function without Temperature control function.</li> </ul>	<ul style="list-style-type: none"> <li>•Combined with Constant temperature Water bath</li> <li>•Cooling for Trap container and Reaction container.</li> <li>•Cooling for Samples</li> </ul>
P.129	150LF		
P.129	400R		
P.129	250DF		
P.129	800R		

•Enable to achieve the max (min) temp. by optional combination and them limited conversely (see left figure). "(Substantial) operational temp. ranges show the product singly and temp. range that can be achieved by optional selection.

(\*1) The operational temp ranges of product singly are shown when the included plastic water bath C-type (heat resistant temp: 70°C) used. (SDN/ EXN-B come with Heat insulation bath B-type and heat resistant temp not problem. 70°C cannot be made due to the heat radiation when hood not used). (\*2) SDminiN cannot be combined with coolers such as Cool pipe. Use it inside a low temp room around 4°C when low temp required. Use antifreeze solution when the temp is below 7°C for other models. Please note not recommended to operate at high temp with antifreeze liquid as it is.

(\*3)SH-10N compatible with heat medium for high temp. Be sure to use it above 100°C (recommended from above 70°C). As other models not compatible with heat medium for high temp, the upper limit 90°C to 100°C even when heat insulation water bath used. •The higher the heater output, the faster the heating speed.

	Settable temp. range	(Realistic)Temperature range			Temp. control accuracy	Heater	Page
		Thermominder only (With the standard accessory Water bath)(*1)	Used with Thermal Insulation water bath and Hood or Lid	Lower limit when immersion cooler added shown in the left (*2).			
	•-20°C to 80°C	•5°C above Room temperature to 70°C	•5°C above Room temp. to 80°C	•Around +10°C	•±0.1°C to	•500W	P.114
	•-20°C to 100°C		•5°C above Room temp. to about 95°C		•±0.1 to 0.3°C	•800W	
	•-20°C to 100°C		•5°C above Room temp. to about 95°C		•±0.02 to 0.08°C	•800W	P.115
	•-20°C to 85°C		•5°C above Room temp. to about 85°C		•±0.1°C to	•500W	
	•-20°C to 95°C		•5°C above Room temp. to about 90°C		•±0.1°C to	•700W	P.116
	•-20°C to 105°C		•5°C above Room temp. to about 100°C		•±0.05°C to	•1kW	
	•-20°C to 180°C		•5°C above Room temp. to 180°C(*3)		•±0.1°C to	•1kW	
	•-20°C to 100°C		•5°C above Room temp. to about 95°C		•±0.1°C to	•1.2kW	

	Shaking method	Temperature range	Temp. control accuracy	Cooling function	Page
	•Reciprocal	•5°C above Room temperature to 70°C	•±0.1 to 0.3°C		P.120
			•±0.02 to 0.08°C		
			•±0.1°C to		
	•Monode	—	—	—	P.121
	•Reciprocal	•-10°C to 50°C	•±0.1°C to	✓	P.122
		•5°C above Room temperature to 180°C	•±0.1°C to		


	Shaking method	Temperature range	Temp. control accuracy	Cooling function	Page
	•Reciprocal	•5°C above Room temperature to 80°C	•±0.0 to 0.1°C		P.124
		•0°C to 50°C	•±0.05 to 0.2°C	✓	

	Shaking method	Temperature range	Temp. control accuracy	Cooling function	Page
	•Reciprocal	※Shaking unit only (temp. control unit and water bath required separately)			P.126

	Temperature range	Cooling capacity	Page
	•-10°C to 30°C	•Approx. 150W	P.128
	•-15°C to 30°C	•Approx. 290W	
	•-30°C to 30°C	•Approx. 370W	
	•-45°C to 30°C	•Approx. 130W	P.129
	•-75°C to 0°C	•Approx. 150W	

#### [Immersion cooler] Cool pipe series

Being easily a low temp. bath by the combination of Thermominder. Recommended when water bath size want to be determined freely.



#### Spring net shaking platform in which various vessels can be mounted.

Stainless steel Spring net shaking platform makes it easy to install both either test tubes and flasks. Moreover, each step height including the bottom plate can be changed according to the vessel sizes (except for Monode shaking type).



#### Selection guide

Constant temperature incubator shaker  
CO<sub>2</sub> incubator shaker

Shaker

Mixer  
Rotator  
Stirrer

Bead beater  
homogenizer  
Ultrasonic homogenizer

Aluminum block bath  
Minimize bath

Water bath  
Shaking water bath  
Immersion cooler

Hybridization oven  
Constant temperature chamber

Centrifugal concentrator  
Cold trap

Freeze dryer






















Electrophoresis and blotting apparatus

Constant temperature circulating water circulating system [Chiller]

Appendix




## Correspondence table of Unit Water bath "Thermminder" series

Released in 1982	1993 - 1996	after 1997
<p>Released in 1982</p> <p><b>Mini-80</b> 10°C to 80°C Popular edition</p> 	<p><b>SM-05</b> 0°C to 100°C Popular edition</p> 	<p><b>SDmini</b> -20°C to 80°C Dedicated for Plastic Water bath</p>
<p>Released in 1980</p> <p><b>Ace-80</b> 0°C to 80°C Temp. control with Electromagnetic Contactor</p> 	<p>Discontinued(SX-10N below is recommended, while it</p>	
<p>Released in 1988</p> <p><b>Jr-100</b> 0°C to 100°C Standard</p> 	<p><b>SJ-10</b> 0°C to 100°C Standard</p> 	
<p>Released in 1988</p> <p><b>Dx-100</b> 0°C to 100°C High Specs</p> 	<p><b>DX-10</b> -20°C to 100°C High Specs, Stirring force adjustment</p> 	
<p>Released in 1988</p> <p><b>Lt-100</b> -50°C to 50°C Enables Low temp. setting, High Specs</p> 	<p>Discontinued(All of Current models can be set from -20 °C.)</p>	
<p>Released in 1988</p> <p><b>G-100</b> 0°C to 100°C Remote setting, High Specs</p> 	<p><b>DG-10</b> -20°C to 100°C Remote setting, High Specs</p> 	
<p>Released in 1988</p> <p><b>H-100</b> 0°C to 180°C High Temp</p> 	<p><b>DH-12</b> 0°C to 180°C High Temp</p> 	<p><b>SH-12</b> 0°C to 180°C High Temp</p> 
<p>Released in 1980</p> <p><b>P-80</b> 10°C to 100°C External Circulation</p> 	<p><b>SP-12</b> 0°C to 100°C External Circulation</p> 	
<p>Released in 1988</p> <p><b>C-630</b> 0°C to 100°C External Circulation, High Specs, Deep water level</p> 		
<p>Released in 1988</p> <p><b>C-650</b> 0°C~+105°C External Circulation, High Specs, Deep water level</p> 		
	<p><b>SD-B</b> 0°C to 100°C Low water level, For attached water bath</p> 	<p><b>EX-B</b> 0°C to 100°C Low water level, High Specs, For attached water bath</p> 

- The models on the most right (model name shown in yellow) are currently being produced. Refer them to replace the old models with new ones.
- The color of bands of background are imaged with the body color of actual product. Refer to them when identify the model you have.

# from First-generation models to Latest models



•The temp. range shown are the "Settable temp. range". Immersion cooler etc. required for temp. control below RT (25°C).

\*Operational temp. range\* and \*Settable temp. range\* might be narrow due to a performance of corresponded optional accessories.

● Protuberances not included in Dimensions. ● Vessels of photo not included.

We contribute to the development of research and industry.  
[ 2019-2020 General Catalog ] **TATEC**

Selection guide

Constant temperature incubator shaker  
OD-MonitorCO<sub>2</sub> incubator shaker

Shaker

Mixer  
Rotator  
StirrerBead beater  
homogenizer  
ultrasonic homogenizerAluminum block bath  
Minimize bathWater bath  
Shaking water bath  
Immersion coolerHybridization oven  
Constant temperature chamberCentrifugal concentrator  
Cold trap

Freeze dryer

Electrophoresis and blotting apparatus

Constant temperature circulating water system [Chiller]

Appendix

# Thermominder SDminiN

**Even water level 45 mm from the bottom of bath can be used. With the included Plastic Water bath easy to use. Optional Thermal insulation Water bath can be combined.**

Optional Combination Examples --> P.119



Model	SDminiN
Temperature range (*1)	5°C above RT(*2) to 80°C
Settable temp. range (*3)	-20°C to 80°C
Temp. control accuracy (*4)	±0.1°C to
Stirring method in Bath	Jet flow
Temperature display	Digitally (Changeable Preset/Current value)
Timer	Buzzer notification for Preset time Operation OFF *Setting range: 1min to 99h59min (Both model)
Other functions	Buzzer notification when Preset temp. reached, Automatic tuning
Minimum Water level	45mm from the bottom of bath
Heater	500W (Time proportional output variable)
Safe devices/ protections	Circuit protector, Fuse, Dry-heating protection with float, Heater protection cover, Sensor error, Short circuit, High/Low temp.-sample protection, Water level alarm, Non-volatile memory error, Automatic tuning error, Alarm setting error
Power failure	Switchable Automatic/Manual recovery
Dimensions inside Bath	180(W)×237(D)×155(H)mm (Occupancy Dim.: 190W x 83D x 147Hmm)
Main unit dimensions	220(W)×127(D)×235(H)mm (When Bath C-type used: 220W x 398D x 262Hmm)
Weight	Approx. 3kg
Power supply/ consumption	AC100V/5.5A 75Wh at 37°C (46Wh when Hood used)(*5)
Standard accessories	1 x Plastic Water bath C-type (*2)

(\*1) Max. temperature might not be reached when optional Hood not used or depending on the usage conditions. As it is dangerous due to steam when using high temp. with tap water, recommended to use our specified Heat medium for High temp.(See 117 page).

(\*2) As Max. Heat-resistant temperature of bath C-type up to 70°C, use optional stainless steel insulated bath F-type above 70°C.

(\*3) Not equipped with cooling function. Used together with our Chiller/Coolnit, Immersion cooler etc. at below RT (25°C).

(\*4) The value under the conditions of RT (25°C), AC 100 V/50Hz, Capacity 6 liters water, Preset temp. 37°C and No heat load.

(\*5) The reference value when Plastic Water bath C-type and optional Heat-resistant Plastic Hood PF-SDM at each preset temp. used under the condition of \*4.


## Features

- Low water level 45mm suitable for small vessels
- Simple operation system with Operation OFF Timer
- Optional Heat-resistant Plastic Hood and Thermal insulation Water bath available

## Applications

- Various incubations such as Enzyme reaction
- Temporarily Incubation of Culture medium
- Inactivation of Serum

## Dimensions SDminiN with Plastic Water bath C-type



## Optional accessories : Heat-resistant Plastic Hood

Heat-resistant Plastic Hood PF-SDM  
Plastic Water bath C  
Combination Examples



Description/Model	Remarks
<b>Heat-resistant Plastic Hood PF-SDM</b>	Heat-resistant temp. 100°C. Low power consumption with reduce evaporation. Dedicated for this SDminiN.

# Thermominder SDN-B/EXN-B

**Even water level 45 mm from the bottom of bath can be used. Thermal insulation tray included. EXN-B equipped with Temp. control accuracy with two decimal places and Timer function.**

Optional Combination Examples --> P.120

## Features

- Low water level 45mm suitable for small vessels
- Simple operation system with Temperature memory
- Comes with Stainless steel Insulation tray (Water bath)


## Applications

- Various incubations such as Enzyme reaction
- Temporarily Incubation of Culture medium
- Inactivation of Serum



EXN-B with Timer function. 3 steps simple program operation with temperature memory in combination.

## Dimensions (Common)



Model	SDN-B	EXN-B
Temperature range (*1)	5°C above RT to 70°C	
Settable temp. range (*2)	-20°C to +100°C	
Temp. control accuracy (*3)	±0.1°C to ±0.3°C	±0.02°C to ±0.08°C
Stirring method in Bath	Jet flow	Jet flow (4-level control)
Temperature display	Digitally (Changeable Preset/Current value)	
Temperature memory	3pc	
Timer	-	Buzzer notification for Preset time Operation OFF Each Setting range: 1min to 99h59min
Other functions	Buzzer notification when Preset temp. reached Automatic tuning	Buzzer notification when preset temp. reached Simple program (3 steps) Safety device output, Temp. input/output (*4) Automatic tuning
Minimum Water level	45mm from the bottom of bath	
Heater	800W	
Safe devices/ protections	Circuit protector, Fuse, Dry-heating protection with float, Sensor error, Short circuit, High/Low temp. sample protection, Water level alarm, Non-volatile memory error, Automatic tuning error, Alarm setting error	
Power failure	Switchable Automatic/Manual recovery	
Dimensions inside Bath	230(W)×350(D)×140(H)mm	
Main unit dimensions	280(W)×109(D)×215(H)mm (When Thermal insulation tray B-type used: 280W x 475D x 240Hmm)	
Weight	Approx. 3kg (6kg: Thermal insulation tray B-type included)	
Power supply/ consumption	AC100V/8.5A, 105Wh at 37°C	
Standard accessories	1 x Thermal insulation bath B-type (*5)	

## Optional accessories : Heat-resistant Plastic Hood

Heat-resistant Plastic Hood  
PF-B  
Combination Examples



Shaking Water bath  
Personal-11  
Combination Examples



Description/Model	Remarks
<b>Heat-resistant Plastic Hood PF-B</b>	Low power consumption with reduce evaporation. Designed for these SDN-B and EXN-B.
<b>Shaking Water bath Personal-11</b>	The water bath equipped with shaking function. See 120page for details.

(\*1) Max. temp. might not be reached when optional Hood not used or depending on the usage conditions. As it is dangerous due to steam when using high temp. with tap water, recommended to use our specified Heat medium for High temp. (See 117 page)

(\*2) Not equipped with cooling function. Used together with our Chiller/Coolnit, Immersion cooler etc. at below RT (25°C). Use our specified Antifreeze fluid (Thermal medium for Low temp. See 117page) when preset temp. below 7°C. The components life might be shortened when those used at the upper limit temp. that can be set. That of EXN-B is actual measured value due to (0.01 unit).

(\*3) The value under the conditions of RT (25°C), AC100V/50Hz, Capacity 6 liters water, Preset temp. 37°C and No heat load with Thermal insulation tray B-type.

(\*4) Optional Signal cable CA-671 required to output. Ask us for Temperature Input.

(\*5) Thermal insulation bath B-type sold singly.



# Thermominder SM-05N/SJ-07N/SX-10N/SH-10N/SP-12N

**Five models: Economy, Standard, High Specs, High Temp. and External Circulation. Select the most suitable model according to your purpose. Each Power consumption example is described for each models (See below).**

Optional Combination Examples --> P.119 Shaking Water bath "Personal-11 SM Set" --> P.120



**Economy** SM-05N

**Standard** SJ-07N

**High Specs** SX-10N

**High Temp** SH-10N

**External Circulation** SP-12N

With Heat resistant  
Plastic Hood  
(sold separately)

## Features

- Comes with Plastic Water bath. Heat insulation Water bath available as option.
- Combined with optional cooler to be Low temp. Water bath.
- Corresponds to High temp. [SH-10N], External circulation [SP-12N]

## Applications

- Various incubations such as Enzyme reaction
- Temporarily Incubation of Culture medium and Serum
- Circulation to Capillary of Evaporator [SP-12N]

Type	Economy	Standard	High Specifications	High Temperature	External Circulation
Model	SM-05N	SJ-07N	SX-10N	SH-10N	SP-12N
Temperature range (*1)	5°C above RT to 85°C	5°C above RT to 95°C	5°C above RT to 105°C	5°C above RT to 180°C	5°C above RT to 100°C
Settable temp. range (*2)	-20°C to 85°C	-20°C to 95°C	-20°C to 105°C	-20°C to 180°C	-20°C to 100°C
Temp. control accuracy (*3)	±0.1°C to		±0.05°C to	±0.1°C to	
External Circulation volume	-				Max. 6.2 L/min (*4)
Stirring method in Bath	Jet flow (Weak)	Jet flow (Controllable)		Jet flow	Jet flow (Controllable)
Temperature memory	-		1pc	-	
Timer	Buzzer notification for Preset time, Operation OFF, Temp. memory interlock (= Temp. transition timer, only for SX-10N) Setting range: 1min to 99h59min (Each model)				
Other functions	Buzzer notification when preset temp. reached, Automatic tuning, Safety device output (SM-05N excluded. Alarm out cable AOC-2 required to output.)				
Min. Water level	55mm from the bottom of bath		80mm from the bottom of bath	70% to 80% of Bath water capacity maintained	
Heater	500W (Time proportional output variable)	700W (Time proportional output variable)	1000W (Time proportional output variable)	1000W (Time proportional output variable)	1200W (Time proportional output variable)
Safe devices/ protections	Circuit protector, Fuse, Dry-heating protection with float, Heater protection cover, Sensor error, Short circuit, High/Low temp. sample protection, Water level alarm, Non-volatile memory error, Automatic tuning error, Alarm setting error				
Power failure	Switchable Automatic/Manual recovery				
Dimensions inside Bath	130(W)×66(D)×135(H)mm			130(W)×75(D)×137(H)mm	130(W)×85(D)×145(H)mm
Main unit dimensions	130(W)×135(D)×304(H)mm			130(W)×165(D)×302(H)mm	130(W)×164(D)×315(H)mm
Weight	Approx. 3.4kg			Approx. 5kg	
Power supply/ consumption	AC100V/5.5A, 41Wh at 37°C (*5)	AC100V/7.5A, 41Wh at 37°C (*5)	AC100V/10.5A, 41Wh at 37°C (*5)	AC100V/11A, 81Wh at 37°C (*6)	AC100V/13A, 167Wh at 37°C (*6)
Standard accessories	1 x Microtube Floater, 1 x Plastic Water bath C-type (*7)				

(\*1) Max. temp. might not be reached when optional Hood not used or depending on the usage conditions. Recommended to use our specified Heat medium for High temp (See 117, 118 page) when using tap water at high temp. as dangerous due to the steam.

(\*2) Not equipped with cooling function. Used together with our Chiller/Coolnit, Immersion cooler etc. at below RT (25°C). Use our specified Antifreeze fluid (Heat medium for Low temp. See 118 page) when preset temp. below 7°C. The components life might be shortened when those used at the upper limit temp. that can be set.

(\*3) The value under the conditions of RT (25°C), AC100V/50Hz, Capacity 6 liters water, Preset temp. 37°C and No heat load. That of SX-10N is an actual measured value due to "0.01 unit".

(\*4) The value varies depending on the inner dia, length of hose and hydraulic head.

(\*5) The reference value when Stainless steel Thermal insulation Water bath F-type and Optional Hood PF-SDM used under the condition of \*3. The power consumption increases by approx. 20% to 50% when Optional Hood not used.

(\*6) The reference value when Stainless steel Thermal insulation Water bath F-type used under the condition of \*3.

(\*7) As Heat-resistant temp. of Plastic Water bath C-type is 70°C, use Optional Stainless steel Thermal insulation Water bath.



# Optional accessories ①

Optional Combination Examples --&gt;P.119

## The combination of Thermominders with Water baths and Hood/Lids


Water bath		Used with Water bath					Applicable Hood/Lids		Used with Water bath and Hood/Lids				
Description/Model	Min. Inner Dim. and Capacity (Water level 75%)	SM-05N	SJ-07N	SX-10N	SH-10N	SP-12N	Model/Product name		SM-05N	SJ-07N	SX-10N	SH-10N	SP-12N
Plastic Water bath A -type (*1)	333 x 533 x 200Hmin Approx. 28L	(*)2	✓	✓	✓	✓	—		—	—	—	—	—
Plastic Water bath B -type (*1)	295 x 450 x 160Hmin Approx. 17L	✓	✓	✓	✓	✓	—		—	—	—	—	—
Plastic Water bath C -type (*1)	180 x 320 x 155Hmin Approx. 7L	✓	✓	✓	✓	✓	Heat-resistant Plastic Hood PF-C		✓	✓	✓		
Stainless steel Thermal insulation Water bath A -type	300x 400 x 200Hmin Approx. 18L	✓	✓	✓	✓	✓	Stainless steel-made Flat lid A/B (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	✓	✓	Stainless steel-made Top lid A/B (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	✓	✓	Stainless steel-made Flat lid A/B (SH)					✓	
		✓	✓	✓	✓	✓	Stainless steel-made Top lid A/B (SH)					✓	
Stainless steel Thermal insulation Water bath B -type	300 x 400 x 150Hmin Approx. 14L	✓	✓	✓	✓	✓	Stainless steel-made Flat lid A/B (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	✓	✓	Stainless steel-made Top lid A/B (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	✓	✓	Stainless steel-made Flat lid A/B (SH)					✓	
		✓	✓	✓	✓	✓	Stainless steel-made Top lid A/B (SH)					✓	
Stainless steel Thermal insulation Water bath D	355 x 600 x 155Hmin Approx. 25L	✓	✓	✓	✓	✓	Stainless steel-made Flat lid D (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	✓	✓	Stainless steel-made Top lid D (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	(*)3	✓	Stainless steel-made Flat lid D (SH)					✓	
		✓	✓	✓	✓	✓	Stainless steel-made Top lid D (SH)					✓	
Stainless steel Thermal insulation Water bath E -type	300 x 500 x 155Hmin Approx. 17L	✓	✓	✓	✓	✓	Stainless steel-made Flat lid E (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	✓	✓	Stainless steel-made Top lid E (SM/SJ/SX)		✓	✓	✓		
		✓	✓	✓	✓	✓	Stainless steel-made Flat lid E (SH)					✓	
		✓	✓	✓	✓	✓	Stainless steel-made Top lid E (SH)					✓	
Stainless steel Thermal insulation Water bath F -type (*4)	192 x 330 x 155Hmin Approx. 7L	✓	✓	✓	✓	✓	Heat-resistant Plastic Hood PF-C		✓	✓	✓		

(\*1) Heat-resistant temperature 70°C (\*2) Temp. reaching time to 37°C or more might be slow due to the heater output capacity to water bath volume. (\*3) Using together Optional stirring unit recommended (130 page) to stir the sample fully in the water bath. (\*4) Comes with Drain cock (Outer dia. 14mm for -10°C to 80°C). Other drain cocks required when used in other temp. range.


•Stainless steel-made Flat Lid/Top lid made to order. •Lid of water bath for SN-12 customized to order. Ask us for details.



### Temperature rising time (25°C --> 70°C)




### SM-05N/SJ-07N/SX-10N/SH-10N/SP-12N



# Optional accessories ②(SM-05N / SJ-07N / SX-10N / SH-10N / SP-12N)

"ThermominderSM-05N/SJ-07N/SX-10N/SH-10N/SP-12N"-->P.116

## Spring net stand

	Model	Dimensions	Applicable Water bath
	A2	270×270×140Hmm	Stainless steel Thermal insulation Water bath A/B, Plastic Water bath A/B
	B2	280×320×140Hmm	
	C	140×240×140Hmm	Stainless steel Thermal insulation Water bathA/B/F (*Stainless steel Thermal insulation Water bathF can't be used with Type-c and SP-12N sets.) Plastic Water bath A/B/C Thermal insulation TrayB/For SDN-B/EXN-B,See page 129 for details.)

## Immersion cooler

Enables temp. control below RT (25°C). by throwing Cool Pipe into water bath.

Description/Model	Cool Pipe 80LF/ 150LF
Applications	Makes Thermominder a Low temp. constant bath.
Main Specs	Temperature range:-10/-15°C to +30°C Cooling capacity :about 150/290W(When Liquid temp and Electrical frequency is 10°C and 50Hz.)

●80LF shown ●See page 142 for details.

## Pipe clamp

Fixes cooling pipe, pipe of Immersion cooler to water bath.

Product	Pipe clamp
Applications	Fixes cooling pipe, pipe of Immersion cooler to water bath.
Clamping width	Adjustable for 5.35 mm or 25.60 mm by replacing the pressing plate according to the thickness of water bath.

## Cooling Pipe

Enables temp. control below RT (25°C). by installing it in water bath to circulate cooling water.

Description/Model	Cooling Pipe A / B
Applications	Temp. control below RT (25°C) of Thermominder
Main Specs	A:For maintaining around 5°C B:For maintaining RT (25°C)

●Installation example of Water bath and Thermominder.

## Energy saving effect of Lid

Enable power consumption saving by covering the lid onto Unit water bath.



Temp.(°C)	Power Consumption (Wh) for SM/SJ/SX	
	Without Lid	With Lid
37	80	41
55	184	109
70	343	182

### USER'S VOICE

Power consumption reduces by nearly 40% just by covering the lid onto Unit.

## Heat medium for Low temp./Heat medium for High temp.

Use it if necessary when Immersion cooler used together or an operation at high temp.

Product / Model	Remarks
<b>Showbrine blue</b>	Heat medium for Low temp.(Antifreeze) 20kg, One (1) Can. Recommended for temp. below 7°C.
<b>Silicone oil MA-50</b>	Heat medium for High temp.18kg, One (1) Can. Kinetic viscosity50mm <sup>2</sup> /s(at 25°C), Focuses Temp. accuracy, Recommended for temp. above 70°C
<b>Silicone oil MA-100</b>	Heat medium for High temp.18kg, One (1) Can. Kinetic viscosity100mm <sup>2</sup> /s(at 25°C), Focuses Low evaporation, Recommended for temp. above 70°C

●See page 199 for details on Heat medium for Low temp. and page 200 for details on Heat medium for High temp.

## Physical properties of Silicone oil

Silicone oil is excellent in thermal and oxidation stability and is flame retardant with high flash points. As the vapor pressure very low and the amount of evaporation loss small it does not contaminate the work environment and suitable for long-time operation at high temp. Has feature that can be easily being warm and cool it as its specific heat is about One third of water.


### USER'S VOICE

Thermal conductivity about One forth compared to water but Eight times as much as air so it can be expected being more effective than air bath.



## Automatic water supply Unit : Level Keeper UB-2

Supplies water automatically to the preset water level when water level drops.



Water level can be set from 40 to 75 mm from upper surface of bath by adjusting the float.

Use a pressure-resistant hose and secure it with the included hose band so that the hose does not come off under water pressure.

Model	UB-2
Max. Operable differential pressure	Below 1.0MPa (10kg/cm <sup>2</sup> )
Heat-resistant temp.	-10 to +70°C (Water supply temp. Up to +40°C)
Control method	Water supply: Solenoid valve, Water level: Float
Dimensions inside Bath	70 x35mm(Clampable width within 30mm)
Water level adjustment width/Water supply port	40 to 75mm(Upper surface of bath), Outer dimensions10mm/14mm
Power supply	AC100V / 0.1A
Standard accessories	Hose band x1

# Optional Combination Examples of Thermominder

Please select a combination of Thermominder, Water bath, Lid, Cooler, Heat medium etc. according to the target temperature range from Low temp. to High temp.


"ThermominderSDminiN" -->P.114 "ThermominderSM-05N/SJ-07N/SX-10N/SH-10N/SP-12N" -->P.116 "Cool Pipe150LF"-->P.128



I want energy saving and compact one as much as possible as I often use at 37°C!!

## Combination Examples (1)


Configuration
①Unit Water bath / Thermominder / SDminiN
②Plastic Water bath C <i>*Included in SDminiN</i>
③Heat-resistant Plastic Hood / PF-C



I want to use an one for Food quality inspection ( 44.5°C ) !!

## Combination Examples (2)



Configuration
①Unit Water bath / Thermominder / SJ-07N
②Stainless steel Thermal insulation Water bath A Approx. 18L (Water level 75%)
③Stainless steel-made Top lid A/B



I want to use as much capacity in the bath as possible while keeping water temp. around 0°C!!

## Combination Examples (3)


Configuration
①Unit Water bath / Thermominder / SX-10N
②Immersion cooler / Cool Pipe / 150LF
③Stainless steel Thermal insulation Water bath B Approx. 14L (Water level 75%)
④Pipe clamp
⑤Heat medium for Low temp. /Showbrine blue



I want to automatically supply water to prevent the bath to boil dry when operating at 37°C for a long time!!

## Combination Examples (4)

Configuration
①Unit Water bath / Thermominder / SM-05N
②Automatic Water supply Unit Level Keeper / UB-2
③Stainless steel Thermal insulation Water bath B Approx. 14L (Water level 75%)



Selection  
guide

Constant  
temperature  
incubator  
shaker  
CO<sub>2</sub>-Monitor

CO<sub>2</sub> incubator  
shaker

Shaker

Mixer  
Rotator  
Stirrer

Bead beater  
homogenizer  
Ultrasonic  
homogenizer

Aluminum  
block bath  
Minimize bath

Water bath  
Staking water bath  
Immersion cooler

Hybridization oven  
Constant temperature  
chamber

Centrifugal  
concentrator  
Cold trap

Freeze dryer

Electrophoresis  
and  
Blotting apparatus

Constant  
temperature  
water circulating  
system [Chiller]

Appendix

# Personal-11 SDN/EXN/SM Set

**Versatile Bench-top Water bath equipped with Reciprocal shaking. Unique Monode shaking platform for shaking L-shaped test tubes available as an option.**

"Thermominder SDN/EXN-B" --> P.115 "Thermominder SM-05N" --> P.116

※The lid cannot be used when Monode Kit combined.

※The lid cannot be used when Monode Kit combined.



Personal-11 SDN Set



Personal-11 EXN Set Combined with Monode Kit (Option)



Personal-11 SM Set

## Features

- Thermominder and Personal (Shaking water bath) come as a set.
- Easy to remove included shaking platform and inner tank.
- Reciprocal shaking, Shaking width adjustment, Monode shaking (Option).

## Applications

- Small-scale Culture of Microbe such as E. coli
- Various incubations such as Enzyme reaction
- Hybridization

Model	Personal-11 SDN Set Thermominder SDN comes as a set	Personal-11 EXN Set Thermominder EXN comes as a set	Personal-11 SM Set Thermominder SM comes as a set (*1)
Thermominder	Temperature range (*2)	5°C above RT to 70°C (Preset/Current temp. displayed simultaneously in SDN/EXN.Switch. Switchable Preset/Current temp display in SM.)	
	Settable temp. range (*3)	-20°C to 100°C	
	Temp. control accuracy (*4)	±0.1°C to 0.3°C	±0.02°C to 0.08°C
	Temperature memory	3	-
	Timer (*5)	-	Buzzer notification for Preset time, Operation OFF/Operation ON
	Other functions	Buzzer notification when preset temp. reached, Automatic tuning	Buzzer notification when preset temp. reached, Automatic tuning
	Heater	800W	500W
	Safe devices/ protections	Circuit protector, Fuse, Dry-heating protection with float, Heater protection cover (SM-05N), Sensor error, Short circuit, High/Low temp. Sample protection, Water level alarm, Non-volatile memory error, Automatic tuning error, Alarm setting error, Switchable Automatic/Manual recovery when power failure	
	Power supply	AC100V/8.5A	
	Standard accessories	-	
Personal	Shaking method	Reciprocal shaking, Shaking speed: 20 to 160r/min (displayed digitally), Shaking width: 10 to 40mm (Stepless variable, Default 30mm)	
	Platform dimensions	220×310mm	
	Other functions	1 x Elapsed time indicator (0.1 to 999.9h, with Automatic reset), 1 x Service outlet for Thermominder	
	Bath inside dim./volume	235(W)×430(D)×140(H)mm, Approx. 11L (80% Water level) (*8)	
	Power supply	AC100V/0.5A (9A: with SDN)	AC100V/0.5A (6A: with SM-05)
	Standard accessories	1 x Spring net Shaking Platform, 1 x Clamp Spacer	
Dimensions/weight		300(W)×495(D)×336(H)mm, Approx. 19kg	300(W)×527(D)×420(H)mm, Approx. 19.4kg

(\*1) Can be combined with SJ-07N and SX-10N.

(\*2) Max. temp. might not be reached when optional Hood not used or depending on the usage conditions.

(\*3) Max. usable temp. in Personal-11 is 70°C. The components life might be shortened when the unit used above 70°C. Both Thermominder and Personal not equipped with cooling function. Used together with Personal Lt-10F (122 page) at below RT (25°C).

(\*4) The value under the conditions of RT (25°C), AC100V/50Hz, 80% Water level, Preset temp. 37°C and No heat load. Only EXN Measured value.

(\*5) Each setting range is 1min to 99h59min.

(\*6) Optional Signal cable CA-671 required to output. Ask us for Temperature Input.

(\*7) This also serves as a packing box whereas not used together with Personal-11 Set.


(\*8) The proper volume is 7L to 9L when shaking with the capacity (max. number) such as Centrifuge tubes and max. 500mL of Erlenmeyer flasks at 120 to 160r/min.



## Optional accessories

If you want larger capacity (Num. of vessels) of Monode Shaking.-->P.124 If you want Monode shaking in Air bath-->P.134, 135

Product configuration~Each part easily detachable~



### Capacity of Vessels Example

Spring net shaking platform (Spring pitch15mm)	φ11mm Test tube (vertical)	187
	Centrifuge tube (30°C tilted) 50ml	12
	Erlenmeyer flask 100mL	8
	Erlenmeyer flask 200mL	6
	Erlenmeyer flask 250mL	5
	Erlenmeyer flask 300mL	5
Monode kit (Corresponds to SDN/EXN set)	Erlenmeyer flask 500mL	3
	Optional L-shaped Test tube (See below)	12

### L-shaped Test tube and Monod Shaking


Monod shaking is the most effective shaking motion in culture using L-shaped Test tubes. We offer L-shaped Test tubes with the dimensions on the left and Shaking platform suitable for it.



### Option for Personal-11 SDN/EXN Set



Heat-resistant Plastic Hood PF-P Example for use



Monod kit MD-1218

### Option for Personal-11 SM Set



Combination example of Stainless steel Roof lid and Asbestos Timer unit.

### Related products: WTB-Shaker




Further compact-size "WTB-ShakerUnit" that is CFC-free and can cool. --> P.121

Description/Model	Remarks
Heat-resistant Plastic Hood / PF-P	Suppresses evaporation and Reduce power consumption. It can not be used when Monod kit MD-1218 below installed. Used by replacing it with Monode Shaking platform for 12pcs x L-shaped Test tube and Spring net Shaking platform of Personal-11.
Monode kit / MD-1218(*)	Used by replacing it with Monod Shaking platform for 12pcs x L-shaped Test tube and Spring net Shaking platform of Personal-11.
L-shaped Test tube (incl. 10)	φ18 x 120x70mm


(\*) Suppresses evaporation and Reduce power consumption.

Description/Model	Remarks
Stainless steel-made Top lid	Suppresses evaporation and Reduce power consumption.
Automatic Water supply Unit Level Keeper / UB-2	Supplies water when water level in the bath drops below the preset.
Cooling pipe B-type	For Cooling water circulation. Enables it use below RT .
Asbestos Timer unit B-type	For automatically turning Shaking ON and OFF. For Preparation of samples for X-ray analysis of Asbestos.

### Personal-11 SDN/EXN Set, MD-1218



### Personal-11 SM set



● Protuberances not included in Dimensions. ● Vessels of photo not included.

We contribute to the development of research and industry.  
[ 2019-2020 General Catalog ] **TATEC**

Selection  
guide

Constant  
temperature  
incubator shaker  
OD-Monitor

CO<sub>2</sub> incubator  
shaker

Shaker

Mixer  
Rotator  
Stirrer

Bead beater  
homogenizer  
ultrasonic  
homogenizer

Aluminum  
block bath  
Minimize bath

Water bath  
Shaking water bath  
Immersion cooler

Hybridization oven  
Constant temperature  
chamber

Centrifugal  
concentrator  
Cold trap

Freeze dryer

Electrophoresis  
and  
Blotting apparatus

Constant  
temperature  
water circulating  
system [Chiller]

Appendix

# Personal Lt-10F-SX Set/Personal H-10-SH Set

**Versatile Bench-top Water bath equipped with Reciprocal shaking. 2 models for Low temp. and for High temp. For Cultivation of Particular Microbe and Preparation of Analysis samples.**

"Thermominder SX-10N/SH-10N" --> P.116



Personal H-10-SH set

Personal Lt-10F-SX Set


## Features

- Benchtop Shaking water bath with -10°C to 80°C.
- Reciprocal shaking, Shaking width adjustable
- Silicone oil can be used for High temp. type

## Applications

- Cultivation of Psychrophilic bacteria and Thermophilic bacteria
- Various incubations such as Enzyme reaction
- Preparation of X-ray analysis samples of Asbestos (Option)

### Falling temperature time (20°C --> -10°C)




Personal Lt-10F-SX Set

Ambient temperature: 25°C  
Liquid volume: 11L Methanol (\*)  
Result: Approx. 60 minutes

(\*) Not recommended to use as antifreeze as it is dangerous due to flammable.

### Rising temperature time (30°C --> 70°C)



Personal Lt-10F-SH Set

Ambient temperature: 25°C  
Liquid volume: 11L Water (\*)  
Result: Approx. 60 minutes

Model		Personal Lt-10F-SX Set Personal Lt-10F and Thermominder SDN come as a set	Personal H-10-SH Set Personal H-10 and Thermominder EXN come as a set
Thermominder	Temperature range (*1)	-10°C to 50°C	5°C above RT to 180°C
	Settable temp. range (*2)	-20°C to 105°C (5°C to 105°C in SX-10N used singly)	-20°C to 180°C
	Temp. control accuracy (*3)	±0.1°C to (±0.05°C in SX-10N used singly)	±0.1 to
	Temperature memory	1	-
	Timer (*4)	Buzzer notification for Preset time, Operation OFF, Temperature transition	Buzzer notification for Preset time, Operation OFF/Operation OFF
	Other functions	Buzzer notification when preset temp. reached, Automatic tuning, Safety device output (Alarm out cable AOC-2 required to output.)	
	Heater	1000W (Time proportional output variable)	
	Safe devices/protections	Circuit protector, Fuse, Dry-heating protection with float, Heater protection cover, Sensor error, Short circuit, High/Low temp. Sample protection, Water level alarm, Non-volatile memory error, Automatic tuning error, Alarm setting error, Switchable Automatic/Manual recovery	
	Power supply	AC100V/10.5A	AC100V/11A
	Standard accessories	1 x Microtube Floater, 1 x Plastic Water bath C-type (*5)	
Personal	Shaking motion	Reciprocal shaking	
	Shaking speed/width	20 to 160r/min (Digitally display), 10 to 40mm (Stepless variable, Default 30mm)	20 to 160r/min (Digitally display), 10 to 30mm (Stepless variable, Default 30mm)
	Platform dimensions	220x310mm	
	Other functions	1 x Elapsed time indicator (0.1 to 999.9h, with Automatic reset), 1 x Service outlet for Thermominder	
	Compressor	75W	-
	Bath inside dim./volume	235(W)×430(D)×140(H)mm, Approx. 11L (80% Water level, Approx. 10L: Combined with Thermominder or Platform) (*6) In Personal Lt-10F-SX Set the water bath cannot be detached and equipped with drain hole (unit rear). In Personal H-10-SH Set the water bath can be detached. and W/D drain.	
	Power supply	AC100V/4.5A (15A: with SX-10N)	AC100V/0.5A (11.5A: with SH-10N)
	Standard accessories	1 x Dedicated Spring net Shaking platform	
Dimensions/weight		381(W)×545(D)×559(H)mm, Approx. 30kg	336(W)×575(D)×457(H)mm, Approx. 22kg

(\*1) Max. temp. might not be reached when optional Hood not used or depending on the usage conditions.

(\*2) Use antifreeze when below 7°C and heat medium for high temp. when above 70°C. See the right page for our specified heat medium.

(\*3) The value under the conditions of RT (25°C), AC100V/50Hz, 80% Water level, Preset temp. 37°C and No heat load.

(\*4) Each setting range is 1min to 99h59min.

(\*5) This also serves as a packing box whereas not used together with Personal-11 Set.

(\*6) The proper volume is 7L to 9L when shaking with the capacity (max. number) such as Centrifuge tubes and max. 500mL of Erlenmeyer flasks at 120 to 160r/min.

## Optional accessories

### Each Detachable part



Personal Lt-10F is equipped with a drain hole as the water bath inside it cannot be removed due to its construction (Water bath inside Personal Lt-10F and Drain cock on the back are shown in left figure). Water bath inside Personal H-10 can be removed.

### Capacity of Vessels Example

Spring net shaking platform (Spring pitch 15mm)	φ11mm Test tube (vertical)	187
	Centrifuge tube (30°C tilted) 50ml	12
	Erlenmeyer flask 100mL	8
	Erlenmeyer flask 200mL	6
	Erlenmeyer flask 250mL	5
	Erlenmeyer flask 300mL	5
	Erlenmeyer flask 500mL	3


### Optional accessories

#### USER'S VOICE

The lid plays an important role for energy saving when at High or Low temp. operated!!



Personal Lt-10F SX set  
Stainless steel-made Top lid  
+Asbestos Timer unit  
Combination Examples



#### USER'S VOICE

High temp. type also useful for culturing Hyperthermophiles slants (Slant medium).




Personal H-10 SH set  
Stainless steel-made Top lid  
Combination Examples




Personal Lt-10F equipped with Service outlet on the lower left of back and Personal H-10 equipped with that on the lower right of back.

Description/Model	Remarks
Stainless steel-made Top lid	Suppresses evaporation and Reduce power consumption.
Automatic water supply unit : Level Keeper / UB-2	Supplies water when water level in the bath drops below the preset.
Cool Pipe B	For Cooling water circulation. Enables it use below RT (25°C).
Asbestos Timer unit B	For automatically turning Shaking ON and OFF. For Preparation of samples for X-ray analysis of Asbestos (for Lt-10F SX Set)
Heat medium for Low temp.(Antifreeze) / Showbrine blue ( 20kg )	Use it when below 7°C.
Heat medium for High temp. Silicone oil / MA-50 ( 18kg )	Kinetic viscosity 50mm <sup>2</sup> /s(at 25°C), Focuses Temp. accuracy, Recommended for temp. above 70°C
Heat medium for High temp. Silicone oil / MA-100 (18kg)	Kinetic viscosity 100mm <sup>2</sup> /s(at 25°C), Focuses Low evaporation, Recommended for temp. above 70°C

### Personal Lt-10F SX set



### Personal H-10 SH set



# Water bath Shaker MM-10/Cool bath Shaker ML-10F

**Shaking baths with High temp. accuracy that used in various testing and research fields. Corresponds to Ames test with Monod shaking and Program operation as option.**

MM-10



## Features

- Low-temperature type available. Draining easily.
- Reciprocal shaking, Shaking width adjustment, Monod shaking as option
- Possible for Program operation as option

## Applications

- Cultivation of Microbe such as E. coli
- Various incubations such as Enzyme reaction
- Ames test [ML-10F with PU-6, Some Modification required]



ML-10F


Model	MM-10	ML-10F
Temperature range (*1)	5°C above RT to 80°C	0°C to 50°C
Temp. control accuracy (*2)	±0.02°C to 0.1°C	±0.05°C to 0.2°C
Shaking motion/Speed range/Amplitude	Reciprocal shake, 20 to 160r/min, 10 to 40mm (Stepless variable)	
Temperature display	Digitally (Changeable Preset/Current value)	
Platform dimensions	400 x 300mm	
Stirring method in Bath	Jet flow	
Other functions	Temperature checking monitor. Remote temperature setting terminal (0V to 5 V input, Enables Temp. program control with optional Program unit PU-5 etc.). (*3) Drain hole (right side of the unit).	
Heater/Compressor	Heater: 1300W	Heater: 800W, Compressor: 125W
Safe devices/protections	Earth Leakage Circuit Breaker, Sample protection (High temp.), Water level alarm, Sensor error	
Bath inside dim./volume	520(W)×365(D)×225(H)mm, Approx. 25L (60% Water level)	530(W)×365(D)×225(H)mm, Approx. 35L (80% Water level)
Dimensions/Weight	810(W)×476(D)×316(H)mm, Approx. 45kg	983(W)×486(D)×338(H)mm, Approx. 62kg
Power supply	AC100V/1.5A	
Standard accessories	1 x Dedicated Spring net Shaking platform, 1 x Drain hole filter	

(\*1) Max. temp. might not be reached when optional Hood not used or depending on the usage conditions. Use heat medium for high temp. when above 70°C. See the right page for our specified heat medium.

(\*2) The value under the conditions of RT (25°C), AC100V/50Hz, 60 to 80% Water level, Preset temp. 37°C. That is an actual measured value due to "0.01 unit".

(\*3) Can be corresponded to Program unit PU-6 that enables the program control for temp. and shaking by a modification of this unit. Helps streamline Ames test.

## Dimensions





## Optional accessories

### ● Capacity of Vessels in included Spring net shaking platform

Vessels	Capacity
φ16mm Test tube (vertical)	204
Erlenmeyer flask 50mL	24
Erlenmeyer flask 100mL	20
Erlenmeyer flask 200mL	12
Erlenmeyer flask 250mL / 300mL	9
Erlenmeyer flask 500mL	6
Erlenmeyer flask 1L	3
Erlenmeyer flask 2L	2
Sakaguchi flask 500mL	6


### ● Mountable number of Clamps (Option)



Above 300mL Clamps comes with Spring and Above 500 mL comes with Octagonal rubber sheet.

Vessels	Model	Number
Erlenmeyer flask	50mL	CF-0050
	100mL	CF-0100
	200mL	CF-0200
	250mL	CF-0250
	300mL	CF-0300
	500mL	CF-0500
	1L	CF-1000
	2L	CF-2000
Sakaguchi flask	500mL	SF-0500

### ● Other Optional accessories



Monode kit MD-3018



MD-3018 Example for use



Program Unit PU-6



Program Unit PU-5

Description/Model	Remarks
Stainless steel-made Top lid for M series	Suppresses evaporation and Reduce power consumption.
Monode kit MD-3018	Monode Shaking platform below of 30pcs x L-shaped. Enables to adjust the angle of shaking by replacing it with Spring net Shaking platform.
L-shaped Test tube ( incl. 10 )	φ18 x120x70mm
Program Unit PU-5	Enables Program control of Temp.
Program Unit PU-6	Enables Program control of Temp. and Shaking (Some processing required separately for use in combination).
Heat medium for Low temp.(Antifreeze) / Showbrine blue (20kg)	Use it when below 7°C.
Heat medium for High temp. Silicone oil MA-50	Kinetic viscosity 50mm <sup>2</sup> /s(at 25°C), Focuses Temp. accuracy, Recommended for temp. above 70°C
Heat medium for High temp. Silicone oil MA-100	Kinetic viscosity 100mm <sup>2</sup> /s(at 25°C), Focuses Low evaporation, Recommended for temp. above 70°C

## Application examples in Ames test

### MM-10/ML-10F optimum for Preculture of Microbes.

Since Temp. control accuracy  $\pm 0.02$  to  $0.1^{\circ}\text{C}$  (ML-10F  $\pm 0.05$  to  $0.2^{\circ}\text{C}$ ) it enables Preculture of Microbes that meets GLP standard (= Keeps the temp.  $37^{\circ}\text{C}$  within  $\pm 0.5^{\circ}\text{C}$ ).

### Increases efficiency using with Program Unit!

Enables Automatic execution of Shaking ON-OFF and Temp. transition from storage temp ( $4^{\circ}\text{C}$ ) to the culture finished. e.g. It can be set before going home so that you can conduct the experiment immediately after the next morning.

#### USER'S VOICE

The combined PU-6 and ML-10F with some processing uses for Ames test conveniently.




ML-10F



PU-6

#### USER'S VOICE

Various other programs can be apparently set up.



### Energy-saving Constant temp. Chamber recommended for Main Culture.

Use "Invitro box iB-130" (page 138) for "Sprinkling it on Min Glucose Agar plate medium  $37^{\circ}\text{C}$  after Preincubation of Specimen liquid and Microbes by mixing" and "Incubation for 48 hours". This unit also contributes to energy saving in laboratory facilities by Three-position temp. control method.



#### USER'S VOICE

Since air volume for air circulation in the chamber can be arbitrarily adjusted recommended to weaken air volume if drying of agar plate is concerned.



## Plus Shaker EP-1

**Combined with Thermomixer becomes a Shaking Water bath at low cost. Possible for combined with Plastic Water bath A/B-type and Stainless steel Heat Thermal insulation Water bath E type.**

"Thermomixer SM-05N/SJ-07N/SX-10N/SH-10N/SP-12N" --> P.116



Model	EP-1
Shaking motion	Reciprocal shaking
Shaking speed/width	20 to 160r/min, 0 to 30mm (Stepless variable)
Platform dimensions	275 x 351mm
Applicable Water tank	Plastic Water bath A/B-type Stainless steel Heat Thermal insulation Water bath E type.
Dimensions inside Bath	280(W)×385(D)×140(H)mm
Clampable thickness	Below 50mm
Dimensions/Weight	330(W)×575(D)×250(H)mm, Approx. 14kg
Power supply	AC100V/1A

•Use this unit below 80°C (Up to 70°C when used together with Plastic Water bath A/B-type)

### Capacity

Vessel	Size	Spring net Shaking platform	When Clamps (Sold separately) used
Test tube	φ11mm	221pcs	-
	100mL	12pcs	13pcs
	200mL	8pcs	11pcs
Erlenmeyer flasks	250mL	6pcs	8pcs
	300mL	6pcs	8pcs
	500mL	4pcs	6pcs

### Optional parts: Holding leg, Clamps, Combination Shaking Water bath


Model	Remarks
<b>Holding leg for Plus Shaker</b>	Used when not attached to water bath, 1pc
<b>Clamp CF-0100</b>	For 100mL Erlenmeyer flasks, 1pc
<b>Clamp CF-0200</b>	For 200mL Erlenmeyer flasks, 1pc
<b>Clamp CF-0250</b>	For 250mL Erlenmeyer flasks, 1pc
<b>Clamp CF-0300</b>	For 300mL Erlenmeyer flasks, 1pc
<b>Clamp CF-0500</b>	For 500mL Erlenmeyer flasks, 1pc
<b>Plastic Water bath A-type</b>	Inner dim. 333 x 533 x 200Hmm, Up to 70°C
<b>Plastic Water bath B-type</b>	Inner dim. 295 x 450 x 160Hmm, Up to 70°C
<b>Stainless steel Heat Thermal insulation Water bath E type</b>	Inner dim. 300 x 500 x 155Hmm

### Features


- Combined with Water bath for use
- Can be used singly with Optional holding leg.
- Comes with Spring net shaking platform.
- Can be used with Clamps together.

### Applications

- Cultivation of Microbe such as E. coli
- Various incubations such as Enzyme reaction
- Hybridization




### Dimensions



"ThermomixerSM-05N/SJ-07N/SX-10N/SH-10N"-->P.116



EP-1 + SM-05N + Plastic Water bath A	EP-1 + SJ-07N + Plastic Water bath A
Temperature range: 5°C above RT to 70°C	Temperature range: 5°C above RT to 70°C



EP-1 + SX-10N + Plastic Water bath A	EP-1 + SH-10N + Stainless steel Thermal insulation Water bath E
Temperature range: 5°C above RT to 70°C	Temperature range: 5°C above RT to 80°C

# Program Unit PU-5/PU-6

**Enables Temperature transition and Program operation of Shaking ON/OFF with our Water baths, Incubator shaker and Chillers for Open circuit in combination.**

"Bioshaker BR-300LF/3000LF Series" --> P.024 to 027/032 to 033 "Cool bath Shaker ML-10F" --> P.124 "Coolnit CL Series" --> P.171

## Features

- PU-5 for Temperature Program
- PU-6 for Temperature and Shaking ON/OFF Program
- Enables Remote control (wired) for the target product

## Applications

- Program operation for Shaking water bath ML-10F
- Program operation for Incubator shaker BR-300LF
- Program operation for Chillers for Open circuit CL-80R etc.

## Conforming products for Program Unit



"Bioshaker" BR-300LF/3000LF Series --> P.032



"Chillers for Open circuit CL-80R" --> P.171




"Cool bath Shaker ML-10F" --> P.124

## USER'S VOICE

Even if our other units without program function this unit enables them to do Program operation. Although it is wired it can be operated from a distance.



## Dimensions



● Protuberances not included in Dimensions. ● Vessels of photo not included.




## Example of Temperature program pattern

PU-6 enables Automatic operation (\*) combining Shaking ON/OFF with such as Temperature program below.


1. Repeats the reaction temperature and storage temperature.

Trapezoidal pattern




2. Temperature condition as Day and Night assumed.

Trapezoidal pattern



3. Stepwise temperature rise and fall and holding.

Multistage trapezoidal pattern




(\*)Some processing is required separately for connection to perform automatic operation of Shaking ON/OFF with PU-6. Ask us for details.

Model	PU-5	PU-6
Temperature program	Number of Storage pattern	2
	Number of Storage segment	8 segments/patterns
	Max. Number of segment	16
	Number of Pattern connection	2
	Number of Repeat	1 to 999 or Infinite
	Preset time range	0 to 99h59min
	Functions	Weight zone, Holding, Step
Time signal	Output	DC 0V to 5V
	Number of Storage pattern	-
	Number of Storage	1
	ON-OFF: Twice (2) per each	
Dimensions	Preset time range	0 to 99h59min
	Output	AC100V/Max. 15A (Resistance load)
	Dimensions	223(W)×181(D)×159(H)mm
Standard accessories	1 x Pt Temperature Sensor, 1 x Connecting cable	

# Cool Pipe 80LF/150LF/400R

**Just throwing the tip of cooling pipe in the object for use. For making Thermominder a Low temp. constant bath and Cooling for trapping vessels for Concentrator.**

Unit water bath "Thermominder" in combination --> P.119 Related products such as Temperature control units --> P.130




Model	80LF	150LF	400R
Temperature range (*1)	-10°C to 30°C	-15°C to 30°C	-30°C to 30°C
Cooling capacity (*2)	Approx. 150W	Approx. 290W	Approx. 370W
Condenser (Air-cooled) output	80W	150W	400W
Cooling Pipe Structure	Stainless steel Flexible tube		
Pipe Immersing part dim.	φ16mm×330mm	φ16mm×500mm	φ16mm×1000mm
Pipe Thermal insulation length	1m		
Unit Dimensions (W x D x H)	270×260×350mm	310×410×320mm	370×530×330mm
Weight	Approx. 21kg	Approx. 28kg	Approx. 36kg
Power supply	AC100V/2.5A	AC100V/4A	AC100V/6A

(\*1) Not equipped with Temperature control function. Use Thermominder etc. together with if necessary. The Min. temp. might not be reached depending on the ambient temperature.

(\*2) The value under the conditions of ambient temp. 25°C and liquid temp. 10°C at 50 Hz.

- When using a temperature controller the temp controller that controls the heater recommended to protect the refrigerator.
- Use a heater with capacity larger than cooling capacity of this product when using together Thermominder etc.
- Use an antifreeze that does not corrode Stainless steel, Chemically and Thermally stable and a viscosity of below 30mm<sup>2</sup>/s (specific gravity 1.0) in within operational temp. range
- Do not bend Cooling pipe extremely (Min bending radius 50 mm for fixed bending). Might crack and cause gas leakage if bent it forcibly. Do not immerse the heat insulation part of Cooling pipe (part that black insulation material is wound) in liquid.
- Cannot be used with Seawater.

## Dimensions



## Features


- Just throwing the tip of cooling pipe in the object for use.
- Cooling pipe made of stainless steel and has movable flexibly.
- Only Cooling function without Temperature control function.

## Applications


- Combined with Constant temperature Water bath
- Cooling for trapping vessels and Reaction container.
- Cooling for Samples

## Cooling curve and Cooling capacity curve

### Cooling curve



### Cooling capacity curve



## Optional accessories/Related products

See 130 page for details on products below

Model	Descriptions
<b>Temperature control Unit TU-100N</b>	Controls Immersion Heater: ±0.1°C to 0.5°C
<b>Temperature control Unit TU-200N</b>	Controls Cooling pipe and/or Immersion Heater: ±2.0°C
<b>Immersion Heater</b>	5 types available by material and capacity (W).
<b>Circulating pump Unit JP-40</b>	Circulates water in the bath to outside of the bath.
<b>Stirring Unit</b>	Stirring for small-scale water inside the bath.
<b>Heat medium for Low temp. Showbrine Blue</b>	20kg/Can; for below 7°C.



# Cool Pipe 250DF/800R

**Usability same as 80LF/150LF/400R while these enable Ultra Low temperature. Coolability of 250DF is to -45°C and that of 800R is to -75°C.**


Unit water bath "Thermominder" in combination --> P.119 Related products such as Temperature control units --> P.130

## Features

- Unitary refrigeration that to -45°C
- Cascade refrigeration that to -75°C
- Other features are compatible with 80LF/150LF/400R.


## Applications

- Combined with Constant temperature Water bath
- Cooling for trapping vessels and Reaction container.
- Cooling for Samples 800R




## Cooling curve and Cooling capacity curve

### Cooling curve



### Cooling capacity curve



Methanol was used for the experiment purpose so not recommended to use it in actual.

## Optional accessories/Related products

See 130 page for details on products below

Model	Description
<b>Temperature control Unit TU-100N</b>	Controls Immersion Heater: $\pm 0.1^\circ\text{C}$ to $0.5^\circ\text{C}$
<b>Temperature control Unit TU-200N</b>	Controls Cooling pipe and/or Immersion Heater: $\pm 2.0^\circ\text{C}$
<b>Immersion Heater</b>	5 Models available by material and capacity (W).
<b>Circulating pump Unit JP-40</b>	Circulates water in the bath to outside of the bath.
<b>Stirring Unit</b>	Stirring for small-scale water inside the bath.
<b>Heat medium for Low temp. Showbrine Blue</b>	20kg/Can; for below $7^\circ\text{C}$ . Note: This cannot be used at Min temp of both 250DF/800R due to its freezing point (Concentration 80% at around $-40^\circ\text{C}$ ).

Model	250DF	800R
Temperature range (*1)	-45°C to 30°C	-75°C to 0°C
Cooling capacity	Approx. 130W (*2)	Approx. 150W (*3)
Condenser (Air-cooled) output	250W	2 x 400W
Cooling Pipe structure	Stainless steel Flexible tube	
Pipe Immersing part dim.	$\phi 16\text{mm} \times 330\text{mm}$	$\phi 34\text{mm} \times 220\text{mm}$
Pipe Thermal insulation length	1m	1.2m
Unit Dimensions (W x D x H)/Weight	310x410x320mm, Approx. 30kg	420x500x560mm, Approx. 70kg
Power supply	AC100V/6A	AC100V/17A (*4)

(\*1) Not equipped with Temperature control function. Use Thermominder etc. together with if necessary. The Min. temp. might not be reached depending on the ambient temperature.

(\*2) The value under the conditions of ambient temp.  $25^\circ\text{C}$  and liquid temp.  $10^\circ\text{C}$  at 50 Hz.

(\*3) The value under the conditions of ambient temp.  $25^\circ\text{C}$  and liquid temp.  $-40^\circ\text{C}$  at 50 Hz. Do not use this product (Upper limit  $0^\circ\text{C}$ ) when setting temp. exceeding  $0^\circ\text{C}$  using Thermominder etc. together.

(\*4) Breaker connection of AC100V and 20A recommended. Other than Single phase 100V might cause failure.

• When using a temperature controller the temp controller that controls the heater recommended to protect the refrigerator.


• Use a heater with capacity larger than cooling capacity of this product when using together Thermominder etc.

• Use an antifreeze that does not corrode Stainless steel, Chemically and Thermally stable and a viscosity of below  $30\text{mm}^2/\text{s}$  (specific gravity 1.0) in within operational temp. range. Our antifreeze (Showbrine Blue) cannot be used at Min. temp. of both 250DF/800R so ask us for details.

• Do not bend Cooling pipe extremely (Min bending radius 50 mm for fixed bending). Might crack and cause gas leakage if bent it forcibly. Do not immerse the heat insulation part of Cooling pipe (part that black insulation material is wound) in liquid.


• Cannot be used with Seawater.

## Dimensions



# Temperature control Unit TU-100N/TU-200N

**Supports an equipment without temperature control function such as Immersion Cooler. Makes it a Low temp. water bath if use it together with Cool pipe.**



## Features

- Dedicated for Heating control. Precision temperature control [100N]
- Controls Heater and Refrigerator alternately. Rough Temp. control [200 N]
- Enables also Operation for Stirring Unit and Circulation Unit.

## Applications

- Precision control for Immersion Heater [100N]
- Roughly controls Immersion Heater and Cool pipe [200N]
- Keeps samples warm and cool inside a constant temperature bath.

Model	TU-100N	TU-200N
Temperature range	-100°C to 100°C	
Temp. control accuracy (*)	±0.1°C to 0.5°C	±2.0°C
Temperature controller	Type K Thermocouple (PID Control)	Type K Thermocouple (ON-OFF Control)
Temperature control	Heating	Heating/Cooling alternately
Temp. setting method	UP/DOWN Key	
Temperature display	Digitally	
Outlet for Heater	Max. 1200W	Max. 1200W (alternately with Compressor)
Outlet for Compressor	-	Max. 600W (alternately with Heater)
Outlet for Stirring unit	Max. 50W (Circulating Pump unit can also be connected)	
Safe devices/protections	Earth Leakage Circuit Breaker, High temperature	
Unit Dimensions (W x D x H)	180x235x90mm	
Power supply	AC100V/Max. 15A	
Standard accessories	Titanium protective tube temperature sensor	

(\*) The value under the conditions that the proper combination and arrangement taken at ambient temp. 25°C. Temp accuracy might be changed depending on ambient temp., preset temp., water bath capacity, with/without or strength of stirring, heat medium, heater, sensor and position of stirrer.

### Option: Immersion Heater Throwing this into the water bath for use



Copper made

Product
Immersion Heater (500W, Copper made)
Immersion Heater (1400W, Copper made)



Titanium made

Product
Immersion Heater (300W, Titanium made)
Immersion Heater (500W, Titanium made)
Immersion Heater (1000W, Titanium made)

•For heating sea water, titanium made immersion heater is recommended.

### Option: Stirring unit Stirring for small-scale water inside the bath.



Product	Stirring Unit
Motor	AC Motor (3W)
Clampable thickness	Maximum 50mm
Dimensions inside Bath	72 x 72 x 180Hmm
Unit Dimensions	72 x 138 x 255Hmm
Power supply	AC100V/0.5A

### Option: Circulating pump Unit Circulates water in the bath to outside of the bath.



Model	JP-40
Motor	AC Motor (40W)
Clampable thickness	Maximum 35mm
Circulating nozzle outer dia.	φ10mm
Circulating capacity (50/60Hz)	Max. Pump head: 4.5/5.7m Max. Flow rate: 10/11.5mL
Dimensions inside Bath	123 x 108 x 150Hmm
Unit Dimensions	123 x 221 x 289Hmm
Power supply	AC100V/1A