

## Application data sheet #04

Constant temperature  
incubator shaker

# DWMax M·BR-034P

## DNA extraction from cultivated yeast by shaking with deep well plates



### M·BR-034P is used for DNA extraction from cultivated yeast by shaking with deep well plates.

#### Overview

Thanks to the newly developed NewMax drive mechanism, the DWMax M·BR-034P has become a constant temperature incubator shaker equipped with a mixing capability by shaking, as well as uniformity and stability that exceeds the Maximizer M·BR-022UP, a similar product that is also currently sold. This Product can uniformly and strongly mix micro tubes and a 96-hole deep well. It is capable of exerting culture efficiency that is sufficient for microorganisms having a large cell body that tend to sink like yeast in a deep well which is generally unsuitable for shake-cultivating.


This document reports on the results of the following experiment. M·BR-034P was used to cultivate yeast (*Saccharomyces cerevisiae* S288C) in a deep well. A DNA extraction kit (Nippon Gene) was then used to extract the DNA from cultivated yeast. As a result, we were able to sufficiently collect the DNA as compared to a culture using a test-tube type culture tube.

#### Equipment used, reagents, etc.

- Constant temperature incubator shaker: DWMax M·BR-034P (TAITEC), BioShaker BR-23FP-MR (TAITEC)
- Deep well: 96-well deep well plate 2.0 mL, rectangular hole (AxyGen Scientific)
- Reagent kit that can extract DNA from yeast: ISOPLANT (Nippon Gene)
- Yeast: *S. cerevisiae* S288C
- Culture medium: YPD
- Centrifugal machine: Can support 12,000g (for a micro tube)
- Test tube mixer (Vortex): Delta mixer Se-04 (TAITEC)
- Reagent and equipment for agarose gel electrophoresis; ultraviolet visible spectrophotometric

#### Method

##### 1) Shake-cultivation of yeast using a deep well




2mL of a culture medium inoculated with yeast was dispensed to positions 1 to 5 in the deep well as shown in the illustration on the left, and then shake-cultivated with the DWMax M·BR-034P (1,500 r/min., +30° C, 24 hours).

As a subject for comparison, 2 mL of an inoculated culture medium was dispensed in a 12 mL disposable culture tube and was then shake-cultivated with a BioShaker BR-23FP-MR (reciprocal 180 r/min., +30° C, 24 hours). This is referred to as (6).

##### 2) DNA was extracted in accordance with the ISOPLANT protocol.

900  $\mu$ L was taken from the above culture solutions (1) to (6), each of which was put into a micro tube and centrifuged, with supernatant removed.



Centrifuging (1) to (6). With visual inspection, the amounts of precipitated fungi were almost the same.

↓  
Each was added with Solution I, and then vortexed.

↓  
Each was added with Solution II, and then vortexed.

↓  
Microtubes (1) to (6) were statically placed in the M-BR-034P and incubated at 50° C for 25 minutes.  
(Protocol specifies 15 minutes for incubation, but as the M-BR-034P in a gas phase was used, incubation was performed for a little longer. Light inverted agitation by hand every five minutes.)

↓  
Each was added with Solution III, and then vortexed.

↓  
Each was incubated in ice for 15 minutes.

↓  
Each was centrifuged at 12,000g for 10 minutes.

↓  
The water phase of each was taken, added with two times that amount of ethanol, and properly mixed and centrifuged at 12,000g for 10 minutes.

↓  
After removing supernatant, sediments were cleansed with 70% ethanol and air-dried.


↓  
The sediment was dissolved into 30 μL TE as a DNA solution, while purity and yield were checked by spectrophotometry and electrophoresis.

### Results and discussion

The yield, purity, and results of the electrophoresis of the extracted DNA are shown below. Compared to yeast cultivated in a culture tube (6), we were able to favorably collect DNA from yeast cultivated in a deep well (1) to (5). In each case, there was a good yield of 1.3 to 2 μg/mg.

	Yield (μg/mg)	A260/280
① (Cultivated in a deep well with M·BR-034P)	1.345	1.855
② ( // )	2.020	1.906
③ ( // )	1.475	1.916
④ ( // )	1.735	1.896
⑤ ( // )	1.470	1.909
⑥ (Cultivated in a culture tube with BR-23FP)	1.515	1.930

1% Agarose S electrophoresis ▶  
M : Marker 6 ( λ /Sty I digest, Nippon Gene)



**Now on sale**

Constant temperature incubator shaker  
DWMax M·BR-034P

**Supports validation!**



Power consumption:  
about **20 Wh**



With a new mechanism, the NewMax drive, high uniformity and stability are achieved in a multi-specimen culture. Four well plates (standard) and up to ninety-six 1.5 mL micro tubes (with an optional rack) can be mounted. Temperature can be adjusted between room temperature -7° C and +60° C with electronic heating and cooling.



### Written and edited by:

TAITEC Corporation  
2693-1 Nishikata, Koshigaya-shi, Saitama-ken, 343-0822, Japan  
TEL: +81-48-986-3228  
FAX: +81-48-988-8363  
E-mail: [overseas@taitec.org](mailto:overseas@taitec.org) Website: <https://e-taitec.com/>

Publication: 2017  
For inquiries about the constant incubator shaker DWMax M·BR-034P and the contents of this leaflet, please contact us as provided on the left.