

# Freeze Crusher $\mu$ T-48

**Powerful crushing of frozen samples with liquid nitrogen. 48 samples can be treated simultaneously. Optimum for extracting hard samples, proteins susceptible to heat denaturation, RNA, etc.**

•Example of various frozen crushed samples including inanimate objects. --> P.106-P.107

## Features

- Crushing of frozen sample in vessels with liquid nitrogen
- 2 mL Microtube or Dedicated metal container are used
- The throughput is 0.2 g to 2 g (Depends on the vessels)

## Applications

- Crushing of Yeast, Mold, Tissue piece of animals and plants
- Crushing of bones, teeth and limbs of small animals
- Crushing of Wire covering and Plastics, Asbestos samples, etc.



Vertically reciprocal shaking



**The procedure for freeze crushing** Work gloves are worn when freezing the sample with liquid nitrogen in the photo. Use gloves that are suitable for handling liquid nitrogen in actual use.



1 Put the specified amount (below) of sample and metal crusher into the vessel.




2 When a vessel requires a holder, attach it in the holder.



3 Soak and freeze in liquid nitrogen until bubbling stops.



4 Quickly attach to the  $\mu$ T-48 main unit.



5 Close the hood and check the time and shaking speed, and then start crushing.



6 After crushing, take out the crusher and proceed to the next process.

## Adapted to various samples with 3 types of vessel holders

Powerful crushing of frozen samples with liquid nitrogen. Living tissues and organs, Hard tissues such as bones, etc. and also some inanimate samples such as rubber and plastic can be crushed. Adapted to Marketed 2 mL tubes and Dedicated stainless steel crushing vessels.

## Optional accessories: Vessel holders



### Product Name / Model / Remarks

#### ① 48pcs-Holder for $\mu$ T-48 TH-0248T

1 pc of Holder (Capacity : 48 pcs of 2.0 mL round bottom Microtubes) and 100 pcs of Metal crushers come as a set. Sample throughput is 0.1 to 0.2 g/1 pc.

#### ② 3pcs-Holder for $\mu$ T-48 TH-0203T

4 pcs of Holders (Capacity : 3 pcs of 2.0 mL round bottom Microtubes), 24 pcs of Metal crushers and Rack come as a set. 4 holders can be attached to  $\mu$ T-48 (max. 12 tubes of 2.0 mL), which is superior to TH-0248T in heat retention.

#### ③ Stainless steel-made strong crushing vessel TH-SPT

Crushing vessel 4 pcs, Dedicated crusher and Rack come as a set. Suitable for samples that cannot be crushed by Microtube with Metal crusher. Larger amount of samples (1.0 to 2.0 g/1 pc) can be crushed than that of Microtubes.

## Recommended shaking speed in each crushing sample

- Stainless steel strong crushing vessel: Up to 1000 r/min
- Metal crusher: Up to 1200 r/min
- Crushing beads: Up to 1600 r/min

A shaking speed that is higher than the above speeds may cause breakage of tubes and vessels. Therefore, make sure to observe the shaking speeds above.

## USER'S VOICE

Very useful when extracting substances susceptible to denaturation and degradation of RNA and proteins.



Model	<b><math>\mu</math>T-48</b>
Crushing method	Vertically reciprocal shaking
Shaking speed	0 to 2500 r/min (*1)
Capacity	2.0 mL Microtube: Max. 48 pcs (*2)(*3) Stainless steel powerful crushing vessel : Max. 4 pcs (*3)
Timer	1 to 999 seconds
Safety devices/functions	Holder attachment detection switch, Cover opening detection switch
Dimensions (WxDxH)	220 x 310 x 315 mm
Weight	Approx. 10.0 kg
Power supply	AC100V/1A (Need a step-down transformer)

(\*1) Around 1200 r/min should be necessary and sufficient condition to crush the sample in actual use.

(\*2) Eppendorf "Safe-Lock Tube 2 mL" is recommended. (\*3) Microtube and Stainless powerful crushing vessel are available as an option.

NEW

Constant-temperature incubator/shaker OD Monitor

For cell culture related products

Shaker

Mixer Rotator Stirrer

Bead beater Ultrasonic homogenizer

Aluminum block Bath Mini-size Bath

Water bath Shaking Water bath Immersion cooler

Hybridization Incubator Constant-temperature Chambers

Centrifugal Concentrators Cold Trap

Freeze dryers

Substrate Electrophoresis apparatus Blotting device for hybridization

Constant-temperature water-circulating system [Chilled]

Appendix

# Example ① Freeze crushing of various samples including inanimate samples

**Embrittlement by freezing enables strong crushing.**  
**The freeze crushing with  $\mu$ T-48 is also suitable for Obligatory anaerobe samples.**




•Freeze Crusher  $\mu$ T-48 --> P.105

## Test results


- Freezing method Immerse the vessels with the sample and crusher into liquid nitrogen (2.0 mL: Vessel holder) and then freeze them.
- Rushing time 30 sec (Additional 30 sec if not completely crushed)
- Judgment whether sample is crushed Whether powder forms or nearly forms (Cut samples into any size that can be put in vessels).
- Vessels Safe-Lock tube 2.0 mL ...Marketed product (Made by Eppendorf)  
 Metal crusher.....Included in Optional 48 pcs-holder for  $\mu$ T-48 (used in this experiment)  
 Stainless steel-made strong crush vessels .....Optional parts (Dedicated crusher is included.)

<h3>Chicken thigh</h3> <p>Vessels: Safe-Lock tube 2.0 mL                  Sample volume: 0.1 g                  Shaking speed: 1200 r/min                  Crushed with: Metal crusher</p>	<h3>Human hair</h3> <p>Vessels: Safe-Lock tube 2.0 mL                  Sample volume: 0.1 g                  Shaking speed: 1200 r/min                  Crushed with: Metal crusher</p>	<h3>Human nails</h3> <p>Vessels: Safe-Lock tube 2.0 mL                  Sample volume: 0.2 g                  Shaking speed: 1200 r/min                  Crushed with: Metal crusher</p>
<h3>Mouse skin (with body hair)</h3> <p>Vessels: Safe-Lock tube 2.0 mL                  Sample volume: 0.2 g                  Shaking speed: 1200 r/min                  Crushed with: Metal crusher</p>	<h3>Mouse heart</h3> <p>Vessels: Safe-Lock tube 2.0 mL                  Sample volume: 0.2 g                  Shaking speed: 1200 r/min                  Crushed with: Metal crusher</p>	<h3>Mouse tail</h3> <p>Vessels: Stainless steel-made strong crush vessel                  Sample volume: 1 g                  Shaking speed: 1000 r/min                  Crushed with: Dedicated crusher</p>
<h3>Hypocotyl of Radish</h3> <p>Vessels: Safe-Lock tube 2.0 mL                  Sample volume: 0.2 g                  Shaking speed: 1200 r/min                  Crushed with: Metal crusher</p>	<h3>Okra seeds</h3> <p>Vessels: Safe-Lock tube 2.0 mL                  Sample volume: 2 pcs                  Shaking speed: 1200 r/min                  Crushed with: Metal crusher</p>	<h3>Hard rubber (Polychloroprene)</h3> <p>Vessels: Stainless steel-made strong crush vessel                  Sample volume: 2 g                  Shaking speed: 1000 r/min                  Crushed with: Dedicated crusher</p>

## Frozen sample/Holder temperature (Reference)



The 48 pcs-holder for 2 mL tube TH-0248T with 48 tubes and Metal crusher that was frozen with liquid nitrogen was attached to the unit. Then, the temperature change of each part while shaking at 1200 r/min was measured. The cryogenic temp. was completely kept for 30 to 60 sec which was required for crushing. The display temp. indicates the temp. stage of the top surface on which the holder is placed.



## Example ② Freeze crushing of various samples including inanimate samples


**Freeze crushing of Plastic samples using Freeze crusher  $\mu$ T-48 with Stainless steel-made strong crush vessel.**

• Freeze Crusher  $\mu$ T-48 --> P.105



### Results and Examination

We tried some crushing of samples such as polystyrene, polypropylene, and polycarbonate that are well known. Each result is as follows.




Polystyrene was able to be powdered completely (③). However, in case of the φ6 mm ball-shaped sample, there were large fragments that remained at a certain rate even after trying with different amounts and crushing times (①②). It seems that the ball-shaped sample remained uncrushed if it is stuck in upward of the crusher.

The result suggests that the shape of the sample is better to be like a chip shape (or a tablet shape) Polypropylene in order to be crushed into fine fragments. However, it was not crushed into powder (①). In order to improve (although the shaking speed limit is exceeded when using a strong crushing container), the shaking speed was performed at 1100 r/min for 5 minutes. It becomes fine, but like a braided piece of cotton (②). At this stage, it becomes difficult to collect unless suspended in a solvent.

Polycarbonate proved to be difficult to crush. Even if the shaking speed was reduced to 1100 r/min or by reducing sample amount, the result of this experiment was that only a small amount of powder was produced and the chip shape remained almost unchanged.

### Freeze crushing procedure when using stainless steel-made strong crushing vessel

An adjustable plier is useful for taking the Stainless steel-made strong crushing vessel (referred to as crushing vessel) in and out of the liquid nitrogen. Make sure to wear leather gloves when touching the frozen crushing vessel or the vessel holder that has become cold by contacting the frozen vessel. Make sure to ventilate the room well when using liquid nitrogen as there is risk to get Anoxia unknowingly because vaporized liquid nitrogen can become a huge volume of nitrogen gas.



\*1. Desirable to use the minimum-sized polystyrene foam container that the required number crushing vessels can be immersed to minimize the amount of liquid nitrogen used.

\*2. The processing capacity of the crushing vessel is 1.2 g per 1 pc while it is better to make it to 0.5 g per 1 pc for plastic samples (Up to 1 g polystyrene can easily be crushed by freezing).

\*3. Wait for at least 2 minutes after the boiling is settled out to freeze the sample in the crushing vessel sufficiently.

\*4. Place at least two "frozen" Stainless steel-made strong crushing vessels for balance and secure fixation. Because metal shrinks when frozen, both vessels must be frozen to ensure a firm fixation.

\*5. In December 2017, the rack was changed to a new type that does not require thumbscrews to secure vessels.

\*6. The crushed sample may stick to the crusher so tap it with the inner wall of the vessel to drop it.

\*7. If the crushing is insufficient, return the crusher to the unit to freeze it and shake again.