

NEW


Guide for aluminum blocks for different type of vessels

Aluminum Blocks for Microtubes

DTU-Mini CTU-Mini	B-0648	B-0835		B-1120A	B-1120B	B-1515
DTU-R CTU-R	B-0696/(B-0648 x 2)	B-0835 x 2		B-1140/(B-1120A x 2)	B-1120B x 2	B-1515 x 2
DTU-1B / 1C	 <p>B-0648 and B-0696 can also be used with 8-row Microtubes.</p>		AL-0836	AL-8136	AL-1136	AL-1136B
DTU-2B / 2C			AL-0836 x 2	AL-8136 x 2	AL-1136 x 2	AL-1136B x 2
Applicable vessels	① 0.1 mL · ② 0.2 mL	③ 0.5 mL	③ 0.5 mL + ④ 1.5 mL	④ 1.5 mL	⑤ 2.0 mL	④ 1.5 mL · ⑥ 5.0 mL




Aluminum block for Disposable centrifuge tubes




Aluminum blocks for Micro Vials



DTU/CTU-Mini	B-1220	AL-0030V	
DTU/CTU-R	AL-1220 x 2	 <p>For the following three Deep blocks, "Raising block" enables easily checking inside vials.</p>	
DTU-1B/1C	AL-1336V	AL-1525V	AL-2412V
DTU-2B/2C	AL-1336V x 2	AL-1525V x 2	AL-2412V x 2

⑩ 15 mL	⑪ 50 mL	⑪ 50 mL + ⑩ 15 mL
B-1612L		
B-1612L x 2		
AL-1616F	AL-2909F	AL-5011F
AL-1616F x 2	AL-2909F x 2	AL-5011F x 2

Aluminum blocks for Test tubes/Glass centrifuge tubes




Guide for Custom-made aluminum blocks

Even for vessels that cannot be used with standard aluminum blocks, we can craft Custom-made aluminum blocks with various numbers, diameters, depths, and bottom-shaped holes.

Example of Custom-made blocks


Aluminum block for 1 pc 500 mL square media bottle

Hole diameter: 75 mm × 75 mm
Hole depth: 60 mm
Hole bottom shape: Flat bottom
Capacity: 1 pc




Aluminum block with screw lid 2 mL Microtubes fit until its cap.

Hole diameter: 11 mm (Min. dia. in hole 10.5 mm)
Hole depth: 34 mm
Hole bottom shape: Drill hole
Capacity: 25 pcs



Aluminum block for square cells for the spectrophotometer

Hole shape: 12.5 mm × 12.5 mm
Hole depth: 35 mm
Hole bottom shape: Flat bottom
Capacity: 20 pcs



We also designed various blocks other than these. Please contact us when you desire custom-made blocks.

NEW

Constant
temperature
incubator shaker
OD Monitor

For cell culture
related products

Shaker

Mixer
Rotator
Stirrer

Bead beater
homogenizer
Ultrasonic
homogenizer

Aluminum
block Bath
Mini-size Bath

Water bath
Shaking Water bath
Immersion cooler

Hybridization
Incubator
Constant temperature
Chambers

Centrifugal
Concentrators
Cold Trap

Freeze dryers

Submarine
Electrophoresis apparatus
Blotting device for
hybridization

Constant-temperature
water circulating
system [Chiller]

Appendix